

# **EXD HME MicroCT Data Acquisition, Processing and Data Request Overview**

Isaac Seetho, William D. Brown, Harry E. Martz, Jr.  
Lawrence Livermore National Laboratory  
Livermore, CA 94551

Work performed on the  
Science & Technology Directorate of the  
Department of Homeland Security  
Statement of Work  
HSHQPM-10-X-00005 P00007

December 6, 2016  
LLNL-TR-714017



This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or Lawrence Livermore National Security, LLC, and shall not be used for advertising or product endorsement purposes.

This work performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344.

# **EXD HME MicroCT Data Acquisition, Processing and Data Request Overview**

Isaac Seetho, William D. Brown, Harry E. Martz, Jr.  
Lawrence Livermore National Laboratory  
Livermore, CA 94551

## **Executive Summary**

This document is a short summary of the steps required for MicroCT evaluation of a specimen. This includes data acquisition through image analysis, for the EXD HME program [1]. Expected outputs for each stage are provided. Data shall be shipped to LLNL as described herein.

## Data Acquisition

Data acquisition procedure is outlined in [2]. Data acquisition should be preceded daily by acquisition of the following full-panel uncollimated system calibration files for each system configuration to be used in scanning:

- Dark current detector output (*drk* file), with the source off.
- Light current detector output (*lit* file), with the source on at experiment-nominal kV and mA values.
- Mid current detector output (*mid* file), with the source on at nominal kV value and 2/3 of experiment-nominal mA values.

Acquisition with LLNL's DRCT software generates the following slit-collimated files:

- Gain-corrected radiographs for each view (\*.sdt and .spr)
- A *RawMode* folder with pre-gain correction detector output (\**R*\*.sdt and .spr), where *R* indicates a radiograph that has not had gain correction applied.
- Following the test plan for data collection, this is done for 3 system configurations:
  - 160kV 2-slit with AlCu filters (Exp1)
  - 100kV 2-slit with Al filter (Exp2)
  - 160kV 1-slit with AlCu filters (Exp1\_Open)

## Image Reconstruction

Image reconstruction is to be performed using Livermore Tomography Tools v1.3 (32b), using script files generated from a template provided by LLNL [3].

These script files:

- Perform bad pixel correction
- Generate attenuation radiographs (*attenRad*\*.sdt and .spr)
- Generate cropped attenuation radiographs (*attenRad*\*.sdt and .spr)
- Generate sinograms (*sinos*\*.sdt and .spr)
- Apply ring removal (if specified) and polynomial beam hardening compensation (for the 100 kV configuration)
- Generate reconstructed images (*reco*bj\*.sdt and .spr, *rz*\*.sdt and .spr)

During the reconstruction, LTT will generate a subfolder for each operation performed on the data, saving files into that subfolder. An example set of folders for a script that does bad pixel correction, generates attenuation radiographs, crops the radiographs, generates sinograms, performs polynomial beam hardening compensation, and then performs filtered backprojection reconstruction is presented below in Figure 1.




Figure 1: Example Reconstruction directory layout for a script performing the steps as described above. For a script with ring removal before pBHC, an additional folder would appear in the list.

For an acquired scan, the set of reconstructed images should include:

- 11 slices each in the upper and lower slits for two-slit experiments (*recobj\*.sdt* and *.spr*)
- 11 slices in the slit region for open experiments (*recobj\*.sdt* and *.spr*)
- Full-bottle reconstructions of the specimen for open experiments (*rz\*.sdt* and *.spr*)

## Data Reduction

LLNL is interested in the output of MicroCT Analysis v2.0, rev6 [4], which processes reconstructed images and pulls copper strip quality assurance data from gain-corrected images. The executable will generate image segmentation files as well as summary Excel spreadsheets (*\*\_characterization.xls* and *CuStrip\_\*.xls*)

Please report any errors encountered when running these executables to:  
Isaac Seetho (seetho3@llnl.gov).

## Additional Records

In addition to the acquired and reconstructed image files, LLNL is interested in an acquisition notes sheet per scan processed, following the format provided in Appendix A, and a data shipping Traveler document following the format provided in Appendix B.

## Requested Items to be Sent to LLNL (per specimen scanned)

Provided below is a list of files to send to LLNL for specimens scanned under DHS EXD. Files sent by disk shall be accompanied by a traveler document (see Appendix B) providing a complete inventory of disk contents. Disks shipments to LLNL shall also be accompanied by email notification to Steve Cocolis ([cocolis1@llnl.gov](mailto:cocolis1@llnl.gov)) and Harry Martz ([martz2@llnl.gov](mailto:martz2@llnl.gov)) with tracking information.

Send by Disk (to: Steve Cocolis, L-154, Lawrence Livermore National Laboratory, 7000 East Ave., Livermore, CA 94550):

- Files for each system configuration used (Exp1, Exp2, Exp1-Open)
  - Pre-acquisition system calibration files (*drk*, *mid*, *lit*)
  - Raw and gain-corrected DRCT output files (*\*R\*.sdt* and *\*.sdt*)
  - Attenuation radiographs (*attenRad\*.sdt*)
  - Sinograms (*sinos\*.sdt*)
  - Reconstructed images (*recobj\*.sdt*, *rz\*.sdt*)
  - Bad pixel files (*bdpix.txt*)
  - Beam hardening compensation files used (*\*\_ISA\_BHC\_Coeffs.txt*)
  - All parameter files (*\*.sct*) generated during data processing
  - Reconstruction scripts used
- MicroCT Analysis GUI output folders (for both specimen/reference image segmentation and for copper strip segmentation)

Send by Email (to: Isaac Seetho and Jeff Kallman at: [seetho3@llnl.gov](mailto:seetho3@llnl.gov), [kallman1@llnl.gov](mailto:kallman1@llnl.gov)):

- MicroCT Analysis GUI output summary spreadsheets
- Acquisition and reconstruction notes sheet

## References

1. William D. Brown, *TP83-MicroCT Data Acquisition, Reconstruction and Analysis Using the ISA MicroCT System*, LLNL-TR-649192, Lawrence Livermore National Laboratory, January 30, 2014.
2. Cary Pincus, Steven Benson, *MicroCT Data Acquisition Checklist*, LLNL-TR-XXXX, Lawrence Livermore National Laboratory, December 30, 2013.
3. William D. Brown, *Computed Tomography Reconstruction Checklist for Data Acquired on the ISA MicroCT System*, LLNL-TR-649121, Lawrence Livermore National Laboratory, January 23, 2014.
4. I. Seetho, W. Brown, J. Kallman, H. Martz, W. White, *MicroCT: Automated Analysis of CT Reconstructed Data of Home Made Explosive Materials Using the Matlab MicroCT Analysis GUI*, LLNL-TR-503291 rev 1, Lawrence Livermore National Laboratory, December 2016.

## Appendix A: MicroCT Acquisition and Reconstruction Notes

Provided below is a sample acquisition and reconstruction notes sheet. Each field should be filled out for each specimen scan performed.

### KVP and Filter Set:

1. 2-slit, Low Energy, 100 kvp, Aluminum 1.943mm
2. 2-slit, High Energy, 160 kvp, Aluminum 1.943mm, Copper 1.905mm
3. 1-slit, High Energy, 160 kvp, Aluminum 1.943mm, Copper 1.905mm

### Reference Standards:

Teflon Cylinder, 0.5" O.D. (outer diameter)

Water in 5 mL 0.5" O.D. test tube

Delrin Cylinder, 0.5" O.D.

Graphite Cylinder, 0.5" O.D.

Magnesium Cylinder, 0.5" O.D.

Aluminum Cylinder, 0.5" O.D.

### HME Containers:

50 mL 1" O.D. Polypropylene digestive vessel with green cap.

### Reconstruction Parameters:

SOD=1098.43 mm

SDD=1406.48 mm

ODD=308.05 mm

PXcenter=1009

PZcenter=663

Calibration Folders = C:\MCT\Calib\_140916\Cal

HME Code = Z1C4

Batch Code = N/A

Specimen ID Code = 141003\_Z1C4\_1

Prepared By = Smith

Preparation Date = 141002

Prescan Net Weight = 68.03

Prescan Volume = 62.99

Prescan Density = 1.08

Postscan Net Weight = 68.03

Postscan Volume = 62.99

Postscan Density = 1.08

Additional Notes: \*ROR development MCT specimen #1

## Appendix B: Data Shipment Traveler Document

<b>Data Tracking Receipt</b>
------------------------------

From:		Date:	
-------	--	-------	--

### Table of Contents

Compound Name:	
----------------	--

Additional Comments: (Total # of Drives Sent, # of CD/DVDs sent, etc...)	
---	--

Shipped By:		Date Shipped:	
UPS Tracking #		Delivery Date:	
Received By:		Date Received:	