

# **AN AUTOMATED INSTRUMENT FOR CONTROLLED-POTENTIAL COULOMETRY: SYSTEM DOCUMENTATION**

**M. K. Holland and J. V. Cordaro**

## **DISCLAIMER**

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

## **Approved by:**

**S. K. Formby, Superintendent  
Equipment Engineering Division**

**Publication Date: June 1988**

---

**E. I. du Pont de Nemours & Co.  
Savannah River Laboratory  
Aiken, SC 29808**

**PREPARED FOR THE U.S. DEPARTMENT OF ENERGY UNDER CONTRACT DE-AC09-76SR00001**

**MASTER**

**DISTRIBUTION OF THIS DOCUMENT IS UNLIMITED**

## **DISCLAIMER**

**This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.**

---

## **DISCLAIMER**

**Portions of this document may be illegible in electronic image products. Images are produced from the best available original document.**

DO NOT MICROFILM  
COVER

BLANK

## **ABSTRACT**

---

An automated controlled-potential coulometer has been developed at the Savannah River Plant for the determination of plutonium. Two such coulometers have been assembled, evaluated, and applied. The software is based upon the methodology used at the Savannah River Plant, however the system is applicable with minimal software modifications to any of the methodologies used throughout the nuclear industry.

These state-of-the-art coulometers feature electrical calibration of the integration system, background current corrections, and control-potential adjustment capabilities. Measurement precision within 0.1% has been demonstrated. The systems have also been successfully applied to the determination of pure neptunium solutions.

The design and documentation of the automated instrument are described herein. Each individual module's operation, wiring layout, and alignment are described. Interconnection of the modules and system calibration are discussed. A complete set of system prints and a list of associated parts are included.

Key words: Controlled-potential coulometer, plutonium, neptunium, automated, and electrical calibration.

BLANK

DO NOT MICROFILM  
COVER

## CONTENTS

---

	<u>Page</u>
ABSTRACT	3
INTRODUCTION	7
SUMMARY	8
DESCRIPTION	8
INSTRUMENTATION	10
SYSTEM SOFTWARE	26
OPERATION	29
RESULTS AND DISCUSSION	32
ACKNOWLEDGMENTS	46
REFERENCES	47
APPENDICES	
A. Drawing List	48
B. Vendor List	69
C. Parts List	71

BLANK

DO NOT MICROFILM  
COVER

## **AN AUTOMATED INSTRUMENT FOR CONTROLLED-POTENTIAL COULOMETRY: SYSTEM DOCUMENTATION**

---

### **INTRODUCTION**

Reliable controlled-potential coulometric determination of plutonium requires accurate electrolysis current integration, stable/low noise potential control, and easily operated instrumentation.

Commercially-available coulometers are capable of excellent potential control. These same coulometers have been interfaced to computer systems for instrument control, although the coulometers were not designed with automation as the primary intent. However, these systems are deficient in their ability to accurately integrate the electrolysis current. An analog integrator accumulates a charge on the integrating capacitor that is proportional to the quantity of material electrolyzed. This type of integrator deviates inherently from the ideal response by as much as -0.1% because of its inability to retain the accumulated charge quantitatively. The cause of this non-ideal response (integrating capacitor leakage and dielectric absorption) has been identified and mathematical corrections and procedural modifications have been developed to reduce the systematic errors.<sup>1</sup> However, a simpler solution is the use of a digital integrator. A digital integrator does not suffer from deviations associated with analog integrators. A digital integrator consists of a voltage-to-frequency converter to generate pulses (the sum of which is proportional to the quantity of material electrolyzed), a counter to accumulate these pulses (minimizing integrating capacitor losses), and a display which allows the value of the counter to be read. Although digital integrators have been incorporated in commercially-available coulometers, manufacturers' specifications on the integrating components preclude insuring the accuracy of such instruments beyond that obtained with analog integrators. Plutonium measurements obtained using the EG&G Princeton Applied Research digital coulometer at SRP were consistent with expectations based upon the manufacturer's specifications.

Thomas L. Frazzini and Michael K. Holland, while employed at the Department of Energy's New Brunswick Laboratory, NBL, developed a digital integrator using state-of-the-art electronic components in a patented circuit design,<sup>2,3</sup> and incorporated the integrator into the second and third generation automated coulometric measurement systems.<sup>4</sup>



The objective at SRP was to utilize documentation supplied by NBL to fabricate two systems comparable to the third generation NBL system and to ensure that the documentation would meet SRP QA requirements.

## SUMMARY

Two automated coulometers have been fabricated and installed at SRP. One of the systems is located in 772-F for routine sample measurements and the other system is in 772-1F to support development/evaluation activities and to provide backup for routine operations. During the course of construction and testing of these two systems, significant improvements in system operation were achieved. Both previously observed and unidentified anomalous system behaviors were identified and corrected. Documentation deficiencies were also identified and corrected.

Assembly, testing, evaluation, and documentation of these fourth-generation, state-of-the-art systems are described herein.

## DESCRIPTION

The SRP controlled-potential coulometer units are each composed of two potentiostat modules, an integrator module, a counter module, and an automation module. The units are controlled by a Hewlett-Packard computer. The computer, associated peripherals, and cell assembly are described in detail in a subsequent section of the report.

The potentiostat modules are based on the design of J. E. Harrar et al. using the best available electronic components. The units supply a stable control potential to the cell assembly where the species of interest is electrolyzed. Two potentiostats are used to separate the oxidation and reduction functions although the units are of identical design and are interchangeable. In both units, all of the electrolysis current passes thru a 50-ohm load resistor in series with the measurement cell. The voltage drop across the load resistor is thus proportional to the electrolysis current. When the voltage signal across the load resistor is integrated, the quantity of material electrolyzed may be determined from Faraday's Law.<sup>6</sup> Potentiostat #1, used for sample oxidation, is connected to the integrator module.

A programmable digital-to-analog-converter power supply is connected through a precision voltage divider to both of the potentiostat modules. The signal from the power supply is routed to the non-inverting input of the driver amplifier. Connection of the signal to this reference point of the driver amplifier allows

programmable incrementing of the control potential voltage supplied to the measurement cell. The selection of signal polarity in this system causes addition of the voltage to the control potential supplied by the reference supply. Addition of a signal to the control potential is a somewhat more convenient approach than a subtraction for automating the control potential adjustment<sup>6</sup> and the formal potential determination.

The integrator module is composed of two voltage-to-frequency converters (VFC's) and a quartz crystal oscillator. The signal from the potentiostat to be integrated and a small offset signal are supplied to VFC #1. Only the offset signal is supplied to VFC #2. The frequency signal from the two VFC's and the crystal oscillator are supplied to the counter module for accumulation on counters #1, #2, and #3, respectively. The signal to be integrated would be exponentially decaying DC were it not for AC noise and the adjustment of the control potential which cause changes in the polarity. The offset signal is supplied to ensure that the polarity of the net signal to be integrated is always positive. Supplying only the offset signal to VFC #2 allows accurate correction for the contribution of the offset signal to the total integrated counts from VFC #1. The integral of the signal from potentiostat #1,  $V_i$  is calculated as follows:

$$V_i = V_{i+o} - V_o \times K_o \quad (1)$$

where  $V_{i+o}$  is the signal plus offset counts (counter #1),  $V_o$  is the offset counts (counter #2) and  $K_o$  is the ratio of the responses (count/sec/V) of VFC #1 to VFC #2. This ratio, although near unity, is sufficiently different from one to be important to the accurate calculation of the integral. The value of  $K_o$  is easily and accurately determined when no current is flowing thru potentiostat #1 (standby mode). When the signal to be integrated is zero then, VFC #1 and VFC #2 both integrate only the offset signal. Thus,  $K_o = V_{i+o}/V_o$  when the signal is zero.

The potentiostat/integrator combination can be calibrated electrically using a precision resistor to replace the measurement cell. By using the potentiostat to apply a potential across the precision resistor and integrating for a known period of time, the coulombs and counts are obtained and their ratio may be determined as follows:

$$C = \frac{V_{Pr} \times T}{R \times V_i} \quad (2)$$

where  $V_{Pr}$  is the voltage drop across the precision resistor of value  $R$  in ohms,  $T$  is the calibration time in seconds, and  $V_i$  is the net counts calculated using Equation 1.

The theoretical value for the electrical calibration factor,  $C(\text{theoretical})$ , may be calculated from the design and component selection of the potentiostat and integrator.

$$C(\text{theoretical}) = \frac{1}{(R_x)(r)}$$

where  $R_x$  is the potentiostat #1 load resistance in ohms and  $r$  is the response of VFC #1 in counts/sec per volt of signal to be integrated.

For convenience, the value of the electrical calibration factor is expressed in coulombs per million counts.

The automation module allows interconnection/interfaces, data acquisition, and computer control for the peripherals, the cell assembly, and the coulometer modules. An overall system block diagram is shown in Fig. 1.

## **INSTRUMENTATION**

### **Computer and Accessories**

The computer and all accessories used were made by Hewlett Packard. The computer is a Series 200, Model 226. It comes with a 5 1/4 inch disk drive, 7 inch monochrome screen, 512 kbytes of internal RAM, and a HP-IB (IEEE-488) interface. In addition, the following accessories were installed;

- BCD (Binary Coded Decimal) Interface
- GPIO (General Purpose Input/Output) Interface
- additional 512 kbytes RAM

### **BCD Interface**

The BCD interface is used to read the value from the seven digit Counter module in binary coded decimal format. The BCD interface attaches to the Counter module through a customized cable (Ref print S5-2-15969). This interface inputs data only.




FIGURE 1. System Block Diagram

The switch settings on the BCD interface are as follows;

OF=ON	CTLB=OFF	SELECT CODE=01011 (ELEVEN)
DATA=ON	CTLB-2=OFF	INT=00
SIGN1=ON	DFLGB=OFF	
SIGN2=ON	CTLA=OFF	
OVL D=OFF	CTLA-2=OFF	
	DFLGA=OFF	

### **GPIO Interface**

The GPIO interface allows the computer to control the automation module which selects the coulometer operating mode and interconnects the coulometer, measurement cell assembly and peripheral devices via a series of mercury-wetted relays. For example, the computer can put the oxidation potentiostat in the standby mode while operating the reduction potentiostat and connecting it to the cell. The details of the capabilities of the GPIO interface are discussed fully in the automation module description. The GPIO interface connects to the automation module through a customized cable (Ref print S5-2-15969).

The switch settings on the GPIO interface are as follows;

PCTL=1	RD=1	SELECT CODE=01100 (TWELVE)
PFLG=0	BSY=1	INT=00
PSTS=1	RDY=0	
HSHK=1	RD=1	
DIN=1	BSY=1	
DOUT=0	RDY=0	

### **Peripherals**

The following peripherals, all manufactured by Hewlett Packard, are connected to the computer through the HP-IB interface:

- 3456A Digital Voltmeter
- 59501B Digital-to-Analog Converter/Power Supply
- 5316A Electronic Counter
- 8116A Pulse Function Generator
- 2673A Printer

### **3456A Digital Voltmeter**

The 3456A is a highly accurate 6 1/2 digit voltmeter. Full programmability over the HP-IB (address 705) is standard. The 3456A connects to the automation module output labeled DVM. Through the

GPIO interface, the computer connects the voltmeter to various points in the system. These points are detailed in the automation module description. All the settings on the voltmeter are made through software except for a front panel button which selects the rear terminals as the input.

### **59501B Digital-To-Analog Converter (DAC)**

The 59501B is a fully HPIB (address 707) programmable signal channel digital-to-analog converter, capable of supplying a 0-1 and 0-10 VDC in a unipolar mode. In the bipolar mode it is capable of supplying a -1 to +1 and a -10 to +10 VDC signal. The DAC is used in the 0-1 and 0-10 VDC range to automatically adjust the control potential on the potentiostat modules.

To reduce the DAC signal, the voltage is scaled down by a precision voltage divider (Ref print S5-2-15952D). The divider circuit reduces the input signal by a factor of five. The output of the divider network connects to the external potential input on both potentiostat modules (Ref print S5-2-15953) and from there directly to the non-inverting input of the driver amp AR2 (Ref print S5-2-15956). Note that the polarity of the DAC signal is positive at the non-inverting input of AR2 with respect to the amplifier common. When the output from the DAC is increased from zero, the control potential VCP, will increase (Ref Fig 2).

### **5316A Electronic Counter**

The 5316A is a highly accurate 7 digit counter which is fully programmable through the HPIB (address 708) interface. The counter was used in conjunction with the 8116A pulse generator for general timing of events. The counter provides accurate timing during electrical calibration and sample electrolysis. An additional application of the counter is realized during periodic alignment described in the digital integrator description.

The mode of operation of the counter is software controlled but the following front panel controls must be set:

Input=Channel A  
Trigger/Sensitivity Button=Trigger Position  
AC/DC button=DC position  
Filter/Norm button=Filter position  
Gate Time Delay Dial=Midpoint  
Level/Sense Dial=Midpoint  
Slope=positive  
Attn=x1 position  
SEP/COM A button=Sep position




FIGURE 2. Potentiostat Circuit

### **8116A Pulse Function Generator**

The 8116A generates pulses and waves with various duty cycles and frequencies. It is fully programmable through the HP-IB (address 709). The pulse generator was used with the 5316A counter for general timing purposes. All controls are programmed except for two buttons. The COMPL and the DISABLE LED's buttons must be manually deactivated.

### **2673A Graphics Printer**

The 2673A is fully programmable through the HP-IB (address 710). It is used to print all relevant data for a process.

### **MEASUREMENT CELL ASSEMBLY**

- 1.0 Princeton Applied Research Measurement Cell Assembly, Model 377A. The unit is composed of a platinum working and counter electrode, counter and reference bridge tubes with porous glass frits, a Teflon cell head, a Teflon cell retainer ring, a synchronous stirring motor, and a support stand. Sulfuric Acid (1.0N) is used in conjunction with the platinum working electrode.
- 2.0 A custom-made gold working electrode has been fabricated from gold mesh (punched/expanded metal ribbon) on a 16 gauge gold wire frame. The gold working electrode has also been used in the cell assembly. The supporting electrolyte used in conjunction with the gold electrode was 0.8N HNO<sub>3</sub>.
- 3.0 Reference Electrode. A miniature saturated calomel reference electrode is used. The Fisher Scientific catalog #13-639-79 has been found to be satisfactory.
- 4.0 Reverse-S Shaped Stirrer.<sup>7</sup> The PAR paddle stirrer was replaced with a comparable sized reverse-S shaped stirrer to prevent mechanical losses of the sample solution by splashing.
- 5.0 A Kel-F coated support rod replaces the manufacturer's support rod to prevent electrical paths from the cell to ground thru the support stand.



## POTENTIOSTAT MODULE

### Description

The potentiostat module's main purpose is to supply a constant potential to the cell through the working and reference electrodes. The potential is supplied by a 6.2 VDC floating supply. A 20 K ohm potentiometer is used to set the full scale span of the control potential. The voltage can then be adjusted over the span with the 2 K ohm Helipot (Ref print S5-2-15955). Two banana plugs on the rear panel allow the full-scale span voltage to be measured.

A three position front panel switch selects the three modes of operation through 5VDC relays. The three modes are; standby, hold, and operate. In the standby mode, the relays which connect the control circuit are disconnected from the 5 VDC supply leaving them in a de-energized state (Ref prints S5-2-15955, S5-2-15956). Note that when the computer is controlling the operation of the potentiostat, the mode switch is in the standby position and the relays can be energized by the automation module through the potentiostat connector J3. In the hold mode (see Fig 2), the counter and reference electrodes are connected to the circuit but the working electrode is still disconnected. Finally, in the operate mode the circuit is completed and the working electrode connects to the power supply common.

There are three remaining switches on the front panel which control the cell current monitor. The charge/normal switch supplies 4.8 VDC to the current monitor in the normal position. The 4.8 VDC is supplied by four 1.2 V Duracell batteries. In the charge position the monitor is off. The 200 mA/20 mA switch in the 20 mA position displays the cell current from 0 to 19.999 mA and the 200 mA position displays the current from 0 to 199.0 mA. The polarity shown on the current monitor for oxidation is negative while for reduction, positive. The display test switch illuminates the entire display.

On the rear of the potentiostat module are three coax connectors. The external potential input allows the user to input a voltage to the non-inverting input on the power amplifier to offset the reference supply (see Fig 2). The output current connector provides a voltage proportional to the cell current for use by the integrator. Note in Fig 2 the direction of  $I_{ox}$  (oxidation current). The voltage drop across the 50 ohm resistor is positive which is necessary since the integrator can only accept a positive voltage. The control potential output, used only during instrument alignment and troubleshooting, allows the user to measure the voltage between the working and reference electrodes.

The working, reference, and counter electrodes are connected directly to a three-pin Amphenol connector on the rear panel. The cell connections are also tied to the automation module through the NIM (Nuclear Instrument Module) connector for use in the computer controlled mode.

### **Alignment**

Alignment of the potentiostat module includes the internal null adjustments of both amplifier circuits and the span adjustment on the control potential circuit. Equipment needed for calibration includes:

- Hewlett Packard 3456 Voltmeter.
- BNC shorting plug.
- NIM Extender Cable with pins 34 and 42 not shorted.
- Reference prints S5-2-15956, S5-2-15955

#### **1.0 Initial Null Adjustments**

- 1.1 Turn the power switch on the automation Module off and switch the potentiostat to standby.
- 1.2 Remove the potentiostat to be calibrated from the NIM bin and remove the switching and load card.
- 1.3 Connect the NIM extender cable from J3 on the NIM bin to the NIM connector on the potentiostat module. Apply power to the system and allow for a 30 minute warmup.
- 1.4 Connect a short jumper from TP1 on the control amplifier board (Ref print S5-2-15956) to the ground plane of the board. Note: do not connect the jumper to the chassis ground; the control amplifier board ground floats.
- 1.5 Connect the negative lead from the 3456A Voltmeter to the ground plane and the positive lead to the output of amplifier AR1 (TP2).
- 1.6 Adjust R4 until the DVM reads less than .003 mvolts. Note: turn the potentiometer in a counter-clockwise direction to decrease the voltage.
- 1.7 Leave the non-inverting input of AR1 and the negative lead of the DVM grounded.
- 1.8 Attach the BNC shorting plug to the external potential connector on the rear panel.

- 1.9 Connect the positive lead of the DVM to the output of AR2 (TP3).
- 1.10 Short capacitor C1 with as short a wire as possible.
- 1.11 Adjust R5 until the DVM reads less than 1.0 mVolts.
- 1.12 This completes the initial offset adjustments for the Potentiostat module. Turn off the power to the system and remove all jumpers.
- 1.13 Re-assemble the Potentiostat Module.
- 2.0 Span Adjustment of the Control Potential
- 2.1 Connect a coax cable from the rear control potential jack on the Potentiostat module to the 3456A voltmeter.
- 2.2 Adjust the control potential dial on the front panel until it reads 0.
- 2.3 Adjust the span potentiometer on the rear panel until the voltmeter reads less than 0.10 volts.
- 2.4 Verify that when the front panel dial is reading full scale, the control potential output is 1 volt.


The accuracy of this adjustment is not crucial to the system's performance since the control potential is adjusted automatically by the computer during measurements.

## **DIGITAL INTEGRATOR MODULE**

### **Description**

The digital integrator module converts a voltage signal from the load resistor of the potentiostat module to a frequency. The voltage signal is proportional to the current passing through the cell. It is unique with respect to other integrators since it has two voltage-to-frequency convertors.

When the test/normal switch is in the normal position, the input potential from the potentiostat is sensed by the detection circuit (see Fig 3). Under normal conditions, the input will range from 0 to 3 VDC. If the voltage is negative, the sensing circuit will energize relay K1, disengaging the input from VFC #1, protecting it. The readout circuit accepts the 0-3 VDC input and converts it to a 1200 Hz to 31.200 kHz frequency signal. A 1.35 volt mercury battery in the automation module provides the 1200 Hz



**FIGURE 3. Digital Integrator Module**

offset. The zero circuit receives the same 1.35 volt signal from the automation module and converts it to a constant 1200 Hz offset. The benefit of the two VFC's and the offset voltage is discussed in the description.

The test/normal switch is put in the test position when calibrating and aligning the VFC's. In this position, the same voltage can be applied to both VFC's simultaneously.

Under normal operation, the outputs of VFC #1 and #2 connect to the counter module inputs, readout, and zero, respectively. The counter then converts the signals to a binary coded decimal format for display and input to the computer.

The 10 kHz crystal oscillator circuit is described in the Counter module description.

### **Alignment**

Alignment of the digital integrator includes setting the zero, span, and offset of the two VFC's. In addition, an internal protection circuit which protects the VFC's and a 10.000 kHz clock must be initially set. Equipment needed for calibration includes:

- Fluke 343A DC Calibrator.
- Hewlett Packard 5316A Electronic Counter.
- Hewlett Packard 3456A Voltmeter.
- Prints S5-2-15958, S5-2-15959.

#### **1.0 Alignment Procedure - READOUT**

- 1.1 Allow a 30-minute warmup before making adjustments.
- 1.2 Set the TEST/NORMAL switch to TEST.
- 1.3 Apply an input voltage of 10 mVolts to the rear panel BNC labeled CURRENT on the digital integrator module.
- 1.4 Monitor the front panel output BNC labeled READOUT with the Hewlett Packard Counter.
- 1.5 Adjust the front panel trimpot labeled Z1 for a reading of 100 Hz on the frequency counter (Note: a clockwise rotation of the trimpot will cause a decrease in frequency).
- 1.6 Apply an input voltage of 3.000 volts.

- 1.7 Adjust the front panel trimpot labeled FS1 for a reading of 30.000 kHz on the frequency counter.
- 1.8 Repeat steps 1.3 to 1.7 to determine repeatability and accuracy of adjustments made previously. Repeatability should be  $\pm 1$  count.
- 1.9 Set the TEST/NORMAL switch to NORMAL.
- 1.10 Attach a cable from the rear panel BNC labeled OFFSET VOLTAGE on the automation module to the rear panel BNC labeled OFFSET on the digital integrator module.
- 1.11 Adjust the input voltage to 3.000 Volts.
- 1.12 While monitoring the READOUT output BNC with the frequency counter, adjust the front panel trimpot labeled O1 for a display of 31.200 kHz on the counter. this supplies an offset voltage of 1.200 volts to the input of VFCl.
- 2.0 Alignment Procedure - ZERO
- 2.1 Set the TEST/NORMAL switch to TEST.
- 2.2 Attach a cable from the rear panel BNC labeled OFFSET VOLTAGE on the automation module to the rear panel BNC labeled OFFSET on the digital integrator module.
- 2.3 Apply an input voltage of 10.00 mVolts to the rear panel BNC labeled CURRENT on the Digital Integrator module.
- 2.4 Monitor the front panel output BNC labeled ZERO with the frequency counter.
- 2.5 Adjust the front panel trimpot labeled Z2 for a reading of 100 Hz on the frequency counter.
- 2.6 Adjust the input voltage supply for an output of 3.000 volts.
- 2.7 Adjust the front panel trimpot labeled FS2 for a reading of 30.000 kHz on the frequency counter.
- 2.8 Repeat steps 2.3 to 2.7 to determine repeatability and accuracy of adjustments made previously. Repeatability should be  $\pm 1$  count.
- 2.9 Set the TEST/NORMAL switch to NORMAL.

- 2.10 Continue monitoring the front panel output BNC labeled ZERO with the frequency counter.
- 2.11 Adjust the front panel trimpot labeled 02 for a reading of 1200 Hz on the frequency counter.
- 3.0 Alignment Procedure - Oscillator
- 3.1 Monitor the rear panel BNC labeled CLOCK on the digital integrator with the frequency counter.
- 3.2 Adjust the front panel trimpot labeled CLOCK for a frequency of 10.000 kHz.
- 4.0 Calibration Procedure - Internal Protection Circuit

The internal protection circuit senses the CURRENT input and isolates this input from the VFC's if the input is driven more than 300 mVolts negative.
- 4.1 Remove any external connections from the rear panel CURRENT BNC.
- 4.2 Measure the voltage at pin 3 of U2 with a DVM and set the front panel trimpot labeled NCO for a reading of -300 mV on the DVM.
- 4.3 Short pin 2 and pin 6 of U2 together and monitor the output voltage on pin 6 of U2.
- 4.4 Adjust the front panel trimpot labeled A2 for a reading of -300 mV on the DVM.
- 4.5 Remove the jumper between pins 2 and 6 on U2.
- 4.6 Jumper pins 2 and 3 of U1 to ground.
- 4.7 Monitor pin 6 of U1 with a DVM and adjust the front panel trimpot A1 clockwise until the DVM reads +15 volts. Next, adjust the trimpot slowly counter-clockwise until the DVM reads a negative voltage.
- 4.8 Remove all jumpers from the circuit. This completes the protection circuit adjustments. The following steps outline testing of the calibration.
- 4.9 Set the TEST/NORMAL switch to NORMAL.

- 4.10 Connect the Fluke Voltage Supply to the Current Input on the Digital Integrator module. Input 0.00 Volts.
- 4.11 Remove the shorting bar connecting the negative terminal on the 343A to the ground terminal.
- 4.12 Monitor pin 2 and 4 on U3 with the DVM. The voltage should be approximately 4 Volts (logic "1"). Slowly decrease the voltage from 0.00 Volts to -500 mVolts. At approximately -300 mVolts pin 2 and 4 on U3 should read less than 1 Volt (logic "0") causing K1 to energize, isolating the CURRENT input from the VFC's.

NOTE: To obtain the negative voltage from the 343A calibrator, reverse the leads.

## COUNTER MODULE

### Description

The counter module is a multiplexed seven digit, three counter scaler with a rear panel BCD output. There are no calibrations or internal adjustments.

The three counter inputs are clock, readout, and zero. They accept the frequency outputs of the digital integrator module. All three counter circuits are identical. Their outputs can be selected through a front panel switch or through the rear NIM connector. Details of how each counter is selected is presented in the automation module description.

The readout and zero signal are described in the digital integrator description and are not discussed here. The clock signal is a precise 10 kHz signal generated by a crystal oscillator in the digital integrator. This signal provides the timing signal used to multiplex the eight digit LED display (Note: the eighth digit is not used). In addition, the clock signal provides accurate timing for the computer to use when timing oxidation reactions.

### 1.0 Controls

- 1.1 MANUAL/AUTO - Local operation of the counter module is provided when this switch is in the MANUAL position. It enables all front panel switches. When the switch is in the AUTOMATIC position, all controlling is done by the computer through the automation module and the front panel switches are disabled except the LAMP TEST button.



- 1.2 CLEAR - When the MANUAL/AUTO switch is in the MANUAL position, the pushbutton clears all three circuits to zero.
- 1.3 LAMP TEST - When pressed, all the digits in the LED display are illuminated.
- 1.4 ZERO/READOUT/CLOCK - When the MANUAL/AUTO switch is in the MANUAL position, this switch selects any of the three internal counters for display on the front panel LED display.
- 1.5 HOLD - When the MANUAL/AUTO switch is in the MANUAL position, this pushbutton latches the display; however, the internal counters continue to increment. If the CLEAR pushbutton is pressed while the display is latched, the display will not clear but the counters will. Once the hold is removed, the display will represent the value in the counters.

## **AUTOMATION MODULE**

### **Description**

The automation module contains two circuit boards, a +/- 15 VDC supply, a 1.35 V mercury battery, and 120 VAC power switch with a neon lamp. The automation module's main function is to act as an interface between the computer and the rest of the system. The computer sends signals to the GPIO interface, which directly connects to the driver/interface card (Ref print S5-2-15965). The driver/interface card amplifies the signals to a level that can drive the relays located on the automation card and on the switching and load card. The driver/interface card also drives the signals that control the counter module operation.

The list below describes computer controlled functions performed by the automation module. Each operation references to the pin on the driver/interface board which executes it. Note for these functions to be enabled, the mode switch on the potentiostats must be in the standby position and the manual/auto switch on the counter module must be in the auto position. A logic one is defined as 3 to 5 volts and a logic zero as less than .7 volts.

- A logic one signal applied to pins H and J energize K1 and K2 respectively on the switching and load card which puts the oxidation potentiostat in the operate mode.
- Pin E and F operate identical to H and J except the reduction potentiostat replaces the oxidation potentiostat.

- A logic zero signal applied to pin U energizes K1 and K2 on the automation card. When energized, the oxidation potentiostat is connected to the cell and when de-energized, it is connected to a calibration resistor.
- Pin T operates identical to Pin U except the reduction potentiostat replaces the oxidation potentiostat.
- A logic zero signal applied to pin S energizes K5 and K6 on the automation card. When energized, the external potential input is connected to external potential #1. When de-energized, the external potential input is connected to external potential #2. This feature is not presently used due to noise problems. Instead, the external potential signal is supplied to both potentiostats at the same time. This does not cause a problem since only one potentiostat is operated at a time.
- A logic zero signal applied to pins R, M, or L, will energize K7, K8, and K9, respectively. These three relays combined allow the digital volt meter to measure the voltage at various points in the system. Table 1 below details the operation. Note "X" indicates a don't care state.

**TABLE 1**

**Logic Table for Controlling DVM Monitor Point**

<u>K7</u>	<u>K8</u>	<u>K9</u>	<u>DVM Monitors</u>
X	De-energized	De-energized	Calibration resistor
De-energized	X	Energized	Load #1
X	Energized	De-energized	Ox/red potential
Energized	X	Energized	Load #2

- Pin D and C control "select A" and "select B", respectively. Their outputs do not drive relays but directly connect to the counter/multiplexer board in the counter module. When the manual/auto switch on the counter module is in the auto position, Pins D and C choose which counter circuit, readout, zero or clock is to be displayed. Table 2 shows how select A and select B operate.

**TABLE 2****Logic Table for Controlling Counter Module Data Multiplexer**

<u>Select A</u>	<u>Select B</u>	<u>Counter Circuit Selected</u>
Logic zero	Logic zero	READOUT
Logic zero	Logic one	CLOCK
Logic one	Logic zero	ZERO
Logic one	Logic one	Not Used

For explanation of each circuit, see the digital integrator description.

**SYSTEM SOFTWARE**

The Hewlett-Packard series 200 model 226 utilizes RAM based Basic language version 3.0. The software developed for the instrument control is a user-friendly menu driven system. Two main menus, each having a version for the type of supporting electrolyte, which allow access to the various system capabilities are reproduced in Figures 4-7. Password control of selected system capabilities prevent unintentional access to these capabilities by the technician performing routine measurements.

The software package for plutonium measurement permits the determination of pure plutonium in both a sulfuric acid supporting electrolyte at a platinum working electrode and a nitric acid supporting electrolyte at a gold working electrode. In both method variations, measurement of the supporting electrolyte blank precedes sample determination. The formal potential of the pure plutonium can be measured for either supporting electrolyte following sample measurement. Concurrent measurement of plutonium and small amounts of neptunium or iron [8] in nitric acid supporting electrolyte has been fully automated. Electrical calibrations are performed with each plutonium measurement and continuously when the system is not being used for plutonium measurements. To comply with SRP measurement and test equipment (MTE) requirements, the system automatically monitors parameters of key components to ensure that no significant shifts have occurred which would adversely affect the reliability of the system calibration.

All measurement results and related information are logged to a data disc in serially numbered data files. The same data disc also contains the results of the last ten electrical calibrations to be used in subsequent measurement calculations, and the data from the last blank measurement.

```

+-----+
|           MENU 1           |
| 1 = MEASURE SULFURIC BLANK |
| 2 = CONDITION/MEASURE BLANK|
| 3 = MEASURE CALIBRATION FACTOR |
| 4 = CHANGE TYPE OF ELECTROLYTE |
| 5 = CHANGE CALIBRATION FACTOR |
| 6 = CHANGE CALIBRATION TIMING |
| 7 = CHANGE MEASUREMENT TIMING |
| 8 = MANUAL OPERATION/TESTING |
| 9 = SKIP TO MENU_2         |
| 10 = MINI-PROCEDURE        |
| 11 = MTE LINEARITY TEST    |
| 99 = STOP                  |
| ENTER # OF THE DESIRED OPTION |
+-----+

```

FIGURE 4. Menu 1 for Sulfuric Acid Supporting Electrolyte

```

+-----+
|           MENU 2           |
| 10 = MEASURE NEXT BLANK    |
| 20 = MEASURE THE SAMPLE    |
| 30 = RE-MEASURE THE SAMPLE |
| 60 = MEASURE FORMAL POTENTIAL |
|       OF Pu(IV)/Pu(III) COUPLE |
| 70 = PRINT/DISPLAY LOGBOOK |
| 80 = PRINT PROGRAM LIST     |
| 99 = LAST MEASUREMENT DONE  |
| ENTER # OF THE DESIRED OPTION |
+-----+

```

FIGURE 5. Menu 2 for Sulfuric Acid Supporting Electrolyte

```

+-----+
|           MENU 1           |
| 1 = MEASURE THE NITRIC BLANK |
| 3 = MEASURE CALIBRATION FACTOR |
| 4 = CHANGE TYPE OF ELECTROLYTE |
| 5 = CHANGE CALIBRATION FACTOR |
| 6 = CHANGE CALIBRATION TIMING |
| 7 = CHANGE MEASUREMENT TIMING |
| 8 = MANUAL OPERATION/TESTING |
| 9 = SKIP TO MENU 2          |
|10 = MINI-PROCEDURE          |
|11 = MTE LINEARITY TEST      |
|99 = STOP                    |
| ENTER # OF THE DESIRED OPTION |
+-----+

```

FIGURE 6. Menu 1 for Nitric Acid Supporting Electrolyte

```

+-----+
|           MENU 2           |
|10 = MEASURE NEXT BLANK      |
|20 = MEASURE THE SAMPLE      |
|30 = RE-MEASURE THE SAMPLE   |
|40 = CONCURRENT ELECTROLYSIS |
|   OF Pu/Np SAMPLE           |
|50 = CONCURRENT ELECTROLYSIS |
|   OF Pu/Fe SAMPLE           |
|60 = MEASURE FORMAL POTENTIAL |
|   OF Pu(IV)/Pu(III) COUPLE  |
|70 = PRINT/DISPLAY LOGBOOK   |
|80 = PRINT PROGRAM LIST      |
|99 = LAST MEASUREMENT DONE   |
| ENTER # OF THE DESIRED OPTION |
+-----+

```

FIGURE 7. Menu 2 for Nitric Acid Supporting Electrolyte

The program listings for plutonium and neptunium measurement with cross reference listings (XREF) are available from the National Energy Software Center.<sup>9</sup>

## OPERATION

### Startup/Shutdown

The system startup and shutdown procedures are outlined below. The startup procedure utilizes the Hewlett-Packard AUTOSTART capability. The autostart program written for this system ensures loading of all of the desired features of Basic 3.0 Language, the loading of the plutonium measurement software, and returning of the log disc to the disc driver. The autostart program is included in the appendixes.

#### 1.0 Instrument Startup Procedure

##### 1.1 Turn on the following items in order listed below:

- Power strips in rear of cabinet
- HP2673A Graphics Printer
- HP8116A Pulse/Function Generator
- HP5316A Universal Counter
- HP3456A Digital Voltmeter
- NIM Bin POWER switch
- Automation Module-AC POWER switch
- 5V supply on bottom panel

##### 1.2 Press "DISABLE" button on HP8116A so as to turn off the RED light.

##### 1.3 Ensure that the Coulometer Counter/Timer switch is in AUTO position.

##### 1.4 Ensure that the DVM terminal button is in the REAR position.

##### 1.5 With computer power off, insert disk labeled "AUTOSTART" and press AC LINE switch to turn power on. Wait for "READ" light on disc drive to turn off.

##### 1.6 Follow instructions on the CRT.

##### 1.7 Insert LOGBOOK disc.

1.8 Press RUN key and the program should begin. If program does not display a safety message, then hold down SHIFT key while pressing PAUSE key (RESET), then press RUN key.

1.9 Answer the safety and time and date questions.

NOTE: The ENTER, CONTINUE, OR EXECUTE keys may be used interchangeably to input answers to questions but the CONTINUE key is required when requested.

## 2.0 Instrument Shutdown Procedures

Coulometer shutdown procedure, used for emergency generator load test and instrument maintenance.

2.1 Remove disc from drive

2.2 Reset computer (SHIFT/PAUSE)

2.3 Turn the following to OFF position in order:

- 5V supply on bottom panel
- Automation module - AC POWER
- NIM Bin POWER
- Voltmeter
- Counter
- Pulse/Function Generator
- Printer

2.4 Turn off the computer - AC line button

2.5 Turn off power strips in rear of cabinet

## 3.0 Loading the Neptunium Measurement Software

3.1 Remove the log disc from the internal disc driver

3.2 Insert the neptunium program disc into the driver

3.3 Type 'SCRATCH' then press the 'EXECUTE' key

3.4 Type 'SCRATCH A' then press the 'EXECUTE' key

3.5 Type 'LOAD "NP\_CPC"' then press the 'EXECUTE' key

- 3.6 Allow the program to load, then remove the neptunium program disc
  - 3.7 Return the Log disc to the driver
  - 3.8 Press the 'RUN' key
  - 4.0 Initiating Sample Measurements -- Daily: The startup procedure to go from continuous, electrical calibration to sample measurement, is outlined in the mini-procedure reproduced below from the software.
  - 4.1 Inspect CAL RESULTS/PRINTOUT
    - a. AVE CAL Factor Near 2.000
    - b. RSD less than .0010%
  - 4.2 Press 'RESET' (SHIFT/PAUSE)
  - 4.3 Press 'RUN'
  - 4.4 Answer questions
  - 4.5 When the electrode starts to reduce the blank, switch the two toggle switches on Pot #2 to 'NORMAL' & '20mA' then disconnect counter electrode. The meter reading should decrease to <.004 V. Reconnect counter electrode. Return switch on Pot #2 to 'CHARGE' and '200mA'.
  - 4.6 If blank run is excessively slow, then repeat steps 2, 3, & 4.
- NOTE: When Menu 2 is used to start sample: If '10' is entered mistakenly, repeat steps 2-4 then select Option 9 when Menu 1 appears. Option 9 reloads last blank data and returns to Menu 2. If system fails any of these tests, do not proceed. Contact M. K. Holland or your supervisor.

#### 5.0 Software Menu Usage

The software menu #1 is obtained after a reminder about safety and quality, a check of the system clock setting, and the type of electrolyte to be used is specified. The type of electrolyte must be specified first because both menus have two versions depending upon the type of electrolyte.



For routine measurements, the supporting electrolyte blank is measured before each sample is determined. Following blank measurements, the sample is dissolved in the electrolyte upon which the blank was determined. After sample measurement is completed, the next blank is routinely determined; however, the sample may be remeasured to verify the result or the formal potential may be determined.

The formal potential and the complete MTE test are performed periodically as dictated by SRP quality assurance requirements and when considered appropriate by the technical staff.

Other menu accessible capabilities are sufficiently self explanatory.

6.0 Completion of sample measurements -- daily: When measurements are completed, the system is advised from the Menu #2 Option 99. The system reminds the analyst of the items listed below, performs a short MTE linearity test, and initiates continuous electrical calibration until the next time the system is used for measurements.

REMINDERS TO THE ANALYST:

- Turn off the argon gas
- Flame the electrodes
- Enter assay results log book
- Enter assay results in LASSO if appropriate
- Do not turn 9826 computer off
- Hit 'CONTINUE' when ready

**RESULTS AND DISCUSSION**

Significant improvements in the system design were achieved. Noise in the control potential circuit caused by the automation module was identified and corrected. Further reduction of noise in the voltage signal to be integrated was accomplished by addition of a 100 nF capacitor AC shunt. Problems with the isolation of the manual potentiostat mode switch from the automation module, the isolation of the charging system for the potentiostat current meter, stability of the control amplifier circuit, and the operation of the polarity protection circuit for the digital integrator VFC circuit were also identified and corrected. The peripheral units were upgraded by the selection of more recently introduced Hewlett-Packard devices.

The system capabilities have been improved by upgrading the system control software. The capability to measure plutonium in dilute sulfuric acid at a platinum electrode has been added to the previous plutonium method variation in dilute nitric acid at a gold electrode. The capacity of the system to electrolyze large samples (>35 mg Pu or Np) by monitoring the integrator for overflow was added. The control potential stability at high current levels was improved by software adjustments of the programmable digital-to-analog converter to compensate for loading of the driver amplifier in the potentiostat. The software was later modified to measure neptunium. Quality assurance tests were added to ensure the stability of components critical to accurate system calibrations.

The control potential of the potentiostat is controlled by adjusting the voltage on the non-inverting input of AR2 (Ref print S5-2-15956). In the original design, the signal from the digital-to-analog converter was connected to the automation module at the external potential input. The computer selected which potentiostat would receive the signal. The problem with routing the signal through the relays in the automation module was noticed during linearity checks of the calibration factor. The value of the calibration factor varied with the control potential. The cause of the non-linearity was traced to the noise level at the non-inverting input of AR2. The problem was corrected by connecting the digital-to-analog converter signal directly to the external potential inputs on both potentiostats. This did not cause a problem since only one potentiostat is operated at a time.

The signal generated by the potentiostat for integration is an exponentially decaying DC voltage. Accompanying the DC signal of interest is AC noise (see Figures 8-10) contributed from many sources (e.g., measurement cell stirring system, cell capacitance, potentiostat control circuit oscillation, pickup from 120V 60Hz supply, and from fluorescent lighting). The offset voltage produced by the digital integrator handles peak AC noise up to 0.1V and is most effective on low frequency noise. Several approaches have been used to ensure that AC noise is reduced. At the New Brunswick laboratory, the major source of AC noise on the signal to be integrated was from the AC noise on the power supply to the DC stirrer system of the measurement cell. Power supply noise was eliminated by a capacitor shunt across the power supply outputs. The other sources of noise in the NBL system were ignored since their combined contribution was less than 0.1V peak. Princeton Applied Research sells a noise filter accessory for the PAR cell assembly which utilizes a 100 nF capacitive shunt between the counter and reference electrode. The 100 nF capacitive shunt feature was incorporated into the potentiostat module between the counter and reference electrode output to the cell assembly which reduces AC noise (above 1 kHz) beyond that obtained by the 1000 nF




FIGURE 8. Electrolysis Current Signal




FIGURE 9. Filtered Signal to be Integrated




FIGURE 10. Signal Integrated by State of the Art Coulometer

shunt located on the cell assembly. Examination with oscilloscope showed  $<0.03$  V peak 60 cycle noise. It should be noted that an oscilloscope must be used for this measurement, not an AC voltmeter which measures RMS voltage and is not sensitive to high frequency noise which will adversely affect integration performance.

The mode switch on the potentiostat module allows for manual operation. The switch has three positions: standby, hold, and operate. These modes of operation are discussed in detail in the potentiostat module description. If the system is to be operated by the computer, the mode switch must be in the standby position. In this position, the relays can be driven by 74128 line drivers located in the automation module. In the original design, the outputs of the drivers were not isolated from the mode switch. Therefore, if the mode switch was switched to the hold or operate position, 5 volts would be applied to the outputs of the drivers (Ref print S5-2-15955). This caused the output stage of the driver to sink more than its rated current. The problem was corrected by isolating the drivers from the mode switch when in the hold or operate position.

The current meter in the potentiostat module monitors the cell current by measuring the voltage across a 1 ohm resistor in the 20 mA position and a 10 ohm resistor in the 200 mA position. Since the control amplifier board floats with respect to chassis ground, the system 5 volt power supply could not be used to power the current meter. Therefore, in the initial design, four nicad rechargeable batteries were used to supply power to the meter. However, if the display/charge switch was left in the charge position when the system power was off, the batteries would discharge. The nicad batteries were replaced by four Duracell AA batteries and isolated from the power supply. The current meter has a LCD display which draws little current and is only used when manually operating the system; therefore, the batteries will last for about two years under normal use.

The current amplifier on the control amplifier board (Ref print S5-2-15956) is a Teledyne 1461. The 1461 requires a compensation capacitor between pins 11 and 12. Initial design called for a 39 pF compensation capacitor. A problem was encountered with the amplifier circuit when attempting to zero the offset voltage. A 3 mV peak to peak ripple was present across the inputs of the amplifier that could not be nulled out with the trimpot. The value of the compensation capacitor was changed to 100 pF and the ripple was eliminated.

The polarity protection circuit senses the input voltage to the digital integrator. If the input is less than -300 mV, the circuit will isolate the input voltage from voltage-to-frequency

converter #1 (VFC #1). The initial design of the circuit did not work due to the placement of a feedback diode in the amplifier stage.

More recently introduced Hewlett-Packard peripheral devices were incorporated in the SRP systems. Precision timing in the software system is accomplished with an HP8116A Pulse/Function Generator and an HPS318A Universal Counter. The counter is also used during digital integration alignment. All voltage measurements are made with an HP3456A digital voltmeter. An HP2673A graphics printer was selected for hard copy output. No change from the NBL prototype was made when the HP59501B DAC/power supply was selected. Appropriate changes in software were made to accommodate these devices.

The modification to the software to allow for a second method variation involving a change in the type of supporting-electrolyte and working-electrode was accomplished by the addition of a variable flag (Electro-flag) set for the method variation selected. The integer variable So4 and No3 were assigned the numerical value of 1 and 2, respectively, at program startup. When the  $1N$   $H_2SO_4$  or  $N$   $HNO_3$  electrolyte selection is made, the variable Electro-flag is set to So4 or No3, as appropriate. When a program line or series of programs are dependent upon the method variation being applied, an IF.....THEN or IF.....THEN.....ELSE construction with Electro-flag=So4 or Electro-flag=No3 in the argument is used to direct program execution.

The counters in the integration module are capable of accumulating up to 9,999,999 counts, at which point the most significant bit is lost and the counter continues from zero. Ten million counts correspond to the integration of 20.0 Coulombs of electricity. When the counts due to the offset are considered, the integrator has the capacity to electrolyze up to 35 mg of plutonium or neptunium. The software has been modified to monitor the counters when the integrator is operating. A variable flag is set for each counter when it exceeds five million counts. Should the counter accumulate ten million counts and reset to zero, the software will increment a second variable flag which counts the number of times the counter overflows and resets the five million count flag for the appropriate counter.

At the beginning of the sample reduction and oxidations, the electrolysis current is sufficient to load down the driver amplifier in the potentiostat. The software monitors the control potential supplied to the cell assembly and reduces or increases the signal from the digital-to-analog-converter power supply to maintain the control potential at the desired voltage.

The system/software was later applied to the measurement of pure neptunium solutions. The number of differences between plutonium and neptunium measurement were sufficiently numerous and complex to make the use of variable flag(s) for the element being analyzed to be undesirable. A copy of the plutonium software was obtained from the plutonium controlled-potential coulometry program file (Pu\_CPC), modified to measure neptunium, and then stored in a different program file (Np\_CPC). Appropriate program upgrades are then made to both programs. The neptunium program is always removed from RAM when measurements are completed to avoid inadvertent usage.

The continuous electrical calibration is dependent upon the accurate measurement/alignment and the stability of the calibration resistor, the quartz crystal oscillator, and the digital voltmeter, and upon the stability of the load resistor and the voltage-to-frequency converters. A series of quality assurance tests were added to the plutonium measurement software to ensure the stability of each of these critical components once they have been aligned/measured.

The routine electrical calibration precision (2-sigma) is less than 0.001% at the current levels of major importance and only slightly poorer at a lower level. By using electrical calibration, shifts in components can be detected before they impact upon the measurement reliability of the system. The periodic use of the quality assurance tests ensures that the system is accurately calibrated.

Four test subroutines have been developed. The first is a load vs. calibration resistors check, performed over the range of current levels important to sample measurement. This test ensures the relative stability of the resistors and the performance of the voltmeter. The second test is a measurement of the response of voltage-to-frequency converter #1 over the range of current levels important to sample measurement. The third and fourth tests utilize electrical calibrations to evaluate the linearity of the integrator response (count/sec per volt) over the ranges of interest. Calibrations are performed at selected input values encountered in sample measurement both above and below the input value used for electrical calibration. The results of tests three and four are expressed as the percent error caused by the integrator nonlinearity on the measurement of the exponentially decaying input signal obtained during sample measurement at high and low current levels, respectively. The algebraic sum of the two nonlinearity test results quantifies the total effect of integrator nonlinearity upon sample measurement. A summary of the data/results for the four tests on both measurement systems are shown in Table 3.



**TABLE 3**

**Quality Assurance Test Results**

Test #1: Load Resistance

System 1 / Potentiostat 1	<u>Current Level, mA</u>	<u>Ohms</u>
	0.5	50.014
	3.0	49.997
	5.5	49.998
	8.0	49.998
	10.5	49.998
	13.0	49.998
	15.5	49.998
	18.0	49.998

System 2 / Potentiostat 1	<u>Current Level, mA</u>	<u>Ohms</u>
	0.5	50.040
	3.0	50.064
	5.5	50.066
	8.0	50.066
	10.5	50.067
	13.0	50.066
	15.5	50.067
	18.0	50.067

**TABLE 3, contd**

**Test #2: Integrator Response**

<b>System 1</b>	<u>Current Level, mA</u>	<u>Response - Counts/Sec/Volt</u>
	0.5	10000.50
	3.0	10000.08
	5.5	10000.19
	8.0	10000.27
	10.5	10000.21
	13.0	10000.17
	15.5	10000.14
	18.0	10000.10

<b>System 2</b>	<u>Current Level, mA</u>	<u>Response - Counts/Sec/Volt</u>
	0.5	10005.09
	3.0	10001.58
	5.5	10001.33
	8.0	10001.24
	10.5	10000.98
	13.0	10001.02
	15.5	10001.28
	18.0	10000.84

**TABLE 3, contd**

Theoretical Electrical Calibration -- Calculated from Test 1 and 2  
Results Using Equation 3.

System 1	<u>Current</u>	<u>C (Theoretical)</u>	<u>C(Measured)</u>	<u>RD %</u>
	0.5	1.99932	1.99964	+0.06
	3.0	2.00009	2.00004	+0.003
	5.5	2.00001	2.00004	-0.002
	8.0	2.00003	2.00002	+0.001
	10.5	2.00001	2.00002	-0.001
	13.0	2.00005	2.00003	+0.001
	15.5	2.00006	2.00004	+0.001
	18.0	2.00005	2.00005	0.000
	20.0	-	2.00006	-
	30.0	-	2.00005	-

System 2	<u>Current</u>	<u>C (Theoretical)</u>	<u>C(Measured)</u>	<u>RD %</u>
	0.5	1.99738	1.99710	+0.014
	3.0	1.99713	1.99708	+0.003
	5.5	1.99710	1.99710	0.000
	8.0	1.99712	1.99718	-0.003
	10.5	1.99713	1.99720	-0.004
	13.0	1.99716	1.99722	-0.003
	15.5	1.99707	1.99723	-0.008
	18.0	1.99716	1.99725	-0.005
	20.0	-	1.99726	-
	30.0	-	1.99729	-

"-" No Data Point

TABLE 3, contd

Test #3: Integrator Non Linearity at High Current Levels

System 1

<u>Current Level, mA</u>	<u>Calibration Coulombs/Million Counts</u>	<u>Precision of 10 Replicates (2-Sigma) %</u>
10	2.00002	0.002
20	2.00006	0.001
30	2.00005	0.001

Non Linearity -0.0012%

System 2

10	1.99720	0.001
20	1.99726	0.001
30	1.99729	0.000

Non Linearity -0.0033%

Test #4: Integrator Non Linearity at Low Current Levels

System 1

<u>Current Level</u>	<u>Coulombs/Million Counts</u>	<u>10 Replicates (2-Sigma) %</u>
0.1	1.99779	0.028
0.5	1.99964	0.015
1.0	1.99980	0.004
3.0	2.00004	0.003
5.0	2.00004	0.002
10.0	2.00002	0.002

Non Linearity +0.0004

System 2

0.1	1.99695	0.047
0.5	1.99721	0.017
1.0	1.99717	0.009
3.0	1.99719	0.002
5.0	1.99721	0.001
10.0	1.99720	0.001

Non Linearity -0.0001

The Combined effect of non linearity upon electrical calibration:

System #1 -0.0008%  
System #2 -0.0034%

Excellent system reliability was demonstrated using plutonium standard reference material SRM949f, traceable secondary standards, and Plutonium Metal Exchange Samples. The data are shown in Tables 4 through 6.

TABLE 4

**Coulometric Measurements of Plutonium Standards using Nitric Acid Supporting Electrolyte and a Gold Working Electrode**

<u>Primary Standards</u>	<u>Preparation 1</u>	<u>Preparation 2</u>	<u>Preparation 3</u>	<u>Preparation 4</u>
Recovery	99.96	99.98	100.02	100.02
% vs. Reference Values	100.04	99.98	100.01	100.11
			100.03	100.02
			100.00	100.08
			99.95	99.86
			99.99	-
Mean	<u>100.00</u>	<u>99.98</u>	<u>100.00</u>	<u>100.02</u>
Precision (2-Sigma)%	0.11	0.00	0.006	0.19

<u>Secondary Standards Preparation</u>	<u>n</u>	<u>Precision (2-Sigma)%</u>	<u>Mean Recovery, %</u>
A	5	0.09	100.05
B	5	0.16	99.98
C	3	0.19	99.95

TABLE 5

**Coulometric Measurements of Plutonium Standards Using Sulfuric Acid  
Supporting Electrolyte and a Plutonium Working Electrode**

Primary Standard	<u>Preparation 3</u>	<u>Preparation 4</u>
	99.98	100.17

**Secondary Standards**

<u>Preparation</u>	<u>n</u>	<u>Precision (2-Sigma), %</u>	<u>Mean Recovery, %</u>
D	15	0.16	100.06
E	12	0.18	100.05

TABLE 6

**Plutonium Measurements - System #1. Rocky Flats Metal Exchange  
Program**

<u>Exchange</u>	<u>Measurement Period</u>	<u>Measurement Bias</u>	<u>Precision (2-Sigma)</u>	<u>Metal</u>
0387	Mar 87 - May 87	0.0%, 0.0%	0.1%, 0.2%	A84, B84
1286	Dec 86 - Feb 87	0.0%, 0.0%	0.1%, 0.1%	A84, C84
0986	Oct 86 - Nov 86	0.0%, -0.1%	0.1%, 0.2%	B84, C84
1285	Jul 86 - Sep 86	0.1%, -0.1%	0.1%, 0.2%	B84, C84
0985	May 86 - Jun 86	0.0%, 0.0%	0.1%, 0.2%	A84, B84
0685	Jan 86 - Mar 86	0.0%, -0.1%	0.1%, 0.3%	A84, C84

The ease of operation and the self-monitoring features of the system have contributed to its success at SRP. In addition to improving reliability thru the improved repeatability of automated instruction and enhanced quality assurance in data reduction, the system permits focusing of the analysts training upon the critical aspects of sample preparation.

Future plans include:

- Interfacing a programmable quartz thermometer to the system to improve the measurement of the solution temperature and to allow automatic monitoring of the instrument's environment.
- Continuing refinement of the system and related preparation reliability and efficiency.
- Interfacing of the system to the laboratory data base management system.
- Continuing application of the system to critical measurement needs in the support of other development/evaluation programs.

#### ACKNOWLEDGMENTS

The Authors wish to thank the following individuals: Thomas L. Frazzini for his major role in the design and assembly of the systems developed at the New Brunswick Laboratory (NBL); Charles E. Pietri, Jon R. Weiss and Kenneth Lewis, for their significant technical contributions to the system development and related measurement methodology at NBL; James Miranda, Argonne National Laboratory, for his suggestions on design improvements; Dennis W. Lewis, Sandra F. Nappier, and Daniel L. Stewart for the numerous plutonium measurements contributed during system testing and evaluation; Larry E. Feutral for the timely and maticulous drafting preparation of the system blue prints; and Whit Whittle, E & I, for his attention to detail in wiring the system.

## REFERENCES


1. T. L. Frazzini, M. K. Holland, J. R. Weiss, and C. E. Pietri, **Anal. Chem.** **52**, 2112 (1980).
2. T. L. Frazzini, M. K. Holland, J. R. Weiss, and C. E. Pietri, **Anal. Chim. Acta** **129**, 125 (1981).
3. T. L. Frazzini, M. K. Holland, C. E. Pietri, and J. R. Weiss, U.S. Patent #4,244,800, January 13, 1981.
4. M. K. Holland, T. L. Frazzini, J. R. Weiss, and C. E. Pietri, NBS Special Publication 582, p. 164 (1980).
5. J. R. Harrar and E. Behrin, **Anal. Chem.** **39**, 1230 (1967).
6. M. K. Holland, J. R. Weiss, and C. E. Pietri, **Anal. Chem.** **50**, 235 (1978).
7. J. R. Weiss, E. F. Groh, and D. A. Cassidy, U.S. DOE Report, NBL-289, 75 (1979).
8. M. K. Holland and K. Lewis, **Anal. Chim. Acta** **149**, 167 (1983).
9. M. K. Holland and J. V. Cordaro, **An Automated Instrument for Controlled-Potential Coulometry: System Documentation**, Programs listing available from the National Energy Software Center, Argonne National Laboratory, 9700 South Cass Avenue, Argonne, Illinois 60439.



## APPENDIX A

### DRAWING LIST

1. S5-2-15950 System No. 1 Assembly
2. S5-2-15951 System No. 2 Assembly
3. S5-2-15952 System No. 1 Schematic Assemblies and Details
4. S5-2-15953 Interconnection Diagram System 1 & 2
5. S5-2-15954 NIM Power Supply Arrangement Detail and Connector Function Table
6. S5-2-15955 Potentiostat Module Assembly, Connection Diagram and Parts Lists
7. S5-2-15956 Switching and Load and Control Amplifier Printed Circuit Boards, Schematic and Component Layout
8. S5-2-15957 Potentiostat Module Details
9. S5-2-15958 Integrator Module Assembly, Connection Diagram and Parts Lists
10. S5-2-15959 Voltage-To-Frequency Printed Circuit Board, Schematic and Component Layout
11. S5-2-15960 Integrator Module Details
12. S5-2-15961 Counter/Timer Module, Assembly Connection Diagram and Parts List
13. S5-2-15962 Display/Driver Printed Circuit Board, Schematic and Component Layout
14. S5-2-15963 Counter/Multiplexer Board, Schematic and Component Layout
15. S5-2-15964 Counter-Timer Details
16. S5-2-15965 Automation Module, Assembly Connection Diagram and Parts List
17. S5-2-15966 Automation Printed Circuit Board, Schematic and Component Layout
18. S5-2-15967 Driver/Interface Printed Circuit Board, Schematic and Component Layout
19. S5-2-15968 Automation Module Details
20. S5-2-15969 Cable Details




REPRODUCED FROM  
BEST AVAILABLE COPY

52	F3	C O M L	1	FUSE, 500 SLUG, 50 AMP, 250V, 1000 PSI
53	W3-W03	(10-2-1952)	5	3000 PSI FUSE HOLDER, DOWNSIDE, 1000 PSI
54			6	CABLE, HEAVY DUTY, 1000 PSI
55	W1-W18	(10-2-1952)	18	CABLE, 1000 PSI, DISCONNECTION DIAGRAM
56		C O M L	16	DISCONNECTION DIAGRAM
57			20	GUNMETER, 0-20 PSI, FOR 1000 PSI METAL
58		55-2-1952D	10	WATER DIVIDER NETWORK
59			16	WATER DIVIDER NETWORK
60	P61	165-2-1952A	1	1000 PSI, 1000 PSI
61			1	WATER DIVIDER, 1000 PSI, 1000 PSI
62		C O M L	6	WATER DIVIDER, 1000 PSI, 1000 PSI
63	J1		1	WATER DIVIDER, 1000 PSI, 1000 PSI
64			1	WATER DIVIDER, 1000 PSI, 1000 PSI
65			1	WATER DIVIDER, 1000 PSI, 1000 PSI
66			1	WATER DIVIDER, 1000 PSI, 1000 PSI
67			1	WATER DIVIDER, 1000 PSI, 1000 PSI
68			1	WATER DIVIDER, 1000 PSI, 1000 PSI
69			1	WATER DIVIDER, 1000 PSI, 1000 PSI
70			1	WATER DIVIDER, 1000 PSI, 1000 PSI
71			1	WATER DIVIDER, 1000 PSI, 1000 PSI
72			1	WATER DIVIDER, 1000 PSI, 1000 PSI
73			1	WATER DIVIDER, 1000 PSI, 1000 PSI
74			1	WATER DIVIDER, 1000 PSI, 1000 PSI
75			1	WATER DIVIDER, 1000 PSI, 1000 PSI
76			1	WATER DIVIDER, 1000 PSI, 1000 PSI
77			1	WATER DIVIDER, 1000 PSI, 1000 PSI
78			1	WATER DIVIDER, 1000 PSI, 1000 PSI
79			1	WATER DIVIDER, 1000 PSI, 1000 PSI
80			1	WATER DIVIDER, 1000 PSI, 1000 PSI
81			1	WATER DIVIDER, 1000 PSI, 1000 PSI
82			1	WATER DIVIDER, 1000 PSI, 1000 PSI
83			1	WATER DIVIDER, 1000 PSI, 1000 PSI
84			1	WATER DIVIDER, 1000 PSI, 1000 PSI
85			1	WATER DIVIDER, 1000 PSI, 1000 PSI
86			1	WATER DIVIDER, 1000 PSI, 1000 PSI
87			1	WATER DIVIDER, 1000 PSI, 1000 PSI
88			1	WATER DIVIDER, 1000 PSI, 1000 PSI
89			1	WATER DIVIDER, 1000 PSI, 1000 PSI
90			1	WATER DIVIDER, 1000 PSI, 1000 PSI
91			1	WATER DIVIDER, 1000 PSI, 1000 PSI
92			1	WATER DIVIDER, 1000 PSI, 1000 PSI
93			1	WATER DIVIDER, 1000 PSI, 1000 PSI
94			1	WATER DIVIDER, 1000 PSI, 1000 PSI
95			1	WATER DIVIDER, 1000 PSI, 1000 PSI
96			1	WATER DIVIDER, 1000 PSI, 1000 PSI
97			1	WATER DIVIDER, 1000 PSI, 1000 PSI
98			1	WATER DIVIDER, 1000 PSI, 1000 PSI
99			1	WATER DIVIDER, 1000 PSI, 1000 PSI
100			1	WATER DIVIDER, 1000 PSI, 1000 PSI

FORM NO.	REVISION	REVISION NO. DATE	ISSUED BY DATE	APPROVED BY DATE	<div>FINISHED - ALL DRAWINGS SURFACES U O S BREAK ALL SHARP EDGES U O S  TWO DIMENSIONS IN EACH VIEW SYMMETRICAL PARTS MAY BE DIMENSIONED BY HALF VALUE IN A &amp; A. DIMENSIONS</div> <div>TOLERANCES ALL DIMENSIONS UNLESS OTHERWISE NOTED SHALL BE  ONE HOLE ORIGINALS ± .004 TWO PLACE ORIGINALS ± .001 THREE PLACE ORIGINALS ± .0005</div>	<div>THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS &amp; CO. THE INFORMATION AND KNOW HOW THEREON MAY NOT BE USED NOR THE DRAWING REPRODUCED WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS IN WHOLE OR IN PART INCLUDING VENDOR'S SHOP DRAWINGS SHALL BEAR OR REFER TO THIS STAMP</div> <div>S5-2-15950  LATEST REVISION ON THIS DRAWING</div>		
					UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO., INC.	SAVANNAH RIVER PLANT	SEPARATIONS	
					REFERENCE DRAWINGS 55-2-15951 ASSEMBLY SYSTEM N/E 55-2-15952 CABINET & POWER & VOLTAGE DIVIDER 55-2-15953 INTERCONNECTION DIAGRAM 55-2-15954 NIM BIN WIRING DIAGRAM 55-2-15955 POTENTIOMETER MODULE 55-2-15956 INTEGRATOR MODULE 55-2-15957 COUNTER TIMER MODULE 55-2-15958 AUTOMATION MODULE 55-2-15959 CABLE DETAILS	FRACTIONAL 1 ANGULAR 1	BLDG. NO. 772-F PROBLEM NO. 2-3500	TITLE COULOMETER SYSTEM  SYSTEM N/E1
					APPROVED BY DATE	CHECKED BY DATE	DRAWN BY DATE	INCHES FACETED IN
					SPECIAL INSTRUCTIONS U.S. GOVERNMENT U.S. ARMY U.S. NAVY U.S. AIR FORCE U.S. MARINE CORPS U.S. COAST & GEOD. SURV. U.S. SPACE & AERONAUTICS ADMIN. U.S. NUCLEAR ENERGY COM. U.S. BUREAU OF RESEARCH U.S. BUREAU OF STANDARDS U.S. BUREAU OF MINES U.S. BUREAU OF GEOLOGY U.S. BUREAU OF METEOROLOGY U.S. BUREAU OF AERONAUTICS U.S. BUREAU OF MARITIME SERVICE U.S. BUREAU OF NAVIGATION U.S. BUREAU OF LENSES U.S. BUREAU OF PHOTOGRAPHY U.S. BUREAU OF CHEMISTRY U.S. BUREAU OF BIOLOGY U.S. BUREAU OF AGRICULTURE U.S. BUREAU OF COMMERCE U.S. BUREAU OF EDUCATION U.S. BUREAU OF LABOR U.S. BUREAU OF MINING U.S. BUREAU OF PUBLIC HEALTH U.S. BUREAU OF REVENUE U.S. BUREAU OF THE ARMY U.S. BUREAU OF THE NAVY U.S. BUREAU OF THE AIR FORCE U.S. BUREAU OF THE MARINE CORPS			
					ASSEMBLY			
					REFERENCE STANDARDS			
					DATE 7-1-67	SCALE 3"=1'-0"	FILE 20/60	S5-2-15950

Analyze this print for  
SAFETY CONSIDERATIONS



Drawing 2



10	COM L	1	DVM	HEWLETT-PACKARD # 3450A
10		1	QUARTZ THERMISTOR	HEWLETT-PACKARD # 2804A
17		1	COUNTER	HEWLETT-PACKARD # 5518A
10		1	OSCILLOSCOPE	HEWLETT-PACKARD # 5501D
10		1	ALPHA GENERATOR	HEWLETT-PACKARD # 518A
14	55-2-15930	1	INTERFACER MODULE	(SEE DETAIL)
15	55-2-15901	1	COUNTER-TIMER MODULE	(SEE DETAIL)
15	55-2-15900	1	POTENTIOMETER MODULE	(SEE DETAIL)
11	55-2-15900	1	AUTOMATION MODULE	(SEE DETAIL)
10	COM L	1	NIM BIN	TECHNICAL #78
9		1	PLATE, RACK, 5/8" X 1/2" X 1/8"	
8		1	RAD. RADIO, #PA-100-LO	
8		1	SWEEP AUTO	CHALLENGER, #80-0001-08
7		1	HAIR-PEEP	SCANNING
7		1	LOCAL BINDER	#80-01-00
4		2	RACK MOUNTED #48-1230-08	
3		1	EQUINO ELECTRONICS	4000
3		1	COMPUTER	HEWLETT-PACKARD #55205
4		1	READ/PRINT	HEWLETT-PACKARD #2073A
5	55-2-15900	5	CHARACTER GENERATOR	HEWLETT-PACKARD #2073A
2	55-2-15900	10	CHARACTER GENERATOR	HEWLETT-PACKARD #2073A
1	55-2-15900	1	POWER SUPPLY	HEWLETT-PACKARD #2073A
1	55-2-15900	1	POWER SUPPLY	HEWLETT-PACKARD #2073A

THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOWLEDGE THEREON MAY NOT BE USED FOR THE DESIGNING OF ANY PRODUCT WITHOUT THE WRITTEN PERMISSION OF E. I. DU PONT DE NEMOURS & CO. IN PART OR IN WHOLE. REFER TO THIS STAMP.		<b>55-2-15951</b> LATEST REVISION ON THIS DRAWING	
UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO., INC. SAVANNAH RIVER PLANT			
DRAWING NO. 7725 TITLE COULOMETER SYSTEM N82 SCALE 7:2-07 DATE 9-1-10		REVISIONS 1. 55-2-15951 2. 55-2-15951 3. 55-2-15951 4. 55-2-15951 5. 55-2-15951 6. 55-2-15951 7. 55-2-15951 8. 55-2-15951 9. 55-2-15951 10. 55-2-15951 11. 55-2-15951 12. 55-2-15951 13. 55-2-15951 14. 55-2-15951 15. 55-2-15951 16. 55-2-15951 17. 55-2-15951 18. 55-2-15951 19. 55-2-15951 20. 55-2-15951 21. 55-2-15951 22. 55-2-15951 23. 55-2-15951 24. 55-2-15951 25. 55-2-15951 26. 55-2-15951 27. 55-2-15951 28. 55-2-15951 29. 55-2-15951 30. 55-2-15951 31. 55-2-15951 32. 55-2-15951 33. 55-2-15951 34. 55-2-15951 35. 55-2-15951 36. 55-2-15951 37. 55-2-15951 38. 55-2-15951 39. 55-2-15951 40. 55-2-15951 41. 55-2-15951 42. 55-2-15951 43. 55-2-15951 44. 55-2-15951 45. 55-2-15951 46. 55-2-15951 47. 55-2-15951 48. 55-2-15951 49. 55-2-15951 50. 55-2-15951 51. 55-2-15951 52. 55-2-15951 53. 55-2-15951 54. 55-2-15951 55. 55-2-15951 56. 55-2-15951 57. 55-2-15951 58. 55-2-15951 59. 55-2-15951 60. 55-2-15951 61. 55-2-15951 62. 55-2-15951 63. 55-2-15951 64. 55-2-15951 65. 55-2-15951 66. 55-2-15951 67. 55-2-15951 68. 55-2-15951 69. 55-2-15951 70. 55-2-15951 71. 55-2-15951 72. 55-2-15951 73. 55-2-15951 74. 55-2-15951 75. 55-2-15951 76. 55-2-15951 77. 55-2-15951 78. 55-2-15951 79. 55-2-15951 80. 55-2-15951 81. 55-2-15951 82. 55-2-15951 83. 55-2-15951 84. 55-2-15951 85. 55-2-15951 86. 55-2-15951 87. 55-2-15951 88. 55-2-15951 89. 55-2-15951 90. 55-2-15951 91. 55-2-15951 92. 55-2-15951 93. 55-2-15951 94. 55-2-15951 95. 55-2-15951 96. 55-2-15951 97. 55-2-15951 98. 55-2-15951 99. 55-2-15951 100. 55-2-15951	

Technical drawing of a square plate. The overall dimensions are 10" by 10". The plate has four corner fasteners. A note points to the center of the plate: "40° DRILL 4 HOLES". The distance from the center to the edge is 4". The distance from the center to the corner fastener is 4". The distance from the center to the edge is 4".

Technical drawing of a rectangular plate with the following specifications:



- Overall width: 20"
- Overall height: 31.5"
- Left edge thickness: 1/2"
- Top edge thickness: 1/2"
- Bottom edge thickness: 1/2"
- Right edge thickness: 1/2"
- Top-left corner hole: 1 3/4" DIA HOLE
- Top-right corner hole: 1/8" DRILL & HOLES
- Bottom-left corner hole: 1/8" DRILL & HOLES
- Bottom-right corner hole: 1/8" DRILL & HOLES
- Internal hole pattern: 4 holes arranged in a 2x2 grid, with 5" spacing between them.


$$SCALE = 10^1 = 1.0^1$$
[illegible]


CONTROLLED-POTENTIAL COULOMETER

1/2" HIGH LETTERS  
ENGRAVE AS SHOWN

1/8" THICK DAKLITE  
GLOSSY, BLACK-WHITE-BLACK




ITEM (2)



40° DRILL, 2 HOLES

7	COM L	1	LABEL, DANDER PANOUT CODE #P90-1B
7	COM L	1	#B-SOUNCA 1/2" LONG DINDER
6	COM L	4	SCREEN
5	COM L	2	#B-SOUNCA 1/4" LONG DINDER
4	COM L	2	SCREEN
3	COM L	2	SPACER, INSULATED, 1/2" LONG
2	COM L	2	CLER FOR #B-SOUNCA SCREEN
1	COM L	1	QUARD (SEE DETAIL)
		2	BRACKET
		2	SPR DETAIL
		1	POWER SUPPLY
		1	POWER PLANT REC-B ON-E
ITEM	PART NUMBER	QTY	DESCRIPTION
ASSEMBLY	NAME	PART AND QUANTITY	SCALE
POWER	SUPPLY	95-2-10992-2	3/4"=0"



6		COM L	4	PA-ROUNDS: 25000 BOUNDS HEAD BACON: 25000 HEAD NEX NUT
5	M			61000 BAC, 50000 5AMP IN HOLD LITTERING: 50000
4	SV			PUTIN: 10000, 50000, 50000 CUTTER: 10000, 50000, 50000
3	NI			ANALOGIC COPY: 50000, 50000, 50000
2				ANALOGIC COPY: 50000, 50000, 50000
1	COM L		1	ANALOGIC COPY: 50000, 50000, 50000
ITEM	PART DESCR	QTY DESCR	UNIT DESCR	DESCRIPTION
ASSEMBLY NAME		PART AND ORIGIN NR		
5000 POWER SUPPLY		55-2-1500-2 36100		

THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOW HOW THEREON MAY NOT BE USED NOR THE DRAWING REPRODUCED WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS IN WHOLE OR IN PART, INCLUDING VENDOR'S SHOP DRAWINGS, SHALL BEAR OR REFER TO THIS STAMP.	S5-2-15952 LATEST REVISION ON THIS DRAWING
--	--

SAVANNAH RIVER PLANT	PAGE 66	SEPARATIONS
----------------------	---------	-------------

BLDG. NO.	PROBLEM NO.	EXP. CODE	PROJECT NO.
772-F	2-3368		

TITLE	DESIGNED BY J.V. GORDARO	DRAWN BY L.B. FEJERAL
-------	-----------------------------	--------------------------

COULOMETER SYSTEM

SYSTEM

SYSTEM	NO. 1	W. P. AYRES	6-10-87
		D. V. RALPH	6-10-87

12-5-1961	W. B. JACOBSEN	6-11-61
12-5-1961	W. B. JACOBSEN	6-11-61

SCHEMATIC

ASSEMBLIES  
AND

AND  
DETAILS

DATE 7-2-87	SCALE AS SHOWN	FILE 20 / GO MECH/ELEC	S5-2-15952
----------------	-------------------	------------------------------	------------




RACK 7

POWER SUPPLY 7

The diagram shows a horizontal row of 12 vertical modules, labeled J1 through J12 from left to right. Each module has a rectangular shape with a small circle at the bottom. Modules J3, J9, and J12 contain a diagonal cross (X) inside, while the others are empty.

DETAIL (N) 0  
DIRECTOR  
FUNCTION TABLE



AMP CONNECTOR  
#202815-5  
(ALL PINS FEMALE)  
DETAIL (H)

CONNECTOR TOB N.I.M. POWER SUPPLY						
PIN NUMBER	POTENTIOSTAT #1 MODULE	INTERMATOR MODULE	COUNTER/TIMER MODULE	POTENTIOSTAT #2 MODULE	AUTOMATION MODULE	
	J3	J4	J6	J9	J12	
1						
2						
3						115 HOT (OUT)
4						
5						
6						
7						
8	HOLD #1					HOLD #1
9	OPERATE #1					OPERATE #1
10						
11			HOLD DISPLAY			HOLD DISPLAY
12			CLEAR DISPLAY			CLEAR DISPLAY
13	5V (IN)			5V (IN)		
14				HOLD #2		HOLD #2
15				OPERATE #2		OPERATE #2
16	12V (IN) →					12V (IN)
17	-12V (IN) →					-12V (IN)
18	5V (IN)		5V (IN)	5V (IN)		5V (IN)
19	5V (RETURN)		5V (RETURN)	5V (RETURN)		5V (RETURN)
20						
21	COUNTER #1					COUNTER #1
22	REFERENCE #1					REFERENCE #1
23	WORKING #1					WORKING #1
24				COUNTER #2		COUNTER #2
25				REFERENCE #2		REFERENCE #2
26				WORKING #2		WORKING #2
27						
28	24V (IN) →					24V (IN)
29	-24V (IN) →					-24V (IN)
30			SELECT A	SELECT A		
31			SELECT B	SELECT B		
32						
33	115 VAC (IN) →					115 VAC (IN)
34	COMMON			COMMON		
35		+15V (IN)				+15V (OUT)
36		15V RETURN (IN)				15V RETURN (OUT)
37		-15V (IN)				-15V (OUT)
38						
39						
40						
41	115 VAC (N) IN					115 VAC (N) IN
42	CHASSIS GND					CHASSIS GND


NOTE: 115 HOT IS NOT CONNECTED  
TO PIN 23 ON SYSTEM 1

\* NOTE: 5V (IN) IS NOT CONNECTED TO PINE 19 ON SYSTEM 1


[illegible]



Analyze this print for  
SAFETY CONSIDERATIONS




CONTROL AMPLIFIER BOARD  
SCHEMATIC




COMPONENT LAYOUT

ITEM	QTY	DESCRIPTION
19	1	SOCKET, I.C., 14 PIN, TELCO INSTRUMENTS, #62714-08
18	1	AMPLIFIER, TELETYPE, PHILBRICE, #1661
17	1	HEAT SINK, THERMALLOY, #6007A
16	1	CAPACITOR, 1000PF, WICA, RANDAM, #CM05FA103J08
15	1	RESISTOR, 5.1K, 1/2 WATT, 1% TOL.
14	1	TERMINAL, TURRET, HOLLOW LUG, 1/4 IN. DIA., STOCK # 2000A
13	1	CAPACITOR, 100PF, ELECTROLYTIC, PANASONIC, #EEF-B10V018
12	1	RESISTOR, 3.9K, 1/2 WATT, 5% TOL.
11	1	CAPACITOR, 100PF, WICA, 1% TOL.
10	1	RESISTOR, 100K, 1/2 WATT, 5% TOL.
9	1	AMPLIFIER, OVER-DRUM, #3351
8	1	SOCKET, I.C., 14 PIN, TELCO INSTRUMENTS, #62714-08
7	1	RESISTOR, 100K, 1/2 WATT, 5% TOL.
6	1	TRANSISTOR, 2N2900
5	1	CAPACITOR, 100PF, WICA, 1% TOL.
4	1	CAPACITOR, 100PF, WICA, 1% TOL.
3	1	RESISTOR, 100K, 1/2 WATT, 5% TOL.
2	1	RESISTOR, 5.1K, 1/2 WATT, 1% TOL.
1	1	BOARD, PRINTED CIRCUIT, POTENTIOMETER, PRINTED WIRING SERVICES

Drawing 7



SWITCHING AND LOAD BOARD  
SCHEMATIC



COMPONENT LAYOUT

ITEM	QTY	DESCRIPTION
11	1	COMMON, 5 VOLTS
10	1	COMMON, 5 VOLTS
9	1	COMMON, 5 VOLTS
8	1	COMMON, 5 VOLTS
7	1	COMMON, 5 VOLTS
6	1	COMMON, 5 VOLTS
5	1	COMMON, 5 VOLTS
4	1	COMMON, 5 VOLTS
3	1	COMMON, 5 VOLTS
2	1	COMMON, 5 VOLTS
1	1	COMMON, 5 VOLTS

THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOW HOW THEREIN MAY NOT BE USED FOR THE DESIGN OR CONSTRUCTION OF ANY EQUIPMENT WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS IN WHOLE OR IN PART REFER TO THIS STATE.		S5-2-15956 LATEST REVISION ON THIS DRAWING ②	
UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO., INC. SAVANNAH RIVER PLANT		SEPARATIONS TITLE COULOMETER SYSTEM SWITCHING AND LOAD AND CONTROL AMPLIFIER PRINTED CIRCUIT BOARDS SCHEMATIC AND COMPONENT LAYOUT	
7-2-87 7-2-87 7-2-87		7-2-87 7-2-87 7-2-87	



[illegible]

TECHNICAL DRAWING OF A RECTANGULAR STRUCTURE, POSSIBLY A CONTAINER OR TANK, WITH VARIOUS DETAILS AND DIMENSIONS.

**Labels and Callouts:**

- ①: 600 DETAIL (Top Left Corner)
- ②: 600 DETAIL (Top Right Corner)
- ③: 600 DETAIL (Bottom Left Corner)
- ④: 600 DETAIL (Bottom Right Corner)
- ⑤: 600 DETAIL (Top Right Corner, Alternative View)
- ⑥: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑦: 600 DETAIL (Top Right Corner, Alternative View)
- ⑧: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑨: 600 DETAIL (Top Right Corner, Alternative View)
- ⑩: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑪: 600 DETAIL (Top Right Corner, Alternative View)
- ⑫: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑬: 600 DETAIL (Top Right Corner, Alternative View)
- ⑭: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑮: 600 DETAIL (Top Right Corner, Alternative View)
- ⑯: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑰: 600 DETAIL (Top Right Corner, Alternative View)
- ⑱: 600 DETAIL (Bottom Right Corner, Alternative View)
- ⑲: 600 DETAIL (Top Right Corner, Alternative View)
- ⑳: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉑: 600 DETAIL (Top Right Corner, Alternative View)
- ㉒: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉓: 600 DETAIL (Top Right Corner, Alternative View)
- ㉔: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉕: 600 DETAIL (Top Right Corner, Alternative View)
- ㉖: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉗: 600 DETAIL (Top Right Corner, Alternative View)
- ㉘: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉙: 600 DETAIL (Top Right Corner, Alternative View)
- ㉚: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉛: 600 DETAIL (Top Right Corner, Alternative View)
- ㉜: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉝: 600 DETAIL (Top Right Corner, Alternative View)
- ㉞: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㉟: 600 DETAIL (Top Right Corner, Alternative View)
- ㊱: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊲: 600 DETAIL (Top Right Corner, Alternative View)
- ㊳: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊴: 600 DETAIL (Top Right Corner, Alternative View)
- ㊵: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊶: 600 DETAIL (Top Right Corner, Alternative View)
- ㊷: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊸: 600 DETAIL (Top Right Corner, Alternative View)
- ㊹: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊺: 600 DETAIL (Top Right Corner, Alternative View)
- ㊻: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊼: 600 DETAIL (Top Right Corner, Alternative View)
- ㊽: 600 DETAIL (Bottom Right Corner, Alternative View)
- ㊾: 600 DETAIL (Top Right Corner, Alternative View)
- ㊿: 600 DETAIL (Bottom Right Corner, Alternative View)

**Dimensions:**

- Overall Width: 7'0"
- Overall Height: 1'0"
- Internal Width: 6'0"
- Internal Height: 1'0"
- Top Flange Thickness: 1/8"
- Bottom Flange Thickness: 1/8"
- Side Flange Thickness: 1/8"
- Top Flange Width: 1/8"
- Bottom Flange Width: 1/8"
- Side Flange Width: 1/8"
- Top Flange Length: 1/8"
- Bottom Flange Length: 1/8"
- Side Flange Length: 1/8"
- Top Flange Area: 1/8"
- Bottom Flange Area: 1/8"
- Side Flange Area: 1/8"
- Top Flange Volume: 1/8"
- Bottom Flange Volume: 1/8"
- Side Flange Volume: 1/8"
- Top Flange Weight: 1/8"
- Bottom Flange Weight: 1/8"
- Side Flange Weight: 1/8"
- Top Flange Density: 1/8"
- Bottom Flange Density: 1/8"
- Side Flange Density: 1/8"
- Top Flange Temperature: 1/8"
- Bottom Flange Temperature: 1/8"
- Side Flange Temperature: 1/8"
- Top Flange Pressure: 1/8"
- Bottom Flange Pressure: 1/8"
- Side Flange Pressure: 1/8"
- Top Flange Flow: 1/8"
- Bottom Flange Flow: 1/8"
- Side Flange Flow: 1/8"
- Top Flange Level: 1/8"
- Bottom Flange Level: 1/8"
- Side Flange Level: 1/8"
- Top Flange Position: 1/8"
- Bottom Flange Position: 1/8"
- Side Flange Position: 1/8"
- Top Flange Orientation: 1/8"
- Bottom Flange Orientation: 1/8"
- Side Flange Orientation: 1/8"
- Top Flange Rotation: 1/8"
- Bottom Flange Rotation: 1/8"
- Side Flange Rotation: 1/8"
- Top Flange Translation: 1/8"
- Bottom Flange Translation: 1/8"
- Side Flange Translation: 1/8"
- Top Flange Scaling: 1/8"
- Bottom Flange Scaling: 1/8"
- Side Flange Scaling: 1/8"
- Top Flange Skewing: 1/8"
- Bottom Flange Skewing: 1/8"
- Side Flange Skewing: 1/8"
- Top Flange Distorting: 1/8"
- Bottom Flange Distorting: 1/8"
- Side Flange Distorting: 1/8"
- Top Flange Transforming: 1/8"
- Bottom Flange Transforming: 1/8"
- Side Flange Transforming: 1/8"

**Notes:**

- ALL DIMENSIONS ARE IN INCHES UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO CENTER UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO FACE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO EDGE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO CORNER UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO MIDDLE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO SURFACE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO INTERIOR UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO EXTERIOR UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO INSIDE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO OUTSIDE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO TOP UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO BOTTOM UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO LEFT UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO RIGHT UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO FRONT UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO BACK UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO SIDE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO END UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO HEAD UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO TAIL UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO POINT UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO LINE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO AREA UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO VOLUME UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO WEIGHT UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO DENSITY UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO TEMPERATURE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO PRESSURE UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO FLOW UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO LEVEL UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO POSITION UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO ORIENTATION UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO ROTATION UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO TRANSLATION UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO SCALING UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO SKEWING UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO DISTORTING UNLESS OTHERWISE SPECIFIED.
- ALL DIMENSIONS ARE TO TRANSFORMING UNLESS OTHERWISE SPECIFIED.

[illegible][illegible]

The technical drawing shows three views of a mechanical component:


- Front View (Top):** A horizontal bar with a total length of 2". The distance from the left end to the first hole is 1 1/8". The distance between the two holes is 1/8". The diameter of each hole is 16/32" O.D. with a tolerance of ±0.000". There are additional dimension lines indicating distances of 3/8" and 1/8" from the ends.
- Side View (Bottom):** Shows the profile of the bar with a thickness of 1/8".
- End View (Right):** A circular cross-section with a diameter of 16/32" O.D. and a tolerance of ±0.000". It also shows a feature labeled R6-ROUNDED.

Figure 1 is an elevation view of a rectangular plate. The overall width is 100 inches, indicated by a dimension line at the top. The overall height is 50 inches, indicated by a dimension line on the left. There are four holes, one in each corner. The holes are specified as 10 inches in diameter and 10 inches deep, indicated by a callout: "10" DIA. HOLES 10" D". The distance from the center of each hole to the nearest edge is 25 inches, indicated by dimension lines: 25" from the top and bottom edges to the center of the vertical holes, and 25" from the left and right edges to the center of the horizontal holes. The holes are represented by circles with a center dot.


Hand-drawn diagram of a rectangular structure. The top edge is labeled "10.00". The left edge is labeled "0.00". The right edge is labeled "10.00". The bottom edge is labeled "0.00". A dashed line connects the top-left corner to the bottom-right corner. A label "R4-BOUNCING" is written near the top-left corner. A label "4 PLACES, 10.00" is written near the bottom-left corner. A label "12.00" is written near the bottom-right corner. A label "10.00" is written near the top-right corner.

FORM NO.	REMARKS	DATE OF REVISE	APPROVED BY	DATE	<p><b>FRAGMENTS ALL</b> <b>MASSIVE SURFACES</b> <b>U.O.D</b></p> <p><b>BREAK ALL SHARP</b> <b>EDGES U.O.D</b></p> <p>THE NUMBER OF THIS FORM, (STANDARD) IS 1. THE NUMBER OF THIS FORM, (STANDARD) IS 1. THE NUMBER OF THIS FORM, (STANDARD) IS 1.</p>	<p><b>TOLERANCES</b> <b>USE ALL DIMENSIONS</b> <b>UNLESS OTHERWISE NOTED</b> <b>SHALL BE</b></p> <p>ONE PLACE DECIMAL, 1.00 TWO PLACE DECIMAL, 1.00 THREE PLACE DECIMAL, 1.000</p> <p>FRACTIONAL: 1/16 ANGULAR: 10° 30'</p>	<p>THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS &amp; CO. THE INFORMATION AND KNOW HOW THEREON MAY NOT BE USED FOR THE DRAWING REPRODUCED WITHOUT THE WRITTEN PERMISSION OF DU PONT DE NEMOURS &amp; CO. IN WASHINGTON, D. C. IF IN REPLY TO THIS DRAWING, PLEASE REFER TO THIS DRAWING.</p> <p>UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS &amp; CO. INC.</p>	<p><b>S5-2-15957</b></p> <p><b>LATEST REVISION</b> <b>ON THIS DRAWING</b></p>
					<p><b>BAVANNAH RIVER PLANT</b></p> <p><b>OPERATIONS</b></p>	<p>BLANK NO. <b>772-F</b> PROBLEM NO. <b>2-958</b></p> <p>FILED</p>	<p>PROJECT NO. <b>772-F</b> TARGET NO. <b>2-958</b></p>	
					<p><b>REFERENCE DRAWINGS</b></p> <p><b>051-19999POTENTIOSTAT MODULE 1997</b></p>	<p>TITLE</p> <p><b>COULOMETER SYSTEM</b></p> <p><b>POTENTIOSTAT MODULE</b></p>	<p>REVISIONS</p> <p>REVISION NO. <b>1</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>2</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>3</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>4</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>5</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>6</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>7</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>8</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>9</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>10</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>11</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>12</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>13</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>14</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>15</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>16</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>17</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>18</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>19</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>20</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>21</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>22</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>23</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>24</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>25</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>26</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>27</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>28</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>29</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>30</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>31</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>32</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>33</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>34</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>35</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>36</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>37</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>38</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>39</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>40</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>41</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>42</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>43</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>44</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>45</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>46</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>47</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>48</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>49</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>50</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>51</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>52</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>53</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>54</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>55</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>56</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>57</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>58</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>59</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>60</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>61</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>62</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>63</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>64</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>65</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>66</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>67</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>68</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J. GORDON</b></p> <p>REVISION NO. <b>69</b> REVISION DATE <b>1-12-77</b> REVISION BY <b>W. J</b></p>	



Analyze this print for  
SAFETY CONSIDERATIONS



- 57 -



FRONT VIEW




REAR VIEW

7	PI	COM L	1	CONNECTOR EDGE BOARD, TRW/CINCH #50-80A-80
6	01		1	SWITCH TOGGLEBAT LEVER ALLOY # MTA 5000
5	02		5	CONNECTOR DMC INSULATED MOUNT JAMBUOL #51-010
4	03		2	RESISTOR 100K, 1/4 WATT, 5% TOL
3	04	COM L	4	#4-80UNC-2A X 1/2 LONG FLAT HEAD MACHINE SCREW W/FLY NUT
2	05-15999		1	VOLTAGE-TO-FREQUENCY PC BOARD SUBASSEMBLY
1	06-21096		1	THIS Dwg. IS A "WIRE" MODIFIED (SEE DETAIL)
ITEM	PART	PART NO.	QTY	DESCRIPTION


Drawing 9



DESIGNED BY CHECKED BY APPROVED BY				FINISH / ALL MACHINED SURFACES UNLESS SPECIFIED SHALL BE BREAK ALL SHARP EDGES 0.005 R.		TOLERANCES UNLESS SPECIFIED SHALL BE FRACTIONAL: DECIMAL: ANGULAR:		THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOWLEDGE THEREON MAY NOT BE USED FOR THE DRAWING REPRODUCED WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS IN WHOLE OR IN PART INCLUDING VENDOR'S SHOP DRAWINGS, SHALL BEAR OR REFER TO THIS STAMP.		S5-2-15958 LATEST REVISION ON THIS DRAWING	
UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO., INC. SAVANNAH RIVER PLANT PROJECT NO. 772-R PROBLEM NO. 2-5988										DEPARTMENTS SEPARATIONS	
COULOMETER SYSTEM INTEGRATOR MODULE ASSEMBLY CONNECTION DIAGRAM AND PARTS LIST										DATE 7-1-87 BY NONE TYPED 10/10 MECH/ELC	
REFERENCE STANDARDS										S5-2-15958	

Analyze this print for  
SAFETY CONSIDERATIONS



- 59 -

[illegible]



Drawing 12

[illegible]



13	U2	C O M L	1	INTEGRATED CIRCUIT, P84T400A
11	U5, U6		2	INTEGRATED CIRCUIT, DATA SELECTOR/MULTIPLYER P84T400A
10	U9, U8		3	INTEGRATED CIRCUIT, DECODER/COUNTER/P84T400A
9	NF		4	RESISTOR NETWORK, 100.0
8	NEW TUBE		5	ELION-BRADLEY, MODEL 216 B
7	U11		6	SOCKET, INTEGRATED CIRCUIT 14 PINS, TEXAS INSTRUMENTS, CT 214-05
6	U10		7	INTEGRATED CIRCUIT, QUAD INVERTER, P84T400A
5	U5		8	INTEGRATED CIRCUIT, LOGIC TO 7 SEGMENT DECODER, 14 PINS
4	U5, U10, U11, U12, U13, U14, U15, U16, U17, U18, U19, U20, U21, U22, U23, U24, U25, U26, U27, U28, U29, U30, U31, U32, U33, U34, U35, U36, U37, U38, U39, U40, U41, U42, U43, U44, U45, U46, U47, U48, U49, U50, U51, U52, U53, U54, U55, U56, U57, U58, U59, U60, U61, U62, U63, U64, U65, U66, U67, U68, U69, U70, U71, U72, U73, U74, U75, U76, U77, U78, U79, U80, U81, U82, U83, U84, U85, U86, U87, U88, U89, U90, U91, U92, U93, U94, U95, U96, U97, U98, U99, U100		9	SOCKET, INTEGRATED CIRCUIT 14 PINS, TEXAS INSTRUMENTS, CT 214-05
3	U1		10	INTEGRATED CIRCUIT, DATA SELECTOR/MULTIPLYER P84T400A
2	U1		11	RESISTOR, 100.0, 1/4 WATT, 5% TOL.
1	U1		12	CAPACITOR, 100.0, ELECTROLYTIC, 50V
		C O M L	13	BOMED, PRINTED CIRCUIT MULTIVIBRATOR/DRIVER, PRINTED WIRING SERVICE
ITEM	PART	PART AND	QTY	DESCRIPTION

NOTE:

1. JUMPER IS TO BE ADDED TO THE PRINTED CIRCUIT BOARD FROM PIN 5 ON U6 TO GROUND.

[illegible]






6	U2 THRU U8S	COM L	4	INTEGRATED CIRCUIT, DUAL 401 LINE SELECTOR #NTELSION
			10	SOCKETS I.C., 16 PIN, NEHALES #761048
6	C1 THRU C7		1	CAPACITOR 0.056 CERAMIC
5	U1 THRU U3		2	NEHALES #761048
1	U4	COM L	21	INTEGRATED CIRCUIT 4 BIT COUNTER/LATCH
			2	SOCKETS, I.C., 36 PIN, NEHALES #761048
			1	BOARD PRINTED CIRCUIT
			1	PRINTED WIRING DEVICE
ITEM	PARTS LISTED	PART AND QUANTITY	REQD QTY	DESCRIPTION

[illegible]





Analyze this print for  
SAFETY CONSIDERATIONS



- 64 -

REAR VIEW

FRONT VIEW

(PARTS LIST CONTINUED BELOW)

Drawing 16

ITEM	PART	PART AND	QUANTITY	DESCRIPTION
NO.	DESC.	DRAWING NO.	REQ.	
15	C O M L	1	1	#4-ROUND-22 1/2" LONG PLAT HEAD
16	P S I	1	1	POWER SUPPLY, 25VDC,
17	C O M L	1	1	SOLE #8-18-1/2" DIA
18	C O M L	1	1	DRIVE/INTERFACE RC BOARD
19	C O M L	1	1	CONNECTOR, 200P BOARD,
20	C O M L	1	1	TYPE 1, #20-20A-50
21	C O M L	1	1	BATTERY, 1.55V,
22	C O M L	1	1	0-100, 1.55V,
23	C O M L	1	1	100000 BATTERY, FOR D-100 BATTERY,
24	C O M L	1	1	KEYSTONE, #10
25	C O M L	1	1	AUTOMATION RC BOARD
26	C O M L	1	1	CONNECTOR, 200P BOARD,
27	C O M L	1	1	TYPE 1, #20-20A-50
28	C O M L	1	1	COMPOUND MOUNT PLATE, MODIFIED,
29	C O M L	1	1	REF DETAIL
30	C O M L	1	1	CONNECTOR,
31	C O M L	1	1	AMPENOL, #MS102A15-10P
32	C O M L	1	1	AMPENOL, #MS102A15-10P
33	C O M L	1	1	AMPENOL, #MS102A15-10P
34	C O M L	1	1	AMPENOL, #MS102A15-10P
35	C O M L	1	1	AMPENOL, #MS102A15-10P
36	C O M L	1	1	AMPENOL, #MS102A15-10P
37	C O M L	1	1	AMPENOL, #MS102A15-10P
38	C O M L	1	1	AMPENOL, #MS102A15-10P
39	C O M L	1	1	AMPENOL, #MS102A15-10P
40	C O M L	1	1	AMPENOL, #MS102A15-10P
41	C O M L	1	1	AMPENOL, #MS102A15-10P
42	C O M L	1	1	AMPENOL, #MS102A15-10P
43	C O M L	1	1	AMPENOL, #MS102A15-10P
44	C O M L	1	1	AMPENOL, #MS102A15-10P
45	C O M L	1	1	AMPENOL, #MS102A15-10P
46	C O M L	1	1	AMPENOL, #MS102A15-10P
47	C O M L	1	1	AMPENOL, #MS102A15-10P
48	C O M L	1	1	AMPENOL, #MS102A15-10P
49	C O M L	1	1	AMPENOL, #MS102A15-10P
50	C O M L	1	1	AMPENOL, #MS102A15-10P
51	C O M L	1	1	AMPENOL, #MS102A15-10P
52	C O M L	1	1	AMPENOL, #MS102A15-10P
53	C O M L	1	1	AMPENOL, #MS102A15-10P
54	C O M L	1	1	AMPENOL, #MS102A15-10P
55	C O M L	1	1	AMPENOL, #MS102A15-10P
56	C O M L	1	1	AMPENOL, #MS102A15-10P
57	C O M L	1	1	AMPENOL, #MS102A15-10P
58	C O M L	1	1	AMPENOL, #MS102A15-10P
59	C O M L	1	1	AMPENOL, #MS102A15-10P
60	C O M L	1	1	AMPENOL, #MS102A15-10P
61	C O M L	1	1	AMPENOL, #MS102A15-10P
62	C O M L	1	1	AMPENOL, #MS102A15-10P
63	C O M L	1	1	AMPENOL, #MS102A15-10P
64	C O M L	1	1	AMPENOL, #MS102A15-10P
65	C O M L	1	1	AMPENOL, #MS102A15-10P
66	C O M L	1	1	AMPENOL, #MS102A15-10P
67	C O M L	1	1	AMPENOL, #MS102A15-10P
68	C O M L	1	1	AMPENOL, #MS102A15-10P
69	C O M L	1	1	AMPENOL, #MS102A15-10P
70	C O M L	1	1	AMPENOL, #MS102A15-10P
71	C O M L	1	1	AMPENOL, #MS102A15-10P
72	C O M L	1	1	AMPENOL, #MS102A15-10P
73	C O M L	1	1	AMPENOL, #MS102A15-10P
74	C O M L	1	1	AMPENOL, #MS102A15-10P
75	C O M L	1	1	AMPENOL, #MS102A15-10P
76	C O M L	1	1	AMPENOL, #MS102A15-10P
77	C O M L	1	1	AMPENOL, #MS102A15-10P
78	C O M L	1	1	AMPENOL, #MS102A15-10P
79	C O M L	1	1	AMPENOL, #MS102A15-10P
80	C O M L	1	1	AMPENOL, #MS102A15-10P
81	C O M L	1	1	AMPENOL, #MS102A15-10P
82	C O M L	1	1	AMPENOL, #MS102A15-10P
83	C O M L	1	1	AMPENOL, #MS102A15-10P
84	C O M L	1	1	AMPENOL, #MS102A15-10P
85	C O M L	1	1	AMPENOL, #MS102A15-10P
86	C O M L	1	1	AMPENOL, #MS102A15-10P
87	C O M L	1	1	AMPENOL, #MS102A15-10P
88	C O M L	1	1	AMPENOL, #MS102A15-10P
89	C O M L	1	1	AMPENOL, #MS102A15-10P
90	C O M L	1	1	AMPENOL, #MS102A15-10P
91	C O M L	1	1	AMPENOL, #MS102A15-10P
92	C O M L	1	1	AMPENOL, #MS102A15-10P
93	C O M L	1	1	AMPENOL, #MS102A15-10P
94	C O M L	1	1	AMPENOL, #MS102A15-10P
95	C O M L	1	1	AMPENOL, #MS102A15-10P
96	C O M L	1	1	AMPENOL, #MS102A15-10P
97	C O M L	1	1	AMPENOL, #MS102A15-10P
98	C O M L	1	1	AMPENOL, #MS102A15-10P
99	C O M L	1	1	AMPENOL, #MS102A15-10P
100	C O M L	1	1	AMPENOL, #MS102A15-10P

(PARTS LIST CONTINUED)





FIG. NO.	REVISION	DESIGNED BY	CHECKED BY	APPROVED BY	DATE	PROJECT / ALL	WORKSHEETS	THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOW HOW THEREON MAY NOT BE USED FOR THE DRIVING OF THE PROJECT WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS IN WHOLE OR IN PART INCLUDING VECTOR'S SHOP DRAWINGS SHALL BEAT ON REFER TO THIS STAMP.	S5-2-15965
						UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO. INC.	LATEST REVISION ON THIS DRAWING		
						SAVANNAH RIVER PLANT	SEPARATIONS		
						FIG. NO. 772-F	INSTRUMENT	2-3888	
						65-1-7506 AUTOMATION RC BOARD	65-1-7507 DRIVE/INTERFACE RC BOARD		
						COULOMETER SYSTEM AUTOMATION MODULE ASSEMBLY CONNECTION DIAGRAM AND PARTS LIST			
						DATE 7-1-87	BY NONE	CHK TO GC MCK-ELC	S5-2-15965




4	EL PASO, CULVER	C O M L	7	RELAT. FORM C BYDC C.F. SLABE, # HMMBRIIIII 800
5		K B	2	RELAT. FORM C BYDC C.F. SLABE, # HMMBRIIIII 500
6	B , R2		2	RESISTOR, 100Ω, .0% STABILITY, .005%
7	E		1	VOLIE RESISTORS BOARD, PRINTED CIRCUIT, POTENTIOMETER, PRINTED WIRING SERVICES
ITEM	PARTY ORDER	PART AND DRAWING NO.	QTY REQD	DESCRIPTION

Drawing 17

[illegible]




Analyze this print for  
SAFETY CONSIDERATIONS




**DETAIL (M)**  
(FRONT PANEL, MODIFIED)  
MATERIAL: NIM D1N  
AL


**ITEM (1)** (55-2-15968)



**DETAIL (L)**  
(FRONT PANEL, MODIFIED)  
MATERIAL: NIM D1N  
AL



**DETAIL (K)**  
(FRONT PANEL, MODIFIED)  
MATERIAL: NIM D1N  
AL




**DETAIL (N)**  
(COMPONENT MOUNTING BRACKET)  
MATERIAL: 504-L-637  
REF'D:

**Drawing 19**


<b>PROJ. NO.</b> 55-2-15968		<b>REVISED BY</b> DATE		<b>APPROVED BY</b> DATE		<b>FORMED / ALL</b> MACHINED SURFACES U.S.S. BREAK ALL SHARP EDGES & CORNERS R0.05 U.S.S. THE NUMBER OF THE DRAWING SYMBOLS OF THE MAXIMUM EXCEEDING THE MAXIMUM VALUE IN A 2.00 INCHES		<b>TOLERANCES</b> ON ALL DIMENSIONS UNLESS OTHERWISE SPECIFIED SHALL BE: ONE PLACE DECIMAL: 1.00 TWO PLACE DECIMAL: 0.01 THREE PLACE DECIMAL: 0.001 FRACTIONAL: 1/16 DECIMAL: 0.001		THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND KNOWLEDGE THEREON MAY NOT BE USED FOR THE CLAIMING PURPOSES WITHOUT THE WRITTEN PERMISSION OF DU PONT. ALL REPRODUCTIONS OF THIS DRAWING, INCLUDING VENDOR'S SHOP DRAWINGS, SHALL BEAT OR REFER TO THIS STAMP.		<b>55-2-15968</b> LATEST REVISION ON THIS DRAWING	
UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO., INC.						SAVANNAH RIVER PLANT SEPARATIONS TITLE: COULOMETER SYSTEM AUTOMATION MODULE							
REFERENCE DRAWINGS 55-2-15968 AUTOMATION MODULE						REFERENCE STANDARDS DATE: 7-1-67 REV: 12-1-10 FILE NO. 55-2-15968 WORKMAN: ELC							

Analyze this print for  
SAFETY CONSIDERATIONS


- 68 -




W1  
(DCD CABLE)





W2  
(OPIC INTERFACE CABLE)



W3  
(+5V POWER SUPPLY CABLE)  
NOTE: W3 IS FOR SYSTEM 2 ONLY.  
W3 IS 2-10000 FOR SYSTEM 1'S  
5V POWER SUPPLY



W4  
\* NOTE: CONNECT SHIELD TO AMPHENOL CONNECTOR



CABLE	LENGTH
W5	0'
W6	24' (MAX)
W7	0'
W8	0'
W9	10'
W10	10'
W11	24' (MAX)
W12	0'
W13	0'
W14	0'
W15	0'
W16	0'
W17	0'
W18	0'
W19	0'

\*\* DO NOT EXCEED MAXIMUM LENGTH

Drawing 20

REVISIONS 1. REVISED TABLE 2. REVISED WIRING DIAGRAM	DESIGNED BY DRAWN BY CHECKED BY APPROVED BY	FROM: ALL TO: SAVANNAH RIVER PLANT DATE: 7-1-67 BY: NCB	THIS DRAWING HAS BEEN FURNISHED BY E. I. DU PONT DE NEMOURS & CO. THE INFORMATION AND DATA HEREON ARE NOT TO BE USED FOR THE REPRODUCTION OF ANY PART OF THE INFORMATION OR DATA HEREON WITHOUT THE WRITTEN PERMISSION OF E. I. DU PONT DE NEMOURS & CO.	S5-2-15969 LATEST REVISION ON THIS DRAWING
UNITED STATES DEPARTMENT OF ENERGY E. I. DU PONT DE NEMOURS & CO. INC. SAVANNAH RIVER PLANT PROJECT NO. 772-R PROBLEM NO. 2-3888			SEPARATIONS APPROVED BY APPROVED BY APPROVED BY APPROVED BY	
COULOMETER SYSTEM CABLE DETAILS			DATE: 7-1-67 BY: NCB GO-ELEC S5-2-15969	

## APPENDIX B

### VENDOR LIST

- 1) ACD-Elcor Products, Tel: 203-734-3510  
75R Wooster St.  
P. O. Box 244  
Bethel CT 06801
- 2) Analog Devices, Tel: 617-329-4700  
Two Technology Way  
P. O. Box 280  
Norwood MA 02062
- 3) Arrow Electronics, Tel: 919-725-8711  
938 Burke Street  
Winston Salem, NC 27101
- 4) Augat/Alcoswitch, Tel: 617-821-2372  
Canton, MA 02021
- 5) Burr-Brown, Tel: 602-746-1111  
International Airport Industrial Park  
P. O. Box 11400  
Tucson AR.
- 6) Digi-Key Corporation, Tel: 800-344-4539  
701 Brooks Ave. South  
P. O. Box 677  
Thief River Falls, MN 56701
- 7) Dixie Electronics, 404-722-2055  
P. O. Box 2307  
Augusta GA
- 8) Grayhill Inc., Tel: 312-354-1040  
561 Hillgrove Ave.  
P. O. Box 10373  
Lagrange IL 60525
- 9) Hamilton Avnet, Tel: 404-447-7500  
5825 Suite D  
Peachtree Corners East  
Norcross GA 30092
- 10) Hewlett Packard, Tel: 803-732-0400  
Brookside Park 1 Harbinson Way  
Columbia SC

VENDORS LIST CONTINUED

- 11) Julie Research, Tel: 212-245-2727  
211 West 61 St.  
New York, NY 10023
- 12) Newark Electronics, Tel: 404-448-1300  
6950 Peachtree Industrial Blvd.  
Norcross GA 30071
- 13) Printed Wiring Services, 815-838-0005  
11 Maryknoll Drive  
Lockport, IL 60441
- 14) Radio Shack, Tel: 404-860-5965  
3830 Washington Rd.  
Martinez GA
- 15) Sorenson Company, Tel: 603-668-4500  
676 Island Pond Road  
Manchester, NH 03103
- 16) Teledyne Philbrick, Tel: 617-329-1600  
Allied Dr. @RT128  
Dedham, MA 02026
- 17) Tenneloc Corporation, Tel: 615-483-8405  
601 Oak Ridge Turnpike  
Oak Ridge, TN 37830
- 18) Texmate Inc., Tel: 619-481-7177  
348 South Cedro Ave.  
Solana Beach, CA 92075
- 19) Thermalloy Inc., Tel: 214-243-4321  
2021 W Valley View LN  
P. O. Box 810839  
Dallas TX 75381
- 20) Vectron Labs (C/O Saber Associates) Tel: 800-327-0853  
P. O. Box 1599  
Melborn, FL

# APPENDIX C

## PARTS LIST

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) NIM Bin, Tennenloc No. TB3	1	450.00	17
2) NIM power supply, +/-12, +/-24 VDC , Tennenloc No. TC911.	1	515.00	17
3) Three Wide Blank NIM Bin Panel, Tennenloc No. NIM-BP3	1	25.00	17
4) NIM extender cable, Tennenloc No. NC-EXC-4-S, Add pins 3 to 3, 18 to 18, 19 to 19, 35 to 35, 36 to 36, 37 to 37.	1	110.00	17
5) Shielded Cable, 3 conductor, 20 AWG, Belden No. 9963	100 FT. (min)	62.45	12
6) Connector, Amphenol 3106A-10SL-3S	1	7.55	12
7) Cable Clamp, Amphenol 97-3057-1004	1	2.68	12
8) Power Supply, 5 VDC, 6 Amps Sorenson No. SLC5-6B	1	185.00	15
9) Banana Plug, White, Pomona No. 1325-9	1	.90	12
10) Banana Plug, Green, Pomona No. 1325-5	1	.90	12
11) Banana Plug, Red, Pomona No. 1325-2	1	.90	12
12) Coax Cable, RG 58, 3 inches length plus connector.	3	12.91	12
13) Coax Cable, RG 58, 13 inches length plus connector.	1	13.00	12
14) Coax Cable, RG 58, 5 inches length plus connector.	1	12.91	12
15) Coax shorting plug.	1	6.41	12
16) Contact Sockets for Standard NIM Bin, Tennenloc No. 2100-03-07	11	.45	17



-2-

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
17) HP 9826S Computer with 1 Meg RAM	1	8360.00	10
18) HP 98622A GPIO Interface	1	345.00	10
19) HP 98623A BCD Interface	1	40.00	10
20) HP 98613A Ram Based Basic 3.0	1	653.00	10
21) HP 8116A Pulse/Function Generator	1	3307.00	10
22) HP 5316A Electronic Counter	1	3307.00	10
23) HP 3456A Programmable DVM	1	3619.00	10
24) HP 2673A Graphics Printer	1	1500.00	10
25) HP 59501B DAC/POWER SUPPLY	1	760.00	10
26) Potentiostat Module	2	1042.02	
27) Digital Integrator Module	1	759.86	
28) Counter Module	1	510.15	
29) Automation Module	1	752.43	

TOTAL SYSTEM COST: 27,433.86

Note: Prices based on 1986 catalogs.

POTENTIOSTAT PARTS LIST

1.1 CHASIS

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Switch, ALCOSWITCH Model MTA 206N.	2	3.66	4
2) Switch, ALCOSWITCH Model MTA 106F.	1	3.93	4
3) Switch, Rotary, by Grayhill Model No. 44D30-01-4-AJN	1	25.62	8
4) Potentiometer, 2K Ohm, Helipot Series A 10 Turn, Wire Wound, Linearity .25	1	23.05	3
5) Dial for Potentiometer, Series A	1	11.69	3
6) Potentiometer, 20K Ohm, 10 Turn, Spectrol No. 536-1-1-203	1	6.27	12
7) Connector, TRW Cinch, No. 50-22A-20	1	5.12	3
8) Connector, TRW Cinch, No. 50-30A-30	1	3.78	3
9) Connector, TRW Cinch, No. 50-20A-30M-1	1	6.07	3
10) Battery, 1.5 Volt, Alkaline, Size AA	4	.25	12
11) Holder, Battery, holds four AA Batteries, Radio Shack No. 270-397	1	1.79	12
12) Knob, for 1/4 " shaft Radio Shack No. 274-407	1	2/1.29	12
13) Connector, Amphenol 3102A-10SL-3P	1	5.09	12
14) Connector, BNC, Amphenol 31-010	3	2.89	12
15) Voltage Source, 6.2V, ACD-Elcor, No. AR6.2-10	1	230.00	1
16) Voltage Meter, Texmate No. PM-45X, 200mV Range.	1	150.00	18

17) Blank NIM Module with NIM connector Tenneloc No. NIM-B3W	1	75.00	17
18) Jack, Banana, Red, Newark No. 81N1228	1	.46	12
19) Jack, Banana, Black, Newark No. 81N1229	1	.46	12
20) Capacitor, Polystyrene, 100nF, Mallory No. SXX010	1	1.26	12
21) Control Amplifier Card	1	355.15	
22) Switching and Load Card	1	119.64	

SUB-TOTAL: 1042.02

POTENTIOSTAT MODULE PARTS LIST CONT.

1.2 Control Amplifier Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Amplifier, Teledyne Philbrick No. 1461	1	243.00	16
2) Heat Sink, Thermaloy No. 6007A	1	2.34	19
3) Amplifier, Burr-Brown No. 3521L	1	72.40	5
4) Trim Pot, 5K Ohm, Spectrol No. 70Y502	1	4.52	12
5) Trim Pot, 10K Ohm, Spectrol No. 70Y103	1	4.52	12
6) Transistor, 2N3906	1	.25	12
7) Resistor, 3.9 Ohm, 1/4 Watt, 5%	2	.05	6
8) Resistor, 3.3M Ohm, 1/4 Watt, 5%	1	.05	6
9) Resistor, 820K Ohm, 1/4 Watt, 5%	1	.05	6
10) Resistor, 5.1K Ohm, 1/2 Watt, 1%	1	.67	12
11) Capacitor, 100pF, Mica Sangamo No. CM05FD101J03	1	.34	12
12) Capacitor, 1000pF, Mica Sangamo No. CM05FA102J03	1	.63	12
13) Capacitor, .01uF, ceramic Sprague No. 1C10Z5U103M050B	4	.09	12
14) Capacitor, 1uF, Electrolytic, 50WV Panasonic No. ECE-B1HV010S	4	.18	6
15) Capacitor, 100uF, Electrolytic, 50WV Panasonic No. ECE-B1HV101S	2	.45	6
16) Control Amplifier Circuit Board Printed Wiring Services	1	12/250.00	13
17) Integrated Circuit Socket, 14 pin Texas Instruments C7214-09	1	1.66	7
18) Integrated Circuit Socket, Round, 8 pins, Augat No. 8059-2G5	1	1.81	12

SUB-TOTAL: 355.15

POTENTIOSTAT PARTS LIST CONT.

## 1.3 Switching and Load Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Resistor, 1.0 Ohm, 25 Watts Dale No. RH-25	2	3.53	12
2) Resistor, 8.0 Ohm, 25 Watts Dale No. RH-25	1	3.53	12
3) Resistor, 4.0 Ohm, 25 Watts Dale No. RH-25	1	3.53	12
4) Resistor, 82K Ohm, 1/4 Watts, 5%	1	3.53	6
5) Capacitor, 1000pF, Mica Sangamo No. CM05FA102J03	1	.63	12
6) Diode, Switching, SI, PIV 75V, 10Ma 1N4151 or equiv.	2	.12	
7) Relay, Form C, 5 Volts C. P. Clare No. HGM51111K00	3	15.79	7
8) Relay, 2 Form C, 5 Volts C. P. Clare No. HGW2MT54111G00	1	32.92	7
9) Switching Load Card Circuit Board, Printed Wiring Services.	1	12/250.00	13

SUB-TOTAL: 119.64

DIGITAL INTEGRATOR PARTS LISTS

## 2.1 Chasis

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Switch, 3PDT, Alco Model MTA-306D	1	6.33	3
2) Resistor, 100 Ohm, 1/4 Watt, 5%	2	.05	6
3) Connector, Amphenol 31-010	5	2.89	12
4) Connector, TRW Cinch No. 50-30A-30	1	3.78	3
5) Blank NIM Module with NIM connector Tenneloc No. NIM-B1W	1	65.00	17
6) Voltage-to-Frequency Card	1	670.20	

SUB-TOTAL: 759.86

DIGITAL INTEGRATOR PARTS CONT.

2.2 Voltage-To-Frequency Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Voltage-To-Frequency Convertor Module, Analog Devices No. 458L	2	182.00	2
2) Crystal Oscillator, 10K Hz, Vectron Labs No. C0235-T-3.	1	200	20
3) Voltage Regulator, uA7805CKC	1	1.43	7
4) Operational Amplifier, 741, 8 pin dip	2	.65	7
5) Quad Invertor, SN7404N	2	.45	7
6) Relay, Form C, 5 Volts, C. P. Clare No. HGW51111K00	1	15.79	7
7) Trim Pot, 500 Ohm, Spectrol No. 70Y501	1	4.95	12
8) Trim Pot, 1K Ohm, Spectrol No. 70Y102	2	4.52	12
9) Trim Pot, 10K Ohm, Spectrol No. 70Y103	2	4.52	12
10) Trim Pot, 20K Ohm, Spectrol No. 70Y203	2	4.52	12
11) Trim Pot, 50K Ohm, Spectrol No. 70Y503	2	4.52	12
12) Resistor, 5.1K Ohm, 1/4 Watt, 5%	1	.05	6
13) Resistor, 11K Ohm, 1/4 Watt, 5%	2	.05	6
14) Resistor, 43K Ohm, 1/4 Watt, 5%	1	.05	6
15) Resistor, 180K Ohm, 1/4 Watt, 5%	1	.05	6
16) Resistor, 221K Ohm, 1/4 Watt, 1%	2	.10	6
17) Resistor, 240K Ohm, 1/4 Watt, 5%	1	.05	6

18) Digital Integrator Card, Printed Wiring Service.	1	6/235.00	13
19) I. C. Socket, 8 pin dip Texas Instruments No. C7208-09	2	1.17	12
20) I. C. Socket, 16 pin dip Texas Instruments No. C7216-09	2	1.83	12

SUB-TOTAL: 670.20



COUNTER MODULE PARTS LIST

3.1 Chasis

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Switch, Alcoswitch No. MPA-106F.	2	3.93	3
2) Switch, Alcoswitch No. MTA-206TA.	1	3.66	3
3) Switch, Alcoswitch No. MTA-106D	1	2.68	3
4) Switch, Alcoswitch No. MST-405N	1	7.80	3
5) Connector, TRW Cinch 50-50A-30	1	6.30	3
6) Connector, TRW Cinch 50-44A-30	1	5.12	3
7) Connector, BNC, Amphenol 31-010	3	2.89	12
8) Connector, Micro-Ribbon, Amphenol No. 57-40500	1	2.65	12
9) LED Display, 4 digit, Hewlett Packard No. HP5082-7414	2	3.00	10
10) Blank NIM Module with NIM connector, Tenneloc No. NIM-B2W.	1	65.00	17
11) Resistor, 1K Ohm, 5%, 1/4 Watt	1	.05	6
12) Display/Driver Card	1	64.76	
13) Counter/Multiplexer Card	1	329.60	

SUB-TOTAL: 510.15

COUNTER MODULE PARTS LIST CONT.

3.2 Display/Driver Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Integrated Circuit, No. SN74LS151N Data Selector/ Multiplexer	4	.62	7
2) Integrated Circuit, No. SN74LS138N Decoder/ Demultiplexer.	1	.57	7
3) Integrated Circuit, No. CD4511BE BCD to 7 segment decoder/driver	1	1.43	7
4) Integrated Circuit, No. SN74LS293N 4 - bit binary counter.	2	1.13	7
5) Integrated Circuit, No. SN7404N Quad invertor.	1	.45	7
6) Integrated Circuit, No. SN7430N 8 - input NAND gate.	1	.45	7
7) Capacitor, 47uF, Electrolytic, 50 WV.	1	.36	6
8) Resistor, 1K Ohm, 1/4 Watt, 5%	1	.05	6
9) Resistor Network, 100 Ohm, Allen Bradley No. 316B	1	2.10	12
10) I. C. Socket, 16 pins, Texas Instrumnets C7216-09	8	1.83	12
11) I. C. Socket, 14 pins, Texas Instruments C7214-09	4	1.66	12
12) Display/Driver Circuit Board, Printed Wiring Service.	1	6/200.00	13

SUB-TOTAL: 64.76

COUNTER MODULE PARTS LIST CONT.

3.3 Counter/Multiplexer Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Integrated Circuit, SN74143, 4 - bit Counter/Latch, 7-seg driver	21	8.90	7
2) Integrated Circuit, SN74LS153N, Dual 4-to-1 Line Selector	14	.66	7
3) Capacitor, .01uF, ceramic, Sprague No. 1C10Z5U103M050B	7	.09	12
4) I. C. Socket, 16 pins, Texas Instruments No. C7216-09	14	1.83	12
5) I. C. Socket, 24 pins, Texas Instruments No. C7224-09	21	3.24	12
6) Counter/Multiplexer Card, Printed Wiring Service.	1	6/235.00	13

SUB-TOTAL: 329.60

### AUTOMATION MODULE PARTS LIST

#### 4.1 Chasis

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Power Supply, +/- 15 VDC, Sola No. 85-15-2150	1	227.63	12
2) Battery, Mercury, D size, 1.35 V, Duracell No. RM42R	1	13.20	12
3) Switch, ALCOSWITCH Model MTA-206N	1	3.66	4
4) Connector, BNC, Amphenol 31-010	5	2.89	12
5) Connector, Amphenol 3102A-10SL-3P	1	5.09	12
6) Connector, TRW Cinch No. 50-44A-30	1	5.12	3
7) Connector, TRW Cinch No. 50-36A-30	1	4.35	3
8) Connector, Amphenol 3102A-22-14P	1	9.55	12
9) Lamp, Dialight No. 507-4538-0931-610	1	4.83	12
10) Lamp Holder, Dialight No. 250-8745-14-504	1	3.38	12
11) Battery Holder for size D battery, Keysone No. 175	1	.81	12
12) Fuse Holder, Panel Mounting, Little Fuse No. 342001	1	1.15	12
13) Fuse, 3A, 3 AG	1	.38	12
14) Blank NIM Module with NIM connector, Tenneloc No. NIM-B3W.	1	75.00	17
15) Automation Card	1	317.13	
16) Receiver/Driver Card	1	66.70	

SUB-TOTAL: 752.43

AUTOMATION MODULE PARTS LIST CONT.

4.2 Automation Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Relay, Form C, 5VDC C. P. Clare No. HGWM51111K00	2	15.79	7
2) Relay, 2 Form C, 5VDC C. P. Clare No. HGW2MT54111G00	7	32.92	7
3) Resistor, 100 Ohm, .01%, stability .003% per year, Julie Research No. R33	2	10.89	11
4) Automation Card Printed Circuit Board, Printed Wiring Services.	1	6/200.00	13

SUB-TOTAL: 317.13

4.3 Receiver/Driver Card

<u>DESCRIPTION</u>	<u>QUANTITY</u>	<u>COST, EA</u>	<u>VENDOR</u>
1) Integrated Circuit, SN7404N, Hex inverter.	7	.45	7
2) 50 Ohm Line Driver, SN74128N	5	1.13	7
3) Resistor Network, Allen Bradley No. 408E221331	3	1.26	12
4) Capacitor, Electrolytic, 1uF, 50WV Panasonic No. ECE-B1HV010S	1	.18	6
5) Capacitor, Electrolytic, 10uF, 50WV Panasonic No. ECE-B1HV100S	1	.24	6
6) Capacitor, Electrolytic, 100uF, 50WV Panasonic No. ECE-B1HV101S	1	.45	6
7) Receiver/Driver Circuit Board, Printed Wiring Service	1	6/200.00	13
8) I. C. Socket, 14 pins Texas Instruments C7214-09	12	1.66	12

SUB-TOTAL: 66.70