


Novel Miniature Spectrometer for Remote Chemical Detection

Dr. Andrew C. R. Pipino
Physical and Chemical Properties Division
Chemical Science and Technology Laboratory
National Institute of Standards and Technology
Gaithersburg, Maryland 20899
Phone: (301) 975-2565
Fax: (301) 975-3672
Email: andrew.pipino@nist.gov
EMSP Contract#60231
(Jeff W. Hudgens, former PI)

ABSTRACT

A novel miniature spectrometer is being developed that permits ultra-high sensitivity chemical detection in remote or hazardous environments. The technology employs a miniature total-internal-reflection-ring optical cavity to extend the cavity ring-down concept to condensed matter for the first time. A detailed theoretical model and two levels of prototyping have been successfully completed. The ultimate potential of the technology approaches single molecule detection at ambient temperature.

“Conventional” Cavity Ring-Down Spectroscopy




- ◆ New technique for measurement of optical absorption
- ◆ $(L_{\text{abs}})_{\text{min}} = L_{\text{mirrors}} * \text{minimum detectable } \Delta\tau/\tau$
- ◆ **Gas-phase diagnostics only**


Extension to Condensed Matter

Use Total Internal Reflection


- ◆ Requires $\Theta_{\text{inc}} > \Theta_c = \sin^{-1}(n_2/n_1)$
- ◆ Ultra-high reflectivity
- ◆ Broadband
- ◆ Evanescent wave defines the sampled path length
- ◆ Can be used to probe adlayers, films, or bulk media



Initial Prototype Design




$$\tau(\omega) = \frac{t_r}{L_{\text{mirrors}} + L_{\text{bulk}} + L_{\text{surf}} + L_{\text{pol}} + L_{\text{abs}}}$$




- ◆ First application of CRDS to a condensed phase
- ◆ Intra-cavity Pellin-Broca prism with ultra-smooth surfaces
- ◆ Low-loss (high finesse) for only one polarization component

Initial Prototype Results



- ◆ Surface adsorbed I_2 used for calibration
- ◆ Detection limit: 4% of a monolayer at 625 nm
- ◆ Peak absorption at 520 nm for this moderate absorber

Optimum Cavity Design




$$\tau(\omega) = \frac{t_r}{L_{\text{bulk}} + L_{\text{surf}} + L_{\text{diff}} + L_{\text{coup}} + L_{\text{abs}}}$$

Intrinsic losses

A synergistic combination


- ◆ Ultra-high sensitivity
- ◆ Broadband
- ◆ Wide T and P range
- ◆ Single point sampling
- ◆ Miniaturizable
- ◆ Potential on-chip sensor

Theoretical Model



- ◆ Permits the optimum cavity size to be identified
- ◆ Optimum size arises from a competition between losses

The Experiment




- ◆ Gaps controlled piezoelectrically/monitored interferometrically
- ◆ Detection limit of ~40 ppm of a monolayer of I_2
- ◆ Final form can be packaged w/fiber optics

Other Geometries

For fused-silica cavities (visible and near-IR) :

- ◆ Square geometry is optimal for thin-film diagnostics
- ◆ Octagonal geometry is optimal for bulk aqueous solutions



Other materials exist for other spectral regions:

- ◆ Fluoride and chalcogenide glasses for the mid-IR
- ◆ Sapphire for the UV

Summary of Achievements

- ◆ A fundamentally new technology for chemical sensing has been:
 - ◆ conceived
 - ◆ elucidated theoretically through a general model
 - ◆ realized experimentally (quantitative agreement with theory!)
 - ◆ patented and published
- ◆ The technology has important properties for the EMSP mission:
 - ◆ spectroscopic detection/characterization
 - ◆ sensitivity approaching single molecule detection
 - ◆ miniaturizable
 - ◆ broadband
 - ◆ rugged and chemically inert
 - ◆ wide temperature and pressure range
 - ◆ can be remotely located with optical fiber

Publications and Patents

Evanescence wave cavity ring-down spectroscopy with a total-internal-reflection-ring minicavity,
A. C. R. Pipino et al., Rev. Sci. Instrum. **68** (8), 2978, (1997).

Evanescence wave cavity ring-down spectroscopy as a probe of surface processes,
A. C. R. Pipino et al. , Chem. Phys. Lett. **280**, 104, (1997).

Condensed matter cavity ring-down spectroscopy with miniature optical cavity,
A. C. R. Pipino, in preparation.

Intra-cavity total reflection for high sensitivity measurement of optical properties,
A. C. R. Pipino et al., U. S. Patent Application Serial No. 08/962,170, pending.

Broadband intra-cavity total reflection chemical sensor,
A. C. R. Pipino, U. S. Patent Application Serial No. 08/962,171, accepted.