UCRL 78
ey 69/A

UNIVERSITY OF CALIFORNIA

Radiation Laboratory

TWO-WEEK LOAN COPY

This is a Library Circulating Copy which may be borrowed for two weeks. For a personal retention copy, call Tech. Info. Division, Ext. 5545

BERKELEY, CALIFORNIA

UCRL-78

Special Review of
Declassified Reports

Authorized by USDOE JK Bratton
Unclassified TWX P182206Z May 79

REPORT PROPERLY DECLASSIFIED

AUG 15 1979

Date
Date



Lawrence Radiation Laboratory Library
University of California, Berkeley

UNIVERSITY OF CALIFORNIA RADIATION LABORATORY

Cover Sheet
Do not remove

DECLASSIFIED

This document contains pages and plates of figures.
This is copy of of 76. Series.

Issued to:



Each person who received this document must sign the cover sheet in the space below.

Route to	Noted by	Date	Route to	Noted by	Date
				<i>j</i>	
	***				itas v ang pto y taman to it a yan pana
	***		ng a magalatan ay ang bahan ay a		

:					-
Mariada Gales, alanmaydo - Falisa, escribant					
-					
Barrinen aus et erseta k _{arrin} ako - <mark>mekalan urabarbilana</mark> n					
Marie Marie Colorina de Marie Compte de la Marie de Marie	A strip Scott Annual Property of the Strip				

UNIVERSITY OF CALIFORNIA Radiation Laboratory

Contract No. W-7405-Eng-48



DECLASSIFIED

CLASSIFICATION CANCELLED

MARCING ... CHANGE

SPECTRAL LINES OF CURIUM FROM 3100 A TO 4200 A

John G. Conway and Milton F. Moore

April 7, 1948

Special Review of Declassified Reports

Authorized by USDOE JK Bratton Unclassified TWX P182206Z May 79

REPORT PROPERLY DECLASSIFIED

Authorized D: Classifier Date

8-15-79

8-20-79

Date

Berkeley, California

DECLASSIFIED

UCRL-78 Physics-General



Standard Distribution

	*
Argonne National Laboratory Atomic Energy Commission, Washington Brookhaven National Laboratories Carbide & Carbon Chemicals Corp. (K-25 Area) Carbide & Carbon Chemicals Corp. (Y-12 Area) Oak Ridge National Laboratories General Electric Company Hanford Engineer Works Iowa State College Los Alamos Office of N. Y. Directed Operations	1-10 11-13 14-21 22-23 24-25 26-33 34-37 38-39 40 41-43
Monsanto Chemical Company	45
National Bureau of Standards Patent Advisor	46-47 48
Library Branch; Oak Ridge, (NEPA)	49
Inbrary Branch, Oak Ridge University of Rochester Office of Chicago Directed Operations	50–64 65 66
University of California, Radiation Laboratory Information Division Chemistry Department G. T. Seaborg J. G. Conway	67 – 69 70 – 74 75 76

3.

SPECTRAL LINES OF CURIUM FROM 3100 A TO 4200 A

John G. Conway and Milton F. Moore

April 7, 1948



ABSTRACT

Fifty-four spectral lines of curium are listed together with estimated intensities. A brief description of the manner of taking the spectra and of measuring the lines is given.

SPECTRAL LINES OF CURIUM FROM 3100 A TO 4200 A John G. Conway and Milton F. Moore

April 7, 1948

Experimental

Three separate samples of curium, in all seven micrograms, were submitted for spectrochemical analysis. The samples in acidic solution (HCl) were evaporated on copper electrodes and analyzed as outlined for the copper spark method (1,2). An Applied Research Laboratories spark source was used.

The spectra were taken on a Baird Associates spectrograph of the modified Eagle type which has a dispersion of 5.6 A/mm in the first order.

Spectrum Analysis #1 photographic plates were used and were processed according to the manufacturer's recommendations.

The plates were placed in the projection comparator (also an Applied Research Laboratories Unit) and the lines measured using the scale associated with the projection screen. The copper lines of the spectra were used as a wave length standard.

In all 54 lines which we feel can be attributed to curium were found. These, together with an estimated intensity based on an arbitrary 1 to 10 scale, are listed in Table I.

Due to the high alpha-activity of curium special care had to be exercised in handling the samples. In this case a pyrex bulb with quartz windows attached was made and all the sparking of curium was done in this container. An outlet from this bulb was connected through a glass wool filter system to a water aspirator. Upon completion of the work the bulb was flushed with nitric acid and then discarded. The acid wash and the glass wool were then returned for recovery of the curium.

Spectral Lines of Curium from 3100 A to 4200 A

Table I

Wave Length A	Inten- sity	Wave Length A	Inten- sity	Wave Length A	Inten- sity
3092.45	1	3546.13	5	3846.92	4
3105.83	1	3566.94	5	3860.9*	1
3145.80	1	3570.16	5		1
3196.51	1	3572.57	5	3936.72	2
3210-00	4	3580.79	3	3947.99	3
3220.79	6	3585.14	1	3948.79	1
3230.32	4	3591.41	1	3960.09	3
3238.36	2	3670.4*	1	3974.74	9
3252.69	4	3683.30	4	4003.33	8.
3279.11	2	3692.32	4	4024.42	6
3296.18	4	3734.3*	1	4027.96	1
3304.8*	1	3738. 18	1	4031.67	1
3367.90	1	3739.31	2	4061.6*	1
3374.7*	1	3740.7*	1	4078.95	3
3416.95	4	3749.76	2	4106.7*	1
3426.5*	1	3778.0*	1	4141.7*	1
3472.97	8	3779.6*	1	4208.94	9
3498.58	6	3840.49	3	4219.88	4

^{*}All ware lengths good to 0.05 A except those marked with (*); these are good to 0.5 A.

References

- 1. Fred, M., N. H. Nachtrieb and F. S. Tomkins, J. Opt. Soc. Am., 37, 279 (1947).
- 2. Bachelder, M. C., J. G. Conway, N. H. Nachtrieb and B. S. Wildi, Atomic Energy Commission Declassified Document, MDDC 511, November 6, 1945.



