

Near Real Time Liquid Analysis for Explosive Destructive System

**Victoria VanderNoot, Scott Ferko, Kamlesh
Patel, Jim Van De Vreugde, and Brent Haroldsen**

**Sandia National Laboratories
*Livermore, CA***

***14th International Chemical Weapons Demilitarization Conference, 23-26 May 2011
Interlaken Switzerland***

We are Working to Decrease the Process Time of EDS

One approach is to develop an on-board, near-real-time liquid analyzer (NRT-LA)


- EDS effluent is not drained until analysis confirms that the agent is destroyed
- Manual sample collection, transport and analysis takes an hour or more
- First sample at one hour is almost always non-detect

An NRT-LA will have several benefits

- Allow continuous monitoring of treatment
- Save time spent collecting and analyzing samples
- Reduce operator workload collecting samples
- Improve safety by eliminating handling of samples

Challenges to be addressed

- Functional while EDS vessel is operational/rotating
- Background matrix in EDS may be quite complex
- Detection sensitivity (sufficient to determine treatment threshold)
- Residual HD must be determined, not just hydrolysis or degradation products



Incorporating an On-board Chemical Analysis System into EDS

Mount Capillary electrophoresis-based analysis system on the EDS

Automatically collect and analyze samples at regular intervals (5 to 7 minutes)

- Filter and prepare sample as needed
- Separate using capillary electrophoresis
- Detect peaks using COTS UV detector
- Analyze data to detect and quantify agent peaks

Drain vessel when destruction target is reached

- Effluent remains under engineering control



Collect final sample for laboratory RCRA analysis

- On-line analyzer is for process control only; verification of agent destruction will still require laboratory analysis
- EDS processing can continue while sample is analyzed



Previous Year's Development Status: Demonstrated Proof of Concept with Agent


- Separation methods were developed for mustard and Lewisite class
- Bench top prototype (1st generation)
 - Assembled and tested at ECBC (Dec 09)
 - Demonstrated proof of concept
 - Observed decrease in mustard (HD) peak with corresponding increase of thiodiethanol (TDE) due to hydrolysis




Current Effort is Geared to a Prototype for Testing on the EDS Vessel

Development and testing of a 2nd generation prototype


- Refined the hardware design based on proof of concept testing carried out at ECBC
- Tested with agent at ECBC (late May 2010)
- Determined detection/quantitation limits in more realistic sample matrices
- Demonstrated significant level of integration and automation



Agent Testing with Mustard Shows Good Resolution and Quantitation




*IPA (isopropyl alcohol)



Quantitation Limit of HD in aqueous samples is ~10ppm

Agent Testing with Mustard Shows Sufficient Sensitivity in Neutralent



Quantitation limit (S/N ratio of 5) of HD in 50% MEA is ~25 ppm

• in 90% MEA: 35 ppm



Replicate injections of 25 ppm HD in 50% MEA

Peak Height ($\times 10^{-3}$)	Peak Area ($\times 10^{-3}$)
0.2267	0.6943
0.2807	0.8154
0.2141	0.6003
0.1800	0.3622
$<0.23 \pm 0.04>$	$<0.6 \pm 0.2>$

NRT-LA Method is also Suitable for Lewisite Agent Class

Hydrolysis of Lewisite ($t_{1/2} \sim 1$ minute) is too rapid to detect directly

- Detection is of CVAA (2-chlorovinyl arsenous acid), also a blister agent
- Comparable to current state-of-the-art GC methods that derivatize Lewisite/CVAA and do not discriminate between the two
- Detection of CVAA not impacted by presence of Lewisite neutralent
- Identical hardware & reagents were used




*CVAOA (2-chlorovinyl arsonic acid)

Current Development Status: Prototype Ready for Testing on the EDS Vessel

3rd Generation Prototype

- Mounted on EDS vessel (shakedown testing underway)
- Fluid, data and electrical connections made to EDS vessel


Operational testing on EDS at APGEA scheduled for this summer





3rd Generation Prototype was Designed to Mount onto EDS

Attributes/Features

- Builds on the same analysis platform
- Can collect & analyze multiple samples during first hour of processing time
 - Minimizes manual intervention required to take a sample (PPE)
- Integrated sample collection tube with the door valve manifold
 - Maintaining chemical agent containment and pressure safety
 - Wastes contained & flushed back into treatment vessel between analyses
- Operation During Rotation
 - Orientation independence of analysis platform
 - Orientation sensing to ensure sampling when port is submerged
 - Slip ring electrical connections to supply electrical power
- Automated operation w/ wireless communication interface to computer
- Pressure transducers to monitor system pressure
- Temperature controlled
- Same hardware can be adapted to other agent classes



CAD Illustration of 3rd Generation Prototype Mounted on P2





Additional Potential Applications of a NRT-LA

NRT-LA is a generic analysis tool that is suitable for a wide range of chemical analyses

– Other Chemical agents

- CE is versatile: can tune to any number of agents (nerve agents, explosives)
- pH, viscosity, voltage/polarity, buffer additives
- Capillary length/chip based separations

– Underwater or ocean burial sites

- Current technology is further miniaturizable, suitable for seagoing vessel or submersible operation
 - Remote sampling w/ analysis on nearby vessel/submersible OR analysis done at sample collection site
- Can sample fluids and pre-concentrate agents out of fluid for injection (e.g. sea water)
 - Dilute samples can be concentrated
 - Locate leaking munitions



Acknowledgements



Research funding was funded and directed by the US Army Project Manager for Non-Stockpile Chemical Materiel. We are grateful for the assistance of Kevin Morrissey, Lucille Forrest, Amanda Schenning and Janet Fouse for facilitating agent testing at ECBC.

Team

Victoria VanderNoot, Scott Ferko, Jim He, Jerry Inman, Dan Knight, Ken Patel, Tom Raber, Ron Renzi, Jeanne Stachowiak, Jim Van De Vreugde, Bob Crocker, Joanne Volponi and Brent Haroldsen