




In operando Investigation of SOFC Electrodes Using Synchrotron-based Ambient Pressure X-ray Photoelectron Spectroscopy (AP-XPS) in a Novel Two- Environment Chamber

Anthony McDaniel, Farid El Gabaly, William Chueh, and Kevin McCarty
Sandia National Laboratories, Livermore, CA

Sossina Haile
California Institute of Technology, Pasadena, CA

Michael Grass, Zhi Liu, and Hendrik Bluhm
Lawrence Berkeley National Laboratory, Berkeley, CA











Motivation and approach



- **Perovskite-structured oxides**
 - SOFC cathodes and proton conducting membranes
- **Open questions that impact electrochemical activity, conductivity, and chemical stability**
 - Material stoichiometry
 - A- and B-site substitutions and dopants
 - Phase segregation and substitution-site enrichment
 - Oxygen defects
 - Reactive environment
 - Gas composition, impurities, applied electrical potential
- **In operando observations using soft X-rays**
 - Determine elemental composition of surface and bulk
 - Determine oxidation state of elements on surface and in bulk
 - Spatially resolve surface potential landscape
- **Experimental platform operating at near ambient pressure with XPS, XAS, and ACIS**




AP-XPS and XAS on SOFC-type materials at realistic operating conditions




ABO₃/YSZ/SDC/Pt membrane electrode assembly (MEA)



ALS BL 11.0.2



Nernst potential





- Electrolyte supported MEA
 - High temperature glass seal between YSZ support tube and MEA
 - Electrodes patterned by PLD and metal evaporation through shadow masks
- $P \leq 10$ torr, $T < 800$ °C
- Precision X-Y-Z manipulation using stepping motors for spatially-resolved measurement ($\sim 5\mu\text{m}$)
- FRA/potentiostat used to conduct electrochemical experiments
 - AC-impedance spectroscopy

Oxygen reduction reaction on $\text{Ba}_{0.5}\text{Sr}_{0.5}\text{Co}_{0.8}\text{Fe}_{0.2}\text{O}_3$ perovskite cathode






- Understand why stoichiometry impacts performance
 - Conductivity, electrocatalytic activity, chemical stability
- Understand how system responds to oxygen activity
 - Cation oxidation state
 - Thermochemically driven
 - Electrochemically driven
- Resolve differences between bulk and surface composition
 - Phase segregation
 - Cation valency




Oriented thin films deposited by pulsed laser deposition



Cell Geometry



Anode



Highly reversible thin film SDC/porous Pt anode



Characteristic electron migration length @ OCV $L_C = \sqrt{\sigma_{eon} \tilde{R}^* l} \approx 2,000 \mu\text{m}$

Film Composition: $\text{Ba}_{0.53}\text{Sr}_{0.53}\text{Co}_{0.80}\text{Fe}_{0.20}\text{O}_x$
(determined via EDS calibrated with ICP-OES)




Top view

- Dense geometrically well-defined microelectrodes





Si and S impurities surface segregate



- Higher photon energies probe deeper into bulk
 - Greater escape depth for photoelectrons with higher kinetic energy
- S impurity comes to surface on heating, returns to bulk on cooling
- All impurities (S, Si, Pt) amount to < 35% of surface atomic fraction excluding oxygen






A-site cations (Ba, Sr) behave differently at OCV



- **Ba atomic fraction enriched on surface and exists in a chemical state different than bulk**
 - Sr does not bond differently in the bulk
- **Co and Fe are also depleted in the near surface region**




Oxygen activity consistent with literature reports




- O 1s spectra very complex and involves contributions from surface-oxide impurities
 - Clear trends are evident with O_2 pressure
- Measured ASR is consistent with other PLD microelectrode studies

XAS reveals interesting B-site cation (Co) behavior



- **Partial electron yield (PEY) probes surface states (< 3 nm)**
 - Co oxidation state invariant to O_2 partial pressure at OCV
- **Fluorescence yield (FY) probes bulk states (< 100 nm)**
 - Bulk Co^{+3} concentration decreases with oxygen activity





Conclusions

- **No significant variations in cation chemical behavior at surface observed when MEA is biased**
 - Photoelectron spectra at 300 eV unchanged at 0.6 V bias
 - Still working through complete data set
- **Successfully developed an experimental platform for applying synchrotron-based AP-XPS and XAS to investigating MEA in operando**
- **BSCF cathode**
 - Ba appears to have a different bulk state than surface state at OCV
 - Co has a mixed valency at OCV
 - Surface state remains oxidized as bulk reduces when P_{O_2} decreases
 - Impurities (S, Si) surface segregate
 - Unclear how surface electrochemistry is impacted

Funding and Support



This research was supported by the U. S. Department of Energy under Contract No. DE-AC04-94AL85000 (Sandia) and DE-AC02-05CH11231 (LBNL)

