


WiTec at Sandia:


Pushing a Great Tool Further





Thomas Beechem

Guiding Questions

What do I get?



How do I analyze it?






Could it be better?

?



Our System: Twins




A Typical Week in our Lab


Week of: 11/28-12/2	Monday	Tuesday	Wednesday	Thursday	Friday
Sample	Solvent ZnO Nanowires on Glass	GaAs Transistors	GaAs Transistors	MFO C60 (AM) SiGE-12416F (PM)	Graphene 1ML/H/SiC 110114 #4 and #11
Scan Type	Map/Line	Temp Measurements	Temp Measurements	AM-Work with TEB PM-Linescans	Map 25X25 Micron 75X75 Points
Grating	600 l/mm 550 Center	2400 l/mm 50 Acquisitions Each Power 400 Center	2400 l/mm 50 Acquisitions Each Power 400 Center	AM-600 l/mm 2046 PM-2400 l/mm 550 Center	600 l/mm 2046 Center
Counts	Max 434 Pak	20,000 Counts for 289 Peak	20,000 Counts for 289 Peak	AM-Experiment PM-24,000 on 520 Peak	7,000-8500 Counts for 2D Peak
Objective	20X/0.45 Expt. With 100X/0.9	20X/0.45	20X/0.45	AM-20X/0.45 PM-100X/0.9	100X/0.9
Laser/Fiber	10 mW 100 Micron	10mW 100 micron	10mW 100 micron	AM-10 mW, 100 micron PM-10 mW, 25 micron	20mW 25 Micron

Core Effort: Raman Imaging

Graphene Transfer




Graphene Science




Nanowires and Tubes




Cell Imaging








Phononics





Graphene Device Fabrication




Raman for Thermal Sciences




Raman + Photomapping



How do I analyze? Calibration



Moving Data: Matlab & .spe






Next Steps: Hardware Flexibility


Tunable Spectroscopy

New SPM Modes

High Throughput Couplers




Beam Combiners?



Monochromator Imaging



No Fiber Options




Next Steps: WiTec Control


Remote Control?




Multiple Acquisition: Map



Built-In Raman Calibration



Appearance Control




Next Steps: WiTec Project

Are you only a PC?



Data Bridge



Scripting



Social Users Network



Project Manager

