

Time domain and level-specific investigations of collisions using picosecond two-color resonant four-wave-mixing spectroscopy: OH X²Π collisions in flames

Xiangling Chen and Thomas B. Settersten
Combustion Research Facility
Sandia National Laboratories, Livermore, CA

Funded by
U.S. Department of Energy, Office of Basic
Energy Sciences, Division of Chemical
Sciences, Geosciences, and Biosciences.


Overview

Time-resolved two-color resonant four-wave-mixing spectroscopy (TC-RFWM) provides a powerful means to investigate both inelastic and elastic collisions. This technique is demonstrated using custom-built picosecond laser systems, which produce nearly transform-limited laser pulses. The 50-ps pulses provide adequate time and spectral resolution for rotationally resolved studies of collisions affecting ground-state hydroxyl radicals in an atmospheric-pressure methane-air flame. The combination of double resonance, time-delayed probing, and independent control of the polarization of each of the four fields involved in the wave-mixing process enables rotational-level-specific measurement of the decay of laser-induced population, alignment, and orientation, as well as state-to-state transfer of these three moments.


Picosecond TC-RFWM for collision studies

- Picosecond pump-probe technique
 - direct measurement of collision rates in time domain
 - requires sufficiently **short** laser pulses
 - requires **synchronization** of pump and probe lasers
- Two-color (double resonance)
 - enables rotational-level-specific measurements
 - requires sufficiently **narrowband** lasers
- Coherent spectroscopy ($\chi^{(3)}$ process)
 - sensitive to alignment/orientation
 - immune to background in flames
 - requires sufficiently **high-power** laser pulses

Custom picosecond laser system enables time-domain studies in atmospheric-pressure flames



Experimental parameters dictate relevant wave-mixing physics



Polarization dependence of time-domain TC-RFWM

The contribution to the signal electric field with polarization \vec{e}_4 resulting from each relevant Feynman diagram can be expressed as follows using diagrammatic perturbation theory.

$$[\mathbf{P}^{(3)} \cdot \vec{e}_4] = C \Phi \sum_K F(\vec{e}_1 \vec{e}_2 \vec{e}_3 \vec{e}_4; K) G(J'' J' J^*; K)$$

Spherical tensor analysis is used to express the problem in terms of the three ranks, K , of the state multipoles.

- $K=0$ Population
- $K=1$ Orientation
- $K=2$ Alignment

The constant C contains the signal dependence on **molecular number density, pump and probe field amplitudes, and transition line strengths.**

$$C = \frac{N_i \mathcal{E}_1 \mathcal{E}_2 \mathcal{E}_3}{8\hbar^3} \frac{|\langle J^* || \mu^{(1)} || J' \rangle|^2}{2J' + 1} \frac{|\langle J'' || \mu^{(1)} || J'' \rangle|^2}{2J'' + 1}$$

The phase term for the first two Feynman diagrams on the left is given below. The **second exponential term** describes dephasing resulting from motion of molecules forming the intermediate-state grating $\vec{k}_G = \vec{k}_1 - \vec{k}_2$.

$$\Phi = e^{-i(\omega_4 t - \vec{k}_4 \cdot \vec{r})} e^{-i\vec{v} \cdot [\vec{k}_4(t_4 - t_3) + \vec{k}_G(t_3 - t_2)]} e^{-\Gamma_{J^*}(t_4 - t_3)}$$

The signal dependence on electric-field polarization is contained in F . **Control of polarization** can be used to ensure signal is dependent only on population or orientation or alignment by forcing $F \neq 0$ for only one K .

$$F(\vec{e}_1 \vec{e}_2 \vec{e}_3 \vec{e}_4; K) = \sum_{Q_1 Q_2 Q_3} (-1)^{K-Q_2} (2K+1) (\vec{e}_1)_{-Q_1}^{(1)} (\vec{e}_2)_{Q_1-Q_2}^{(1)} (\vec{e}_3)_{Q_2-Q_3}^{(1)} (\vec{e}_4)_{Q_3}^{(1)} \times \begin{pmatrix} 1 & 1 & K \\ Q_1 & Q_2-Q_1 & -Q_2 \end{pmatrix} \begin{pmatrix} 1 & 1 & K \\ Q_3 & Q_2-Q_3 & -Q_2 \end{pmatrix}$$

G contains dependencies on the spectroscopic transition branches and the decay rates $\Gamma_K^{J'}$ for the K^{th} moment of the intermediate level J' .

$$G(J'' J' J^*; K) = (2J'+1) \left\{ \begin{matrix} 1 & 1 & K \\ J' & J' & J'' \end{matrix} \right\} \left\{ \begin{matrix} 1 & 1 & K \\ J' & J' & J^* \end{matrix} \right\} e^{-\Gamma_K^{J'}(t_3 - t_2)}$$

Polarization/excitation schemes


$|F(\vec{e}_1 \vec{e}_2 \vec{e}_3 \vec{e}_4; K) G(J'' J' J^*; K)| e^{+\Gamma_K^{J'}(t_3 - t_2)}$

	$K=0$	$K=1$	$K=2$
Population			
CC ↑ ↑ PP	1/9.07	0	1/177
CC ↑ ↑ PQ	1/9.07	0	1/177 J'+1
↗ ↗ ↑ ↑ PQ	1/9.14	0	1/91.4 J'+1
Orientation			
CC ↑ - PP	0	1/12.05	1/39,400
CC ↑ - PQ	0	1/12.05 J'+1	1/39,400 J'+1
Alignment			
↗ ↗ ↑ - PQ	0	1/7,890 J'+1	1/60.2 J'+1

Pump "magic angle" ~ 54.9° (phase match: $\beta=5^\circ$, $\gamma=1^\circ$)


Two-beam, polarization spectroscopy configuration

- Simplest set-up
- Automatically phase matched
- Detect component of signal with polarization perpendicular to probe polarization
 - $\vec{e}_4 \perp \vec{e}_3$
 - Susceptible to background from probe
 - Collision rate measurements are always influenced by alignment and/or orientation




Four-beam, grating spectroscopy configuration


- Phase matching is critical
- Complete experimental control over polarization of all three waves driving $\chi^{(3)}$
- Polarization-specific detection
 - population
 - orientation
 - alignment




TC-RFWM/PS is immune to thermal background



Δ-doublet resolution and sub-ns time resolution



Direct measurement of total decay rates in flame



Level-to-level transfer of population/anisotropy

