

2007 IEEE Conference on Technologies for Homeland Security
May 16-17, Woburn, MA

A Review of Microfabricated Preconcentrators for Portable Chemical Analysis Systems

**Ronald P. Manginell, Matthew Moorman, D. Adkins[†], P.
Lewis, D. Porter, R. Hadizadeh, D. Copic, R. J. Simonson, R.
Kottenstette, J. Whiting, S. Radhakrishnan*, R. Cernosek**

Sandia National Labs, Albuquerque, NM

rpmangi@sandia.gov

[†]Defiant Technologies, Inc., Albuquerque, NM





www.defiant-tech.com

***ESI-Group, Hunstville, Alabama**

Sandia is a multiprogram laboratory operated by Sandia Corporation, a Lockheed Martin Company, for the United States Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.

Sandia's MicroChemLab™


A hand-held chemical analysis system that uses three integrated modular components





Preconcentrator accumulates analytes of interest

Gas Chromatograph separates analytes in time

Acoustic Sensors provide sensitive detection




1997






SAW Array

US Patents: 6,171,378, 6,527,835
IEEE Sensors Journal, 6 (3) 784-795, 2006.



Planar MicroFabricated Preconcentrators

- **Low C, high efficiency adsorbent platform**
 - 2000°C/W; 10msec ramp
- **Minimal flow restriction**
 - 5 psig, 200 mL/min
- **Concentrate targets**
- **Reject interferants**
- **Rapid release - a non-mechanical GC injector**
- **Bosch or KOH etched to SiN**




CFD Modeling and DOE

DOE


- Carboxen 1000 and light hydrocarbon
- Varied [C], collection time (t_c), desorption flow (f), temperature (T) and desorption time (t_d)
- GC/FID Agilent 6890
- Statistica - full quadratic

Conclusions

- Peak area, A_p , increases with [C]
- **Maximum in A_p with f**
- Peak width, W , is not influenced by [C]
- W decreases with f
- Max in W with T
 - heated area increase, degradation
- Increased tailing with T
- $Pe \sim$ convection/diffusion increases with f , [C] and decreases with T




R. J. Simonson, et. al., "Optimization of a Microfabricated Planar Preconcentrator,"
Proceedings of the 2nd Joint Conference on Point Detection for Chemical and
Biological Defense, Williamsburg, VA 3/1-5/2004, Manuscript K1.




CFD Modeling (ESI Group) & DOE





- Computational Fluid Dynamics
- 2D – flow and adsorbent scaled
- Simplified thermal model
- Unity sticking
- Calibration on DOE
- Adsorption: $A + s \rightarrow A(s)$
 - $k = 36,500 \text{ s}^{-1}$; $25,300 \text{ s}^{-1}$ from Modified-Wheeler
- Desorption: $A(s) \rightarrow A + s$
 - first-order Arrhenius 30.1 kJ/mol
- Can predict other DOE runs
- Aris-Taylor Diffusion
- Diffusion is a dominant effect
- Did not predict fall off in A_p with f
 - Quadratic or simplicity of model; turbulence not an issue




R.P. Manginell, Sekhar Radharishnan, et. al., "Two-dimensional modeling & simulation of mass transport in microfabricated preconcentrators", accepted IEEE Sensors Journal.



Eye Candy: Thermal pulse can be seen on a downstream detector






3DPCs as a supplement or replacement for the planar PC


Planar PC

1. Low C, high thermal efficiency
2. Fast response, low power
3. Collection limitations




3DPCs

1. Planar PC items 1 & 2 retained
2. Smaller diffusion length, higher area, flow through
3. Pressure balance possible







Improved
collection
performance


More analytes like
volatile organic
compounds (VOC)









Types of 3DPCs: perpendicular flow. Etching thanks to ITC.







3DPC: parallel flow and tortuous path




3DPCs, coatings and target analytes: enhanced collection


- **Spray and drop coating of sols**
 - CW agents, TICs
 - Explosives – usually need tortuous 3DPC
 - Automated spray with tilt
- **Commercial packing in PDMS binder OR using packing stops**
 - Toxic Industrial Chemicals (TICs) and Tri-Halomethanes (THMs)
 - PoropakQ, HayesepA, Carboxen
- **Laser ablation of nanoporous carbon**
 - Conformal coating; TICs




“Perpendicular flow”



“Parallel flow”





US 7,118,712 Non-Planar Chemical Preconcentrator





Comparison of Collectors


Preconcentrator Device Comparison with TIC Analytes



0.6 W vs 3 W







Design	t _{90 to 200C}	Power
Planar	10 msec	100 mW
3D	0.6-1 sec	200-600 mW
Tubular	~minutes	~watts



Benefits of 3D Design

- Improved collection/desorption – higher surface area, better contact, lower dead volume
- Lower level detection, new analytes
- Ease of assembly
- Water, and VOCs can now be addressed
- TIC, THM, CW
- Explosives





DARPA MGA: Performance requirements drive system architecture enhancements

Maintain: Low false alarm rate, $< 1/200,000$

Increase: Analysis speed, Analytical channel capacity


Decrease:

System volume < 20 cc

Limit of detection





Energy consumed per analysis

$< 20 \text{ cm}^3$






DARPA







Smart PC™ combines preconcentration and detection to accelerate and automate detection

- DOD: reduced analysis times with increasing target concentration
- MEMS resonator with a heater/adsorbent weighs the sample & decides when it has collected enough
- Modular fixtures
- Circuit autotunes, autozeros
- Software subtract reference and smooth



SMART PC™



**Reduced need for
trained operators**



USPTO 7,168,298 Mass Sensitive Preconcentrator

Detection with selectivity: 7 times faster at LC50 of Sarin

SPC - GC - SAW 12/16/05
Vapor System 1ppm DMMP
Cooked DKAP on SPC, DKAP on SAW



Next time



Conclusions

- **Microfabrication can be used to make a variety of preconcentrators**
- **Planar devices can be very fast, non-mechanical GC injectors**
- **3D designs have higher capacity and can allow a broader range of analytes and/or lower concentration**
- **3D designs with integrated passive valves are in testing for DARPA**
- **Smart PC allows for fast analysis when concentrations are the highest**