

Femtosecond Pure-Rotational CARS for Gas-Phase Diagnostics

SAND2011-5657C


Sean P. Kearney and Darcie Farrow
Sandia National Laboratories
Albuquerque, New Mexico 87185 USA

Joseph D. Miller and Terrence R. Meyer
Iowa State University
Ames, IA 50011 USA

1. Background


"Typical" Nanosecond CARS

- Measurements performed in the *frequency domain*
- The nanosecond approach has its limitations
- Nonresonant background
- Linewidth/collisional modeling
- Data rates are limited to around 10 Hz
- Limited single-shot precision: broadband Stokes dye laser statistics



Femtosecond CARS (fs-CARS)

- Nonresonant background suppression
- Near collision-free measurements
- Transform-limited fs pulses exhibit low amplitude and phase noise
 - Superior precision
 - Efficient signal generation
- Increased single-shot data rates (1-10 kHz)
- Paradigm shift for gas-phase community: *time-domain vs. frequency domain*




2. Vibrational fs-CARS: Time-Domain and Single-Shot Measurements

- Transient Q-branch ($\Delta J = +1$, $\Delta J = 0$) Raman polarizations reflect *frequency spread dephasing*
- Direct *time-domain* measurements



$$P_{res} \sim \sum_i \Delta N_i \left(\frac{d\sigma}{d\Omega} \right)_i \cos(\omega_i t) \exp(-\Gamma_i t)$$

- Single-shot measurements require modification to the fs-CARS approach





Lucht et al., *Appl. Phys. Lett.* **89**, 251112 (2006)

Chirped Probe Pulse (CPP)



Roy et al., *Opt. Lett.* **34**, 3857 (2009)

Transform-Limited Probe Pulse (TLPP)



Miller et al., *Opt. Lett.* **35**, 2430 (2010)


3. Can Pure-Rotational CARS be Applied in Similar Fashion?

- More precise at "low" (<1000 K) temperatures
- Rotational transitions are closely spaced – more species with a single Raman instrument
- Rotational Raman spectra have wider bandwidth (200 cm⁻¹) than vibrational spectra (20 cm⁻¹)
 - More rapid initial dephasing of the Raman coherence
 - Overall decay is collision instead of dephasing dominated
 - Strong signals at long probe delay – nonresonant background suppression
- "Rephasing" results from ~uniform energy 4B spacing in $\omega_{JJ'}$


$$P_{res}(t) = \sum_k \sum_{J, J'} W_{JJ'}^{(k)} \exp(-i \omega_{JJ'}^{(k)} t) \exp(-\Gamma_{JJ'}^{(k)} t)$$

$$W_{JJ'}^{(k)} = X_k \gamma_k (N_J - N_{J'}) F_{JJ'} b_{JJ'}$$


Time Domain Air at 294 K with Scanned fs Probe



Time-Domain N₂ at 294 K with Scanned fs Probe




Frequency Domain N₂ Spectrum (ns-CARS)




4. Temperature Sensitivity of Rotational Transients


Experiment: Air in Tube Furnace



Calculation N₂ at Early Time




Calculation N₂ Late Time







- Structure of the observed transients is not particularly sensitive to temperature
- Some increased sensitivity is observed at late times and arises from higher order corrections to the Raman frequencies

5. Fourier Transform Spectra

- Better temperature sensitivity observed in the frequency domain.
- Fourier spectra of the observed transients show Stoke/Stokes, anti-Stokes/anti-Stokes and Stokes anti-Stokes quantum beats
- Distinct migration of Fourier spectral energy to lower frequencies at temperature increases
- Next step: direct frequency domain measurements for single-shot data





6. Initial CPP Rotational CARS Results



- Mechanically scanned probe data do not permit practical temperature measurement on a single laser shot
- Probe pulse is chirped by dispersive SF-10 rod
- CPP maps time-domain information into the frequency domain
- Rotational analog of CPP Vibrational CARS technique of Roy et al.
- Results are extremely sensitive to laser field phase and amplitude

7. Early TLPP Measurements in Cold N₂

4-f Pulse Shaper



- 4-f pulse shaper slit removes probe bandwidth and stretches pulse to ~10 ps
- Available probe energy is reduced dramatically relative to CPP approach
- Raman spectra are well resolved and much simpler to interpret
- With 20-ps probe delay, these spectra are free of nonresonant background
- Single-shot spectra in N₂ display very low noise, a key advantage of fs pump/Stokes preparation

8. Summary and Conclusions

- Rotational CARS/CSRS is being developed for single-shot gas-phase measurements. We present some early results here
- The observed transients in N₂ and air are in good agreement with previous observations and are well-modeled
- Time domain transients are less sensitive to temperature than in the frequency domain
- Both CPP and TLPP are being investigated for single-shot measurements in the frequency domain
- CPP permits high probe-laser energies, but modeling of results is very sensitive to laser electric fields
- TLPP spectra are more readily modeled, but probe laser energy may be limited