

A Capture-Based Technique for Enhancing RNA-Seq Analyses of Bacterial Transcriptomes During *In Vitro* and *In Vivo* Infections

Zachary W. Bent
Systems Biology
Sandia National Labs


RNA-Seq

- Sequencing of all RNA transcripts in a sample
 - Coding and non-coding
- Create snap shots of gene expression
 - transcriptomics
- Find regulatory RNAs
 - lncRNA
 - miRNA
- Analyze splice variants
- Large scale SNP analysis

Host-Pathogen Interactions using RNA-Seq



Total RNA
Extraction



Mix of **Host** and **Pathogen** RNA

Prepare RNA-Seq
Libraries


HiSeq 2500



Host-Pathogen
Transcriptomics

RNA transcripts of MΦ infected at an MOI of 10

99.75% Host




- 0.25% of transcripts are bacterial
- 3% of bacterial reads are CDS
- 0.0075% of total reads are bacterial CDS
- 12000 reads per HiSeq lane
- ~83 HiSeq lanes = 1,000,000 reads
- ~\$125,000 per sample


Non-specific methods to enrich bacterial transcripts

- Pre-filtering of bacteria
 - 5 μ M filter
 - Density centrifugation
- Ribo-Zero Epidemiology (Epicentre)
 - Removes host and bacterial rRNA
- DSN or HAC treatment
 - Removes most abundant transcripts
 - Typically host rRNA


Capture-based bacterial transcript enrichment and rRNA depletion




In vitro Y. enterocolitica infections



Growth Media
(LB w/o NaCl)
26° C



Conditioned RPMI
MΦ grown for 2 days
Filter sterilized
5% CO₂
37° C
30min, 1hr, 2hr, 4hr



Infection of MΦ
MOI = 10
37° C
30min, 1hr, 2hr, 4hr

Gentamicin Treated
MOI = 10
1hr of infection
1hr of gentamicin

Transcriptome of
internalized bacteria


37° C and
media
induction

Infection-dependent
genes

Infection-dependent differential gene expression


Gene ID	Name	Function	Fold Change
YE2639	-	Sugar Transporter	107.65
YE4027	-	Epi-inositol hydrolase	63.44
YE2638	-	Sugar Transport	54.79
YE4025	-	Inosose dehydratase	46.21
YE4026	<i>idh</i>	Myo-inositol 2-dehydrogenase	32.73
YE0267	<i>fadA</i>	Fatty acid oxidation	31.52
YE0447	-	Unknown membrane transporter	26.63
YE3198	<i>aglB</i>	6-phospho- α -glucosidase	21.75
YE0309	<i>acs</i>	acetyl-CoA synthetase	20.38
YE4031	-	aldehyde dehydrogenase	18.21

Gene ID	Name	Function	Fold Change
YE0553	<i>scrY</i>	sucrose porin	-24.15
YE3092	-	Unknown	-18.47
YE3297	<i>xni</i>	exonuclease	-15.63
YE0789	-	thiol:disulfide interchange	-14.42
YE2667	-	Unknown	-14.32
YE1740	-	Unknown	-12.59
YE2772	<i>hisl</i>	ATP pyrophosphatase	-12.19
YE0552	<i>scrA</i>	Sucrose transporter	-11.30
YE0195	<i>yigB</i>	Flavin mononucleotide phosphatase	-9.69
YE0781	<i>mrfC</i>	Fimbrial protein	8.12




Infection-dependent differential expression by functional category

(Genes differentially expressed between RPMI and infection)



Expression of the Ysa T3SS



	<i>yspA</i>	<i>D</i>	<i>C</i>	<i>B</i>	<i>scyB</i>	<i>U</i>	<i>T</i>	<i>V</i>	<i>E</i>	<i>6-7</i>	<i>H</i>	<i>1</i>	<i>2</i>	<i>ysrR</i>	<i>S</i>
Control FPKM	5.2	13.6	3.1	1.5	0.8	1.3	1.6	0.5	1.0	2.5	1.7	14.9	15.1	19.4	6.7
RPMI FPKM	41.1	78.4	62.3	54.9	147.2	0.9	0.5	1.0	0	1.3	1.2	16.4	12.4	10.8	2.9
Infection FPKM	97.1	162.2	125.1	119.0	335.6	21.7	7.2	4.3	2.5	16.3	21.7	220.2	144.5	52.1	15.2


In vivo *S. Typhimurium* infections



Livers $2.5 \times 10^2 - 5.5 \times 10^3$ CFU/sample

Spleens $1.3 \times 10^4 - 5.3 \times 10^5$ CFU/sample

In vivo enrichment



- Synergistic effect of combining Ribo-Zero with capture ~2X
- Despite use of Ribo-Zero, significant amount of bacterial rRNA is present

Conclusions

- RNA from infected samples is dominated by host transcripts
- Capture is a highly effective and cost efficient tool to enrich bacterial transcripts from infected samples
- In combination with Ribo-Zero, capture can achieve over 10,000X enrichment
- Capture technique is highly versatile with potential applications in variety of fields
- Inositol metabolism may play an unrecognized role in Yersinia pathogenesis
- The Ysa T3SS system is expressed during the infection of macrophage cells

Acknowledgements

Sandia National Labs - CA

Robert Meagher
Stan Langevin
Kunal Poorey
Anupama Sinha
Deanna Curtis
Steve Branda
Victoria VanderNoot
Mary Tran-Gyamfi



**Sandia
National
Laboratories**

University of Colorado

Corrie Detweiler
Eugenia Silva-Herzog
Dave Brazel
Angelika Krivenko



Sequencing

Vincent J. Coates Genomic Sequencing Laboratory
(UC Berkeley)
Minyong Chung

UC Davis

Glenn Young
Karen LeGrand
Rachelle Hamblin



Funding

Sandia National Labs LDRD Office



Contact:
zbent@sandia.gov