


Vibrational Sum Frequency Generation Spectroscopy Microscope

SAND2016-6836PE

Emmanuel Valenton, B.S. Chemistry, UC Berkeley, May 2016
Christopher Klierer, Org. 8353, Combustion Chemistry Dep.

- Background
- Design
- Current Progress



*Exceptional
service
in the
national
interest*




Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000. SAND NO. 2011-XXXXP

"This work was supported in part by the U.S. Department of Energy, Office of Science, Office of Workforce Development for Teachers and Scientists (WDTS) under the Science Undergraduate Laboratory Internship (SULI) program."

Motivation


- Many chemical and physical phenomena depend on surfaces



- Surface-specific signal is weak compared to bulk-specific signal

- Fosca Di Gabriele et. al. *Nucl. Eng. Des.* **2014**, 280, 69-75
- Korenobu Matsuzaki et. al. *J. Am. Chem. Soc.* **2016**, 138, 7551-7557

Sum Frequency Generation



$$\mathbf{P} = \epsilon_0 \chi^{(1)} \mathbf{E}$$

$$\mathbf{P} = \epsilon_0 (\chi^{(1)} \mathbf{E} + \chi^{(2)} \mathbf{E}^2 + \chi^{(3)} \mathbf{E}^3 + \dots)$$

$$\mathbf{E} = \mathbf{E}_1 \cos \omega_1 t + \mathbf{E}_2 \cos \omega_2 t$$

$$\mathbf{P}^{(2)} = \epsilon_0 \chi^{(2)} (\mathbf{E}_1 \cos \omega_1 t + \mathbf{E}_2 \cos \omega_2 t)^2$$


rearranging as before

$$\mathbf{E}_1^2 + \mathbf{E}_2^2$$

$$\mathbf{E}_1^2 \cos 2\omega_1 t + \mathbf{E}_2^2 \cos 2\omega_2 t$$

$$\frac{1}{2} \mathbf{E}_1 \mathbf{E}_2 \cos(\omega_1 - \omega_2) t$$



$$\frac{1}{2} \mathbf{E}_1 \mathbf{E}_2 \cos(\omega_1 + \omega_2) t$$



- Lambert, A. G., Davies, P. B., Neivandt, D. J. *Appl. Spectrosc. Rev.* **2005**, 40, 103-145
- Shen, Y. R. *J. Phys. Chem. C* **2012**, 116, 15505-15509

Sum Frequency Generation


- Selection rules allow sum frequency generation (SFG) if it takes place in an asymmetric medium



- Sum frequency generation is a non-linear spectroscopic method that is surface-specific
- Marialore Sulpizi. *J. Phys. Chem. Lett.* **2013**, 4, 83-87

Experimental Setup

- Goal: to construct an SFG microscope
- Beamline Setup:




Experimental Setup

- Olympus optical microscope
- Thorlabs xy-stage
- Reflective Optic
- (500 – 20000 nm)



Current Progress




Conclusion

- SFG spectroscopy is a non-linear surface specific method
- Sum frequency generation spectroscopy microscope from optical microscope, xy-stage, and wide wavelength range reflective optic
- LabVIEW control software has been created

Future Work:

- Installing microscope into beam line
- Potential topics: energetic materials degradation and soot formation.



Acknowledgements