

New Approaches to Relating the Macroscopic Behavior of Shales to Pore-Scale Constitutive Properties

Jason Heath¹, Thomas Dewers¹, Charles Bryan¹, Mei Ding², Rex Hjelm²


¹Sandia National Laboratories, ²Los Alamos National Laboratory

GSA Annual Meeting, September 25, 2016


Motivation 1 – Shale Pores

***Characterization of multiscale porosity alone is difficult.
Including coupled processes is even harder...***

Example of large-area porosity characterization for context



Backscatter SEM



Shale Pores

mSEM: 2D imaging; pores visible but organics not-so-much; difficult sample prep.; static sample


Helpful but not enough for in situ coupled processes

mSEM



(Heath, Dewers, Milliken, Mozley, Yoon, *in prep*)

Motivation 2 – Geologic CO₂ Storage



Will initially dry CO₂ dehydrate caprock and exacerbate leakage pathways?

Recent studies investigate intercalation of CO₂ and H₂O in clay and shrink-swell behavior under hydrostatic stress (e.g., Loring et al., 2014)

How does lithostatic stress affect shrink-swell of clay?


Oedometric SANS

Key advances:


- Data collection over 1 to 1000 nm (or to 10s of microns if used with USANS)
- Accommodates pore fluids at high pressure and temperature
- Non-hydrostatic stress state applied to sample

Difficulties and opportunities:


- We are developing oedometer v4.0...
- Sample preparation for *in situ* fluids, pressure, and temperature at beam facilities (LANSCE; NIST with USANS); thin sample
- Data interpretation – you don't "see" pores




Oedometer v1.0



SWy-2




- MAWP: 27.6 Mpa
- Made of Ti
- Designed for geomechanics
- Drawback: 1-inch sample chamber
- Multiple scattering




High Pressure Piping System

Oedometer v2.0




- MAWP: 6.89 MPa
- Al window; steel
- Designed for SANS neutron optics
- Drawback: “penny-shaped” crack in metal...



LQD OEDOMETER ASSEMBLY



ASSEMBLED CONFIGURATION



Oedometer v3.0



- A - Elastic Regime
- B - Initiation of yielding
- C - Near dilational “turn-around”
- D - Approaching sample failure

Oedometric-SANS performed on Mancos Shale


- New oedometer design with sapphire windows acting as pistons for uniaxial loading
- MAWP: 60 MPa axial and pore pressure
- Stress parallel to beam
- Used successfully on SANS and USANS instruments (at NIST)
- Fluid injection for contrast measurements

Example from Oedometer v2.0

What is the change in pore structure of montmorillonite (SWy-2) as a function of non-hydrostatic stress conditions and dissolved water in CO₂?

Approach:

- Measure compaction or swelling with oedometer, coupled to SANS
- Take the same clay sample through a stress path with different pore fluids (dry and wet CO₂) and measure pore structure with SANS




H₂O?

CO₂?


Initially
dry supercritical
CO₂ in pores

CO₂ density and H₂O solubility





Concept: Deliver water into sample using CO₂

Timeline of oedometer stress path





Oedometer stress path and stages






$$\text{Pressure} \times \frac{\text{Area}_{\text{piston}}}{\text{Area}_{\text{ram}}}$$





Unified & Fractal Fitting




Power-law relationship:

$$I(Q) = \frac{A}{Q^n} + B$$


Stage	Avg. pore size (nm)	n	D
I: 1 st comp.	37.9	3.64	2.36
II: 2 nd comp.	37.5	3.66	2.34
III: HP, dry CO ₂	34.4	2.83	3.17
IV: LP, dry CO ₂	37.6	3.60	2.40
V: HP, wet CO ₂	39.5	3.14	2.86
VI: LP, wet CO ₂	34.6	3.71	2.29

- Recoverable changes in fractal dimension relate to elastic compression
- Irreversible changes in D relate to plastic and/or creep strains
- Intercalation not obvious


I




II




III




IV




V



VI





- Effects of consolidation: initial compaction shifts to slightly smaller pores
- Release of pore pressure for dry CO₂ seems reversible: falls back near initial compaction curve – elastic strain
- Release to lower pressure of wet CO₂ shows a change in pore sizes (irreversible strain) – there may be multiphase effects that we have yet to account for

Acknowledgements

- Many thanks to the U.S. Department of Energy, Office of Science, Basic Energy Sciences under Award Number DE-SC0006883 for funding.
- This work has benefited from the use of Low-Q Diffractometer at LANSCE at the Lujan Center at Los Alamos Neutron Science Center, funded by DOE Office of Basic Energy Sciences. Los Alamos National Laboratory is operated by Los Alamos National Security LLC under DOE Contract DE-AC52-06NA25396.
- For contributions to the design of oedometers and pressure systems
 - Steve Bauer, Jiann Su, Greg Flint, SNL
 - Mark Taylor, LANL

Transmission and multiple scattering

