


Exceptional service in the national interest



Diagnostics and Development Efforts on Z


Gregory A. Rochau
Manager – Radiation and Fusion Experiments

National Diagnostics Workshop – December 5-6, 2017




Sandia National Laboratories is a multi-mission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC., a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

Z Supports a broad portfolio of experiment campaigns with diverse requirements



Dynamic Materials Properties (DMP)



Inertial Confinement Fusion (ICF)



Secondary Assessment Tech (SAT)




Radiation Effects Sciences (RES)




Fundamental Science (FS)

The CY18 Z shot schedule emphasizes Science Campaigns with ~30% of the shot days in ICF


2018 Shot Days by Research Area



2018 Shot Days in Science Campaigns




2018 Shot Days in ICF



179 planned shots out of 234 shooting days

The US is studying a form of magnetic direct drive called Magnetized Liner Inertial Fusion (MagLIF).


Pre-Magnetize




Beryllium liner

- $B_z = 10-30$ T
- Inhibit e^- conduction
- Confine α 's

Preheat



Compress



- Laser Energy = 1-4 kJ
- $T_0 \sim 100$'s eV
- Reduce required implosion velocity

- CR ~ 35
- $\rho R \sim 0.003$ g/cm²
- P ~ 5 Gbar
- BR ~ 0.5 MG-cm

The MagLIF stagnation column is much different in size, pressure, and time-scale than a typical ICF hot spot.

Laser Indirect Drive On NIF

1 mm



E_{driver}	~ 1.8 MJ
P_{HS}	~ 300 Gbar
R_{HS}	~ 30 um
Y_{DT}	~ 2E16
Δt	< 100 ps

Laser Direct Drive On OMEGA



E_{driver}	~ 26 kJ
P_{HS}	~ 55 Gbar
R_{HS}	~ 20 um
Y_{DT}	~ 5E13
Δt	< 100 ps

MagLIF On Z





E_{driver}	~ 2 MJ
P_{HS}	~ 1 Gbar
R_{HS}	~ 100 um
L_{HS}	~ 10 mm
Y_{DD}	~ 4E12
Δt	> 1 ns

Multiple diagnostics are used to interpret the Magnetic Direct Drive implosion dynamics and resulting stagnation conditions



Optical

Point VISAR




X-rays

Imaging & Radiography




Neutrons


nTOF (side & bottom)




In & Cu Activation



Power History




Spectroscopy




Spherical-crystal diagnostics are a workhorse capability for MagLIF, but all are presently time-integrated

2-D Backlighting and Self-Emission




1-D Space Resolved Spectra



Spherical-crystal diagnostics are a workhorse capability for MagLIF, but all are presently time-integrated

2-D Backlighting and Self-Emission


- Liner morphology (2-frames/shot)
- Hot column morphology



- How does this morphology evolve?

1-D Space Resolved Spectra


- Plasma Conditions
- Mix



- When does the Fe come?
- How does the window mix propagate?
- What is the radial mix distribution?

Fast-gated single line-of-sight hCMOS sensors are a key part of the strategy to time-resolve the MagLIF stagnation

hybrid CMOS Focal Plane Array




Pixelated 25 μm thick Si diode layer on top of a burst-mode CMOS circuit



Daedalus Specifications (Late FY18)

- 512x1024 25 μm pixels
- 12.8 x 25.6 mm active area
- 1 side abutable (12.8 mm x 51.2 mm)
- 3 frames per pixel
- Interlacing for more frames
- 1 ns per frame
- Sense 1-500 photons/pixel at 10 keV

6.2 keV Radiograph at Mag = 6




>6 keV Spectrum at Mag = 1/3



Fast-gated single line-of-sight hCMOS sensors are a key part of the strategy to time-resolve the MagLIF stagnation



hybrid CMOS Focal Plane Array



Pixelated 25 μm thick Si diode layer on top of a burst-mode CMOS circuit

Daedalus Specifications (Late FY18)


- 512x1024 25 μm pixels
- 12.8 x 25.6 mm active area
- 1 side abutable (12.8 mm x 51.2 mm)
- 3 frames per pixel
- Interlacing for more frames
- 1 ns per frame
- Sense 1-500 photons/pixel at 10 keV




hCMOS cameras are now taking data on all three US ICF facilities: Z, NIF, and Omega

Z

Opacity Spectrometer




Axial Pinhole Imager




NIF

**Gated Crystal Backlighting
(Pulse Dilation SLOS)**




Gated LEH Pinhole Imager




OMEGA

**Gated Pinhole Imager
(Pulse Dilation SLOS)**




+30 ps





+90 ps

A new capability for time-resolved opacity measurements gave first data on Z in September.




Hybrid CMOS (hCMOS)
X-ray camera provides
4 time-gated frames of
absorption spectra





Z-pinch
X-ray
Source

Mg/Fe
opacity sample


A new capability for time-resolved opacity measurements gave first data on Z in September.



Hybrid CMOS (hCMOS)
X-ray camera provides
4 time-gated frames of
absorption spectra



Magnesium K-shell spectra are observed



Z-pinch
X-ray
Source

Mg/Fe
opacity sample

hCMOS is entering a production phase, but supply chain issues have caused some delay in delivery for key applications

We need dozens of sensors

Z

Opacity Spectrometer
Axial Imager
Gated Crystal Backlighting
Gated Crystal Spectra
Gated Crystal Imaging
Blast Wave Shadowgraphy

NIF

SLOS1-CBI
SLOS2-KBO
Gated LEH
TARDIS-t
DMX Imager (LMJ-CEA)

OMEGA

SLOS-TRIXI
SLOS-TRIXI-2

- We have fabricated enough ROIC and diode wafers to produce hundreds of sensors.




- Hybridization has become the key risk





- Preferred method is wafer-to-wafer direct bond interconnect (DBI)
- Previous DBI supplier quit supporting this size wafer in April causing a delay in delivery
- New supplier under evaluation.
- Die-to-die indium bump bonding has been demonstrated at Sandia this summer and expected to provide ~10 sensors this year

Die-To-Die (D2D)




We are continuing to develop new hCMOS sensor technology with an emphasis on more, faster frames


Timeline	FY15	FY16	FY17	FY18	FY19	FY20	FY21
	Hippogriff	Icarus	Icarus-2	Daedalus	Daedalus HE		Horus
Furi / Hippogriff		<ul style="list-style-type: none">▪ 1st full-scale multi-frame sensors▪ 1.5-2ns minimum shutter▪ Optimized for 1-10 keV x-ray detection					
Icarus		<ul style="list-style-type: none">▪ 1st cameras compatible with pulse-dilation▪ 1 ns min. shutter with 2 or 4 (Icarus-2) frames per hemisphere▪ Optimized for soft x-ray, visible, and e- detection					
Daedalus		<ul style="list-style-type: none">▪ 3 frames per hemisphere (≥ 6 frames with interlacing)▪ 1-side abutment for spectroscopy and z-pinch imaging applications▪ Large well for high energy x-rays while maintaining low end sensitivity					
Horus		<ul style="list-style-type: none">▪ 1st 180 nm sensor in Sandia's new CMOS-8 architecture <p>Design goals include</p> <ul style="list-style-type: none">▪ 6 frames per pixel w/ independent quadrants (≥ 12 frames with interlacing)▪ 0.5 ns minimum shutter					
							UXI Sensors

The tritium era on Z has begun, with contained (and uncontained!) capability tested at 0.1% DT

First Z shot with 0.1% DT (August, 2016)



Z Gas Transfer System



Neutron time of flight diagnostics measured a primary DT neutron signal for the first time on Z


Vacuum chamber post-shot characterization found tritium contamination levels were indistinguishable from background

DT on Z will enable more precise diagnosis of the thermonuclear processes occurring during stagnation


	2015	2016	2017	2018	2019	2020
Contained	$D_2, ^3He$	0.1% T	--	--	3% T	3% T
Uncontained	$D_2, ^3He$	$D_2, ^3He$	0.1% T	1% T	1% T	3% T

Physics	Measurement	T Fill %	
		0.1%	1-3%
Behavior of tritium in the Z pulsed power environment	Sampling of tritium contamination, migration		
Thermonuclear scaling (DD vs. DT)	DT yield		
Ion temperature and non-thermal population	Precision nTOF and DT/DD yield ratio		
Liner/fuel mix	DT yield with tritiated gas fill and deuterated liner		
Fuel morphology	Neutron imaging 		
Thermonuclear reaction history	Gamma Reaction History,  Thompson parabola		
Liner/fuel density, non-thermal effects (peak shifts)	Compact Recoil Spectrometer (CRS),  MRS, precision nTOF		


Neutron imaging is finally becoming a reality on Z (at least in 1-D)




1-D Neutron Image



2-D X-ray Image



Relative Brightness
Comparison



- Designed for $Y_{DD} > 3E12$
- Analysis w/ 900 μm binning


X-ray and neutron images aligned to best fit

- need to develop absolute co-alignment




Line VISAR is being developed on Z as a collaboration between LLNL and SNL for assessing current delivery to MDD targets

- Z-pinch exerts pressure on a D₂ gas cell
- Line VISAR resolves pressure vs. radius




Line VISAR is being developed on Z as a collaboration between LLNL and SNL for assessing current delivery to MDD targets

- Z-pinch exerts pressure on a D₂ gas cell
- Line VISAR resolves pressure vs. radius



- Optics and back-end electronics developed and supplied by LLNL
- Beam path, facility mods, and operations provided by SNL



The initial Z Wolter design is complete and we are preparing for commissioning to begin in January


- Wolter optic with coating optimized to image Mo K- α emission region from a z pinch (Ag He- α region in FY18)
- 5 mm x 5 mm x 5 mm Field-of-view (FOV)
- < 100 micron resolution over majority of FOV
- Complex instrument response function will require additional analysis method development

Four-way Collaboration


- Optic designed by LLNL
- Mandrel Fabricated by NASA
- Coatings applied by Harvard
- Z integration by SNL



Pinhole Image




In-chamber Optics mount



Prototype Optic




Modeled PSF



A novel approach to diffraction is under development using crystal optics to capture and relay the diffracted x-rays

- X-ray diffraction is a key capability needed to determine the phase and compression of a sample under high pressure.
- Debris and background pose problems on Z for typical diffraction schemes
- Spherical crystal optics provide a way to relay the signal to a protected detector
- Only practical up to \sim 8 keV
- First test on Z planned for April




New diagnostic capability development on Z is phased through FY20

Diagnostic	Status	FY18				FY19				FY20			
		Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4	Q1	Q2	Q3	Q4
Opacity Spectrometer	Commissioning												
In-chamber MCP PHC	Test & Iteration												
Axial PHC	Test & Iteration												
Neutron Imager	Test & Iteration												
Wolter	* Final Design & Assembly												
Line VISAR	* Final Design & Assembly												
Gated Xtal Backlighting	Final Design												
Diffraction	* Final Design												
Compact Recoil Spectra	Conceptual Design												
Gated Xtal Spectra	Conceptual Design												
Gated Wolter	Not Started												
Z SLOS	Not Started												
Gamma Reaction History	Not Started												
DT-nlImager	Not Started												

Extra

Diagnostic access on Z is limited by the geometry of the pulsed power components



Radial 12°, 0°; Top Axial; Bottom Axial; Chamber

In-chamber diagnostics provide opportunities for improved access and sensitivity

Fiber-Based Diagnostics


- VISAR
- PDV
- Streaked visible spectra

X-ray Spectrometers


- High Energy
~7-70 keV @ $E/\Delta E \sim 500$
- High Resolution
~5-10 keV @ $E/\Delta E \sim 4000$
- MCP-gated
~10-22 keV

X-ray Imagers

- Spherical Crystal Backlighter
- Spherical Crystal Imager




Shock and debris pose major challenges



The Bremsstrahlung background forces long fielding distances and significant line-of-sight shielding for nTOF.

Example nTOF data on Z



Bremsstrahlung background source:

- Originates from power-flow and late-time z-pinch instabilities
- End-point energy > 10 MeV
- Highly variable depending on electrical load design

To mitigate background effects:



- 4" lead shielding
- Collimation
- Fast Mesh PMT's (NSTec) with a large linear current and fast recovery

Diagnostic Operation on Z is heavily influenced by the unique Z environment

- Target chamber must be vented in-between every shot
 - Removal of all in-chamber and axial diagnostics every shot
- Debris environment requires inspection of diagnostic components between each shot
 - Venting of most fixed diagnostics every shot
- Shock environment and replacement of load hardware requires re-alignment of diagnostics every shot
- X-ray Background requires heavy, 4-pi shielding (Typically 1" Tungsten)
- EMP environment requires carefully shielded cables and RF enclosures for all active electronics
 - CCD cameras are challenging to properly shield and operate

Magnetic direct drive (MDD) ICF is studied on the Z Facility

Z Facility




Drive Pressure

$$P = \frac{B^2}{8\rho} = 105 \frac{\pi}{\epsilon} \frac{I_{MA}}{R_{mm}} / 26 \frac{\text{eV}}{\text{J}} \text{ Mbar}$$


Imposed B-fields relax ρr requirements

Contours of positive energy gain at various values of BR



MagLIF experiments on Z have demonstrated the key features of magnetized inertial fusion.


High Convergence Implosion



6+9 keV Emission Image
CR > 40

MagLIF experiments on Z have demonstrated the key features of magnetized inertial fusion.


High Convergence Implosion



6+9 keV Emission Image
CR > 40

B-field Flux Compression


Secondary DT Spectrum



DT Secondary Spectra
BR > 0.35 MG-cm


MagLIF experiments on Z have demonstrated the key features of magnetized inertial fusion.

High Convergence Implosion




6+9 keV Emission Image
CR > 40

B-field Flux Compression



DT Secondary Spectra
BR > 0.35 MG-cm


Thermonuclear Neutrons



$$Y_{DD} = \frac{1}{2} n_D^2 \langle \sigma v \rangle_{DD} V \tau$$

Yield, Volume, Duration
Consistent with DD reactivity

Z Target Diagnostic Capabilities



These broad categories represent hundreds of different diagnostic configurations and fielding locations