


Liquid Electrochemistry in a TEM: Lithium Cycling


Andrew J. Leenheer, Katherine L. Jungjohann, John P. Sullivan, Michael J. Shaw, Kevin R. Zavadil, and C. Thomas Harris
Center for Integrated Nanotechnologies (CINT), Sandia National Labs, Albuquerque, NM

TEM Liquid Cell Design

- A sealed liquid electrochemical cell enables unique *in-situ* experiments on a transmission electron microscope (TEM) in standard volatile electrolytes.
- To link observed morphological changes in electrodes to standard electrochemical cells (e.g. batteries), quantitative current/voltage control is required at fA-pA levels.
- We developed an electrochemical platform optimized for nanoscale measurements and TEM imaging:




Our priorities: narrow electrolyte gap for good resolution, quantitative electrochemistry capability, ability to add active materials, and multiple electrodes for multiple experiments.




- To localize electrochemistry to viewable area, electrodes are masked with ALD Al_2O_3 and patterned using electron-beam lithography.

Ti working/counter/reference electrodes




Nanoparticle assembly also possible:

Dielectrophoresis between starred electrodes resulted in localized immobilization of nanoparticles at electrode tips from a dispersed solution.



Lithium Electrodeposition and Stripping


- Metallic Li is an ideal battery negative electrode, but high-surface-area dendrites cause degradation and safety issues.
- TEM liquid cell enables unprecedented visualization of dendrite initiation conditions and electrodeposition/dissolution dynamics.
- Plate and strip in typical electrolyte (1:1 EC:DMC / 1 M LiPF_6) at typical Li-battery current density: 1, 10, and 25 mA/cm^2 .



Experiment:


- Apply galvanostatic current to 0.26- μm^2 Ti working electrode to induce Li deposition
- Counter / reference electrodes are 750- μm^2 Ti circles
- Take first image halfway through electrodeposition.
- Image periodically through deposition and stripping.
- Low-density Li appears light in BF STEM images.

Cycling at 1 mA/cm^2 : Smooth deposit, high-efficiency cycling.




During electrochemical stripping, dissolution initiates from discrete weak points in surface film rather than uniformly.

Cycling at 10 mA/cm^2 : Needle-like dendrites appear.



Cycling at 25 mA/cm^2 : Larger dendrite growth, poor efficiency.





Electron Beam Effects

- Electron beam induces radiolysis in electrolyte that can interfere with the intended electrochemistry even at minimal dose rates.
- For reactive metal plating, surface films (solid-electrolyte interphase, SEI) critically affect deposition location/morphology.
- Test beam effect: scan in electrolyte near Li at high dose rate, and do plating/stripping with nearly constant imaging.


Imaging dose rate: $10 \text{ e}^- \text{ nm}^{-2} \text{ s}^{-1}$ for 5 s.
Dose rate during test: $3,750 \text{ e}^- \text{ nm}^{-2} \text{ s}^{-1}$.

- Beam creates a SEI surface film upon extended exposure.
- Beam-induced film continues to grow without Li present.
- Scanning beam in dark deposit reduces it back to light-colored Li:



Note: Little beam effect while scanning in electrolyte away from Li.

Cycling at 10 mA/cm^2 while imaging: Uniform, spherical grains.



Compare to deposition morphology without beam exposure.
Beam mediates Li nucleation and growth even with minimal beam current.

Summary

- CINT TEM liquid cell platform development enables visualization of the difficult Li plating process.
- Small exposed electrode area allows quantitative electrochemical measurements, linking observed nanoscale phenomena to bulk electrochemistry.
- Initiation of needle-like Li dendrites visible in TEM and are more pronounced at higher current density.
- Electron beam exposure induces a surface film and modifies the Li plating morphology, causing spherical grains rather than crystalline needles.

Acknowledgements: Special thanks to John Nagan for cleanroom process development.

This work was performed, in part, at the Center for Integrated Nanotechnologies, a U.S. DOE Office of Basic Energy Sciences user facility. Portions of this work were also supported by a Sandia National Lab LDRD project. Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly-owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.

This research was supported as part of the Nanostructures for Electrical Energy Storage, an Energy Frontier Research Center funded by the U.S. Department of Energy Office of Science.