

This paper describes objective technical results and analysis. Any subjective views or opinions that might be expressed in the paper do not necessarily represent the views of the U.S. Department of Energy or the United States Government.


SAND2018-6382C

3D particle location from perspective-shifted plenoptic images



PRESENTED BY

Elise Munz Hall, Daniel R. Guildenbecher, & Brian S. Thurow



Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.




Photo courtesy of Sandia National Labs


Particle tracking for a variety of highly 3D applications

- Fluid flow analysis
- Explosion mitigation
- Measurement of fragment size, shape, velocity

Measure and improve particle location capability

- Quantify/reduce uncertainty
- Reduce computational requirements

3 | Plenoptic Imaging




Camera modified by the insertion of a *microlens array* between the main lens and image sensor

Captures *spatial and angular* information which can be processed to extract 3D information

Refocus and change perspective from a *single snapshot in post processing*

Single compact camera allows for experimental simplification and flexibility



Plenoptic imaging





Image
Sensor


Microlens
Plane

Virtual Image Sensor Plane


Aperture
Plane





5 Refocusing



- Allows a change in the depth of the focal plane
- Resulting images are focused over a range of depths
- Object depth determined based on focus
- Integration = computationally expensive






- Allows a change in perspective within the range of the aperture
- Resulting image as if from a small portion of the aperture
- Object depth can be determined based on apparent motion
- Selection of single pixels = computationally *inexpensive*
- Technique of choice in this work




Perspective-shift algorithm

- Volumetric calibration
- Perspective view particle location
- 3D particle location





- Collect images of dot cards over depth range of interest
- Measure dot positions in perspective views for each image
- Define known dot card positions in object space (X, Y, Z)
- Create polynomial mapping between image sensor coordinates (u, v, s, t) and object space





$$E(s, t) = L(u_0, v_0, P_s(X, Y, Z, u, v), P_t(X, Y, Z, u, v))$$



9 Perspective view particle location



3D particle location





3D particle location








Experimental assessment

-
- Experimental configuration
 - Error definition
 - Accuracy and precision
 - Computational expense



- Simulated static particle field
- Translated to provide known displacement
- Varied nominal magnification
- Large data set allows statistically significant quantitative measurements
- Same data set used in previous depth from focus study

Uncertainty analysis




plot all the measured locations of the selected particle


as a function of the corresponding traverse position

measured depth error

- Average depth displacement error = accuracy
- Standard deviation of depth displacement error = precision



- Smaller average errors with perspective shift
- Narrower confidence intervals



- Smaller standard deviations and narrow confidence intervals with perspective shift
- Perspective shift results more consistent with depth




For a single image, from raw image to 3D particle positions:

by refocusing...



by perspective-shifting...






- Perspective-shifted images show *improved accuracy and precision* as well as dramatically *reduced computational time* over refocusing based methods
- Preliminary results indicate increasing particle density will increase the probability of mismatching particle images between views
- Future work includes:
 - Exploration of alternative matching schemes
 - Use of hybrid refocusing and perspective shift method for depth determination
 - Adaptations to allow/exploit particles with varied size and shape




QUESTIONS?




- Perspective shift allows particle location at greater depths
- Result of location metrics: sharpness vs. location

Sharpness maps from refocusing at center and extreme depths



Perspective views at extreme depth



Depth displacement precision: normalized by theoretical depth resolution

