



Electric Field Induced Second Harmonic Generation Using Ultrafast Lasers for Measurements in Pulsed Nanosecond Plasmas




PRESENTED BY

Benjamin M. Goldberg



Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

Importance of Reduced Electric Field (E/N)



Reuter, S et. al., *J. Phys D* 2018


- In order to understand how we can use the plasma we must have an in depth understanding of fundamental plasma processes**

Electron collisions with neutrals mainly excite internal degrees of freedom

- Leads to non-equilibrium energy distributions

By controlling E/N, we can choose where that energy goes

- Rotational (< 1Td), Vibrational (1-50 Td), Electronic (50 – 100 Td), Ionization (100 Td)



Yuri P Raizer. Gas Discharge Physics


Noninvasive Electric Field Measurement Methods

OES


- Field inferred by monitoring plasma emission
- Requires specific model depending upon gas composition

E-Field CARS

- 4-wave mixing technique for sub-nanosecond measurements
- Resonant technique which has only been demonstrated in H_2 and N_2



A. Sobota et al, PSST 2019




E. V. Barnat et al Appl. Phys. Lett. 2004

LIF Dip


- Able to generate 2D field maps with sensitivity of ~ 10 V/cm
- Pressure limited to ~ 1 Torr due to pressure broadening

Second Harmonic Generation

- Non-resonant method with sub-nanosecond resolution
- Signal scales as E^2
- Sensitivity to electric field vectors



M. Simeni Simeni et al J.Phys. D 2017



B. M. Goldberg et al, Optics Lett 2019

Electric Field Induced Second Harmonic Generation (E-FISH)

In a centrosymmetric medium, second harmonic generation is not possible

- By applying an electric field, that symmetry is destroyed allowing for SHG

- $\mathbf{P}$: Induced polarization
- χ : Electric Susceptibility
- ϵ_0 : Permittivity of vacuum
- $\mathbf{E}$: Electric Field

For SHG, define second order polarization to be: $\mathbf{P}^{(2)}(t) = \chi^{(2)}\epsilon_0\mathbf{E}^2(t)$

Assuming an oscillating electromagnetic field of frequency ω :

- $E(t) = E_0e^{-i\omega t} + E_0^*e^{+i\omega t}$

$$\Rightarrow P^{(2)}(t) = (2\epsilon_0\chi^{(2)}E^*E) + 2\epsilon_0\chi^{(2)}(E_0^2e^{-2i\omega t} + c.c.)$$

E-FISH is described as a third order nonlinear process:

- $P^{(2\omega)} = \frac{3}{2}N\chi^{(3)}(-2\omega, 0, \omega, \omega)E_{Ext}E_{Pump}E_{Pump}$

- $\Rightarrow I^{(2\omega)} = A \cdot N^2(E_{Ext})^2(I_{Pump})^2$

- $P^{(2\omega)}$: Induced polarization at 2ω

- N : Number Density

- $\chi^{(3)}$: Nonlinear Susceptibility (species dependent)

- E_{Ext} : Electric Field to be measured


- E_{Pump} : Electric field of incident laser

Benefits of E-FISH Method:


- Signal scales as N^2
- Signal scales as E^2
- Time resolution determined by pump beam duration
- Spatial resolution determined by beam focusing parameters
- Non-resonant method works in any species and gas mixtures

Dogariu, A., *Physical Review Applied*, 7(2), 024024.

800 nm



400 nm




Early E-FISH Work

Initial similar work completed in 1960's and 70's

- Used nanosecond pulsed lasers, limiting signal generation
- Measured nonlinear susceptibility polarization components
- Measured nonlinear susceptibilities of gasses against known electric fields
- More recently used for detection of Terahertz waves
- We use a more general approach to measure an externally applied field**


Second Harmonic Signal vs $\lambda/2$ Orientation
Measured Signals vs Caprioli's measurement



components in plasmas:

2019

d Flame, 2019



Typical Experimental Schematic

Most experiments utilize either femto- or picosecond lasers

- Fs Laser: 7 mJ per pulse at 1 kHz with 50 fs pulses and ~ 15 nm bandwidth
- Ps Laser: 30 mJ per pulse at 10 Hz with 100 ps pulses and frequency tunable

Detection accomplished with either PMT or CCD depending upon desired parameters


- 1-D sensitivity has been achieved using cylindrical focusing and ICCD

Key experimental difficulties:

- High intensity lasers generate SHG on upstream optics which needs to be filtered
- White light from multiphoton ionization can damage detection equipment
- Cannot rely upon simple spectrometer to separate pump from SHG signal

Discharges studied:

- Bare metal electrodes in 100 Torr room air
- Atmospheric pressure plasma jet with Argon core flow and N₂ co flow
- Dielectric barrier discharges in Fuel and Oxygen containing mixtures





Measurements with Picosecond Laser in 100 Torr Room Air

Used different electrode cross sections to increase local field at anode

- Probe region is center of discharge gap

Electric field follows the applied voltage until breakdown occurs




- Kink in field corresponds to initiation of breakdown
- Field continues to rise as plasma fills the gap
- Field begins to fall once gap is filled with highly conductive plasma



Cylindrical Focusing for 1-D Sensitivity

Alternatively, we can focus into a laser sheet to reduce the net intensity

- Using a camera for detection allows for 1-D measurement sensitivity
- Detect SHG signal using PI-MAX4 ICCD with pixels size of $12.8 \times 12.8 \mu\text{m}$
- Scanned a localized field across laser sheet to determine spatial sensitivity



Resolution is $\sim 50 \mu\text{m}$




Cold Atmospheric Pressure Plasma Jet

Laminar Argon core flow with N₂ co-flow

- Streamer front dimensions are $\sim 30\text{-}100\ \mu\text{m}$
- Plasma travels as an ionization wave, creating spatially and temporally evolving fields
- Measurement volume: $10\ \mu\text{m} \times 220\ \mu\text{m} \times 2\ \text{cm}$

FID GmbH FPG 20-10 NM

- 3.7 kV applied with FWHM of $\sim 30\ \text{ns}$
- Run at 1 kHz synchronized with fs laser
- Time $t=0$ defined to be $\frac{1}{2}$ peak voltage amplitude



3 ns ICCD gate with 1 ns time steps

Red line represents laser probe region

CAP Calibration Routine



CAP does not have a well defined ground plane making calibration difficult

Key concerns:

- Different gas mixture fractions
- Localized field enhancements due to dielectrics
- “Extent” of calibration field compare with field due to plasma jet
- Systematic errors such as Gaussian pump beam intensity

Use cylindrical electrodes using short focal length lens

- Measured and accounted for gas species using Rayleigh scattering
- Removed jet output to remove dielectrics
- Measured confocal parameter to be $220\text{ }\mu\text{m}$ and coherence length to be $\sim 1\text{ cm}$
- Pixel-by-pixel calibration map



Time Resolved Electric Field Measurements

Individual laser shots collected and analyzed in a post processing routine






- Time bins of 500 ps used
- ~100 shots per bin

Ionization wave propagation clearly visible

- Wave speed ~ 0.3 mm/ns, consistent with ICCD images

Peak field ~ 20 kV/cm, much lower than expected

- Likely due to spatial sensitivity as plasma filaments move



3 ns gate with 1 ns time steps

- Averages represent 16 images
- Red line represents laser probe





DBD with Fuel and Oxidizer

Reactor

- Pressure: 60 Torr
- Initial Temperature: 293 K
- Gap distance: 14 mm
- Quartz double dielectric barrier: 1.6 mm
- Argon bath gas with CH_4 and/or O_2


Two pulsers used:

- FID GmbH FPG 20-10 NM
 - Applied Voltage: ~ 7 kV
 - Pulse Duration: ~ 30 ns FWHM
 - Continuous 1 kHz synchronized with fs laser
- FID GmbH FPG 30-50MC4
 - Peak Voltage: ~ 10 kV
 - Pulse Duration: ~ 12 ns FWHM
 - Pulse burst mode: 150 pulses, 30 kHz




1 kHz – Time Resolved Imaging


0.9 Ar + 0.1 CH₄



0.8 Ar + 0.2 O₂



0.7 Ar + 0.1 CH₄ + 0.2 O₂



Averages of 5 shots per frame

Plasmas appear to be uniform, although radial field appears as well

Breakdown occurs on reflected pulse as well

Camera gate: 1 ns

All videos on same intensity scale




1 kHz – Calibration Results

$$I^{(2\omega)} \propto A \cdot N^2 (E_{Ext})^2 (I_{Pump})^2$$

Signal sensitive to $E_{Ext}^2 \rightarrow$ can not determine +/- field


****NOTE**** - PIXIS 512 CCD Used with 500 ms gate

Laser Propagation





1 kHz – Electric Field Results

$0.7 \text{ Ar} + 0.1 \text{ CH}_4 + 0.2 \text{ O}_2$



****NOTE**** - PIXIS 512 CCD Used with 500 ms gate

Laser Propagation



Burst Mode - Time Integrated Imaging

Using high speed sCMOS cameras, we can image each pulse in a burst


Filaments begin forming around pulse 100

- Within a given burst, filaments are stationary
- ****Note**** movement in flow direction due to convection

6.6% CH₄ 13.4% O₂ 80% Ar

End View: Streamwise direction

Pulse Num 0




Side View: Orthogonal direction


Pulse Num 0

Burst Mode - Time Resolved Imaging


93.3% Ar + 6.6% CH₄
Pulse 150



86% Ar 13% O₂
Pulse 150



6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 150



Nanosecond resolution of single pulses




- 1 ns gate, but only 10 Hz camera repetition rate
- Image not from the same burst: uncorrelated shot to shot

Diffuse discharge with just O₂ or CH₄ addition, but significant filamentation with both present




- Filaments stationary within a given burst, but move from burst to burst

Burst Mode – Averaged Electric Field Results




93.3% Ar + 6.6% CH₄
Pulse 150



86% Ar 13% O₂
Pulse 150




6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 150



Single Shot Measurement Capability

Data presented was an average of ~ 100 individual laser shots

- We can look at single shot data for streamer formation!






Conclusions

Electric field induced second harmonic generation is a unique method for non-invasive electric field measurements

- Key benefits include: time response, field vector sensitivity, and species independent measurements
- Measurements have been demonstrated using femtosecond and picosecond lasers

Cylindrical focusing can be used for 1-D spatial sensitivity

- Measurements in a cold atmospheric plasma jet demonstrated an ionization wave travelling with ~ 0.3 mm/ns
- Measurements in DBD demonstrated a uniform diffuse discharge
- Single shot measurements can be used to observe localized streamer formation for plasma instabilities



Current Work – Looking at gas surface interactions using SFG

SHG is a specialized form of **Sum Frequency Generation**

- At the surface interface, inversion symmetry is broken, allowing for SFG
- Used to measure vibrational spectra of molecules at interfaces


Key benefits to SFG:

- Signal generation localized to interfacial molecules between gaseous phase and bulk media
- Can measure the degree of order of interfacial molecules
- Can measure the orientation of the molecule
- Can determine the average tilt angle of the interfacial molecules

$$I_{SFG} \propto \underbrace{\left| \chi_R^{(2)} \right| e^{i\delta}}_{\text{Resonant Contribution}} + \underbrace{\left| \chi_{NR}^{(2)} \right| e^{i\varepsilon}}_{\text{Non-resonant Contribution}}^2$$

$$\chi_R^{(2)} = \frac{B}{\omega_v - \omega_{IR} - i\Gamma}$$

$\chi_R^{(2)}$: Resonant susceptibility
 B : Strength of vibrational mode
 ω_v : Frequency of mode
 ω_{IR} : Probe Frequency
 Γ : Resonant linewidth




Resonant enhancement occurs when $\omega_{IR} = \omega_v$

Experimental Schematic for SFG


Coherent Legend Elite FS laser system

- Use 1 leg to pump a Light Conversion High Energy TOPAS + nDFG
 - frequency tunable pulse from $\sim 2\text{-}10\ \mu\text{m}$ with $\geq 10\ \mu\text{J}$ per pulse
 - 50 fs pulse with $\sim 200\ \text{cm}^{-1}$ bandwidth
- Generate visible light using second harmonic bandwidth compression
 - Split single pulse into two pulses with equal but opposite orientation of chirp
 - Beams recombined on SHG crystal, generating a frequency narrow pulse with $\sim 6\ \text{ps}$ duration

Signal detection done with Newton iStarr ICCD



Initial Results – CH Stretch on Gold Surface



$$I_{SFG} \propto \left| \left| \chi_R^{(2)} \right| e^{i\delta} + \left| \chi_{NR}^{(2)} \right| e^{i\varepsilon} \right|^2$$

Assume:

- $\left| \chi_R^{(2)} \right|$ is given by a Lorentzian
- $\left| \chi_{NR}^{(2)} \right|$ is given by a Gaussian



Used 10 second camera gate with 50 averaged images
 → 500,000 laser shots averaged!!

Nonlinear LSQ fit for following parameters:

- Resonant frequency (constrained by known CH stretches)
- Resonant linewidth (constrained by realistic values)
- Resonant amplitude
- Phase between resonant and non-resonant
- Non-resonant amplitude

Measured and calibrated non-resonant frequency and linewidth using polystyrene calibration film

Initial Results – Measurements on a Steel Wall with a Flame



Acknowledgments

Collaborators Include:

- Professor Richard Miles (Princeton and Texas A&M), Professor Yiguang Ju (Princeton), & Professor Egemen Kolemen (Princeton and Princeton Plasma Physics Laboratory)
- Dr. Stephan Reuter (Polytechnique Montréal) and Dr. Arthur Dogariu (Princeton)
- Aric Rousso (Princeton) and Timothy Chen (Princeton)
- Luka Hansen (Kiel University) and David Feng (Princeton)
- Dr. Christopher Klierer (Sandia National Labs)

Supported in part by MetroLaser, Inc. through a Phase II NASA SBIR

Based upon work supported by U.S. Army Research Laboratory and the U.S. Army Research Office under grant number W911NF1510236

Exxon Mobil, NETL, and the Princeton SEAS innovation fund and ACEE grant


Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.



Thank you for your attention



Questions?




In transparent media, the index of refraction increases as a function of frequency

For second harmonic generation:


- $\Delta k = 2 \cdot \frac{n_{pump}\omega_{pump}}{c} - \frac{n_{SHG}\omega_{SHG}}{c} \neq 0$
- $I_{SHG} \propto \text{sinc}^2\left(\frac{\pi\Delta kL}{2}\right)$

If measurement length is greater than $2/\Delta k$, efficiency goes down as second harmonic goes out of phase with driving polarization

$$Z_R = \frac{\pi w_0^2}{\lambda} = \frac{4\lambda}{\pi} \left(\frac{f}{D}\right)^2$$



$$Z_R I_{2\omega}^2 \propto \frac{\pi W_0^2}{\lambda} \sin^2 \left(\frac{4\lambda L}{\pi L_{coh} D} \right)^2$$



Focal Length	Rayleigh Range
150 mm	0.894 mm
300 mm	3.575 mm
500 mm	9.93 mm
1000 mm	39.73 mm

Benefits of E-FISH Method

Quadratic Dependence upon Electric Field

- Sensitive down to 100's of V/cm

Time Resolution

- Fs laser acts as a δ function compared with ns changes in fields

Measurement length


- Determined by beam focusing parameters allowing for sub-mm spatial resolution

Species Independence

- Non-Resonant technique that can be used in any gaseous mixture

Field Vector Sensitivity

- Nonlinear susceptibility polarization components indicate the ability to measure field vector components





Dogariu, A., *Physical Review Applied*, 7(2), 024024.

APPJ does not have well defined ground


Calibration completed between cylindrical electrodes

$$E(z) = \frac{2 \cdot V_{Applied}}{d \cdot \cosh^{-1} l/2a} \cdot \frac{1}{\left[1 + 4 \frac{(z-z_0)^2}{d^2}\right]}$$




Burst Mode - Time Resolved Imaging


6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 25



6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 100



6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 150



Nanosecond resolution of single pulses




- 1 ns gate, but only 10 Hz camera repetition rate
- Image not from the same burst: uncorrelated shot to shot

Shows evidence of significant filaments by pulse 100




- Non-uniformity moves from burst to burst

Burst Mode – Averaged Electric Field Results




6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 25



6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 100




6.6% CH₄ 13.4% O₂ 80% Ar
Pulse 150



1 kHz – Electric Field Results

****NOTE**** - PIXIS 512 CCD Used with 500 ms gate

Laser Propagation



$0.7 \text{ Ar} + 0.1 \text{ CH}_4 + 0.2 \text{ O}_2$

