




# In situ ion irradiation transient grating spectroscopy (I<sup>3</sup>TGS)




Cody A. Dennett<sup>1</sup>, Khalid Hattar<sup>2</sup>, Michael P. Short<sup>1</sup>


<sup>1</sup>Department of Nuclear Science and Engineering, Massachusetts Institute of Technology, <sup>2</sup>Center for Integrated Nanotechnologies, Sandia National Laboratories

## Transient Grating

TGS offers laser-based, non-contact, non-destructive evaluation of acoustic, elastic, and thermal properties over micron length scales




- Impulse periodic surface excitation from pulsed laser causes transient changes in temperature and displacement



- Probing laser monitors the decay and oscillation of these excitations


- Surface acoustic wave (SAW) frequencies are easily extracted and provide a link to elastic mechanical properties




- Fitting to a thermal equilibration model provides a direct measure of thermal diffusivity

## I<sup>3</sup>TGS Beamline @ Sandia


Laser diagnostic installed on a dedicated beamline on the 6 MV tandem ion accelerator at SNL



top view of beamline with laser enclosures side view of control system



- Combining a TGS diagnostic with an ion beamline allows for real-time monitoring of material property changes during irradiation



I<sup>3</sup>TGS system schematic for laser and ion beams

- Properties are measured with orders-of-magnitude higher resolution in applied dose than traditional post-exposure examination


- Microstructure effects are related to property evolution through theory, prior testing, or coordinated modelling

## High Temperature Self-Ion Irradiation


Initial experiments demonstrate void swelling can be monitored via changes in SAW speed

material	Single crystal {001} Ni
ion species	Ni <sup>5+</sup>
ion energy	31 MeV
temperature	550°C
spot size	1.6 mm dia.
beam current	33.7 nA (avg.)
acoustic polarization	13° from <100>
grating spacing	4.54 μm

- Ion irradiation conditions set to induce void swelling
- Length scale between the damage layer depth and measurement sensitivity is matched by tuning the ion beam energy




- Ion beam current, sample temperature, SAW speed, and thermal diffusivity monitored continuously during day-long exposures



- The observed increase, followed by decrease in SAW speed is characteristic of void swelling in single crystal FCC metals\*

- No change in diffusivity observed – damage localized beneath thermal penetration

\*Dennett et al., Acta. Mater. 145 (2018)




post exposure cross-sectional electron micrographs confirm large faceted voids

## Capabilities Under Development

More user-friendly alignment and data acquisition tools – coming online as a user facility soon!

- Upgraded ion beam diagnostics and imaging for improved coincidence positioning
- Platen-based sample transfer for easy exchanges
- Azimuthal sample rotation for alignment and exploration of higher-order elastic behavior



new high temperature sample stage with manual azimuthal rotation

## Acknowledgement

Facility Description Dennett et al., NIMB 440 (2019)

Contact MIT Cody A. Dennett – cdennett@mit.edu Sandia Khalid Hattar – khattar@sandia.gov

This work was performed, in part, at the Center for Integrated Nanotechnologies, an Office of Science User Facility operated for the U.S. Department of Energy (DOE) Office of Science. Sandia National Laboratories is a multi-mission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International, Inc., for the U.S. DOE's National Nuclear Security Administration under contract DE-NA-0003525. The views expressed in the article do not necessarily represent the views of the U.S. DOE or the United States Government.