





Population-Based Comprehensive Health Monitoring For Combating Infectious Disease and Bioterror Threats

Conrad D. James, Ph.D.
Biosensors and Nanomaterials Dept.
Sandia National Laboratories



Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000.



Acknowledgements



**Sandia
National
Laboratories**

Peptide Design

Cathy Branda, Team Lead
Peter Anderson
Joe Schoeniger

Microfabrication

Conrad James, Team Lead

David Santos
Jason Harper
Adrian Schiess
Michael Brumbach
David Wheeler
Mark Kinnan

Thayne Edwards
Komandoor Achyuthan
Tony Ohlhausen
Ed Heller
Jaclyn Murton
Katlin Schroeder

Systems Analysis

Julie Fruetel, Team Lead
Margot Kimura
Dave Reichmuth
Trisha Miller

Paula Krauter
Todd West



Prof. Neal Woodbury
Prof. Stephen Johnston



Innovations In Medicine:

Stephen Johnston

Neal Woodbury	Lucas Restrepo	Valeriy Domenyuk
Chris Diehnelt	Nidhi Gupta	John Lainson
Phillip Stafford	Zhan-Gong Zhao	Meredith Blynn
Alex Carpenter	Donnie Shepard	Kari Kotlarczyk
Bart Legutki	Andrey Loskutov	Kristen Seifert
Penny Gwynne	Jinglin Fu	Loren Howell
Douglas Daniel	Rebecca Halperin	



Tom Slezak, LLNL Team Lead
Beth Vitalis
Kevin McLoughlin
Adam Zemla
Jimmy Thissen
Shea Gardner
Ana Paula De Oliveira Sales
Dan Faisso
Tom Edmunds
Carl Melius
Richard Wheeler



Bill W. Colston
David Smith
Kathy Sykes
Robert Kumelis

John Galgiani: U. of Arizona
Hoda Anton-Culver: UC Irvine
Sam Hanash: Fred Hutchinson Cancer Center
Adi Gazdar: UT Southwestern
Adrienne Scheck: Mayo Clinic
Dawn E. Jaroszewski: Mayo Clinic
Matthew Grieving: Nextval, Inc.
George Poste: Complex Adaptive Systems Initiative





Infectious disease scenario

An outbreak of symptoms among a large population after an event (dust storm, etc.)



Massive Dust
Storm Hits
Phoenix, AZ
July 5th, 2011

<http://www.youtube.com/watch?v=8W4Cx44XKZ4>




Coccidioides immitis spherule
with endospores



Coccidiomycosis:

- About 10,000 reported cases annually
- Mostly mild cases, can be life threatening
- Flu-like symptoms

Need: a biodetector for rapid detection of *ill* and *near-ill* individuals



Immunosignaturing (invented by S.A. Johnston and P. Stafford):


- *Produce a longitudinal profile of individuals to identify the ill and predict individuals that will become ill*
- Only microliters of blood needed; Easy - sample just diluted
- Can be self administered, minimally invasive
- Convenient - can be mailed on filter paper

Legutki, J. B., D. M. Magee, P. Stafford, S.A. Johnston (2010). Vaccine 28, 4529-4537
Kukreja, M., S.A. Johnston, P. Stafford (2012). J Proteomics Bioinform (6), 10.4172/jpb.S6-001
Stafford, P., R. Halperin, J.B. et al., (2012). Molecular & Cellular Proteomics 11(4), 10.1074/mcp.M111.011593


HealthTell™

ASU BIODESIGN
INSTITUTE
ARIZONA STATE UNIVERSITY

Challenges in biodetection and biosurveillance




Blood dilution problem...



Hori and Gambhir, *Sci Transl Med*, 2011

The immune system detects & amplifies...



1000s of Abs produced/min;
 $\sim 10^8$ Abs in serum

B-cells respond quickly/specifically to changes in health by producing Abs


Sequence B- and/or T-cell variation in blood or...

Profile Abs over time

Immunosignaturing with peptide microarrays



- Standard fabrication instruments
- 330K peptide per array (assay)
- 312 assays per wafer





High density → High information content → High volume → Low cost

Detection of Valley fever



- 10K peptides on original array, 100 peptides selected
- 90 blinded samples from patients at the clinic
- Zero false positives (100% specificity)
- Zero false negatives (100% sensitivity)

All Patients with Valley Fever Presented with Zero CF Titers, but were later shown to have the disease






Discriminating cancers with the peptide array technology

Example array with 8 micron spots



Technology maturation drivers

- Manufacturability – product readiness for the market
- Robustness – consistent performance across platforms
- Quality assurance – accuracy of synthesized peptides



Technology maturation: U.S. National Laboratories


Peptide design & performance




Process development



Characterization



Sample testing & assay optimization



Exploratory production



Low/mid-level production




Mature, optimize, and innovate



Sandia
National
Laboratories

Pilot peptide array manufacturing




Microarray fabrication process- photolithography and reaction chemistry




SiO₂ on Si




Aminosilane




Boc'd linker




Photoresist coating



Light exposure




Photoresist removal




L Coupling




Repeat - expose and couple AAs



YGGFL 5mer build






Chemical deprotection




2nd Ab with fluorophore

5mer peptide fabrication with control deletions




- Vehicle for testing/optimizing new peptide coupling chemistry
- Positive and negative controls




Stability of coupling reagents in a high-volume production environment

Examine reagent stability under
storage conditions; potential need
for *in situ* mixing





in situ ellipsometry for non-destructive sample assessment

Film analysis for sample-to-sample variability and longitudinal tracking of chemical processing




sample	film thickness (Å)	condition
1	155.25	oxide
2	157.21	oxide
3	157.24	oxide
4	154.78	oxide
5	160.83	silane
6	162.17	silane
7	163.93	silane
8	164.51	silane




Quantitative measurement of amine-group surface density

Surface chemistry of peptide arrays:



UV/VIS/NIR Spectrometry:




Affinity and avidity – impact on peptide array performance

Quantitative measurement of
molecular surface density:

$$\Gamma = 1/2 \left[\frac{A_\lambda N_A}{\epsilon_\lambda} \right] \text{ (molecules} \cdot \text{cm}^{-2}\text{)}$$





Surface	Molecules/cm ² [10 ¹²]	Surface amines [%]
Total free amines	9.12	100
Coupled Met-Boc	7.25	79.5
Neutralized	1.11	12.2
Unreacted free amines	0.76	8.4
Deprotected amines	4.54	49.8



Many factors to explore to optimize Ab binding to
peptide array: density, linker length, etc.



Peptide deletion control with photolithographic confinement

Photolithographic confinement of the aminosilane linker - reduces variation in array element size, improving array repeatability, reduces occurrence of deletions




Spectroscopic analytical techniques for quality assurance and control



X-ray photoelectron spectroscopy (XPS):






- XPS provides elemental and chemical surface analysis; ratiometric analysis to handle SNRs



Analysis of the substrate chemical state after peptide array processes





Quantitative quality control for peptide fabrication




Identification of peptide signatures with mass spectrometry

- Time-of-flight secondary ion mass spectroscopy (TOF-SIMS) analyzes material *ion fragments* from ion bombardment; high mass resolution (ppt)




Gly

Chemical Formula: CH_4N^+
 m/z: 30.03 (100.0%), 31.04 (1.1%)



Chemical Formula: $\text{C}_2\text{H}_6\text{NO}^+$
 m/z: 60.04 (100.0%), 61.05 (2.3%)
 G-Immonium Ion




Chemical Formula: $\text{C}_5\text{H}_{12}\text{N}^+$
 L-Immonium Ion:

Henry and Bertrand, *Surface & Interface Analysis* **2009**, 41, 105
 Sole-Domenech, et al., *Analytical Chemistry* **2010**, 82, 1964




Identification of elemental & molecular signatures of amino acids







Tracking of peptide synthesis with TOF-SIMS






<http://pepbank.mgh.harvard.edu/interactions/details/57761>



TOF-SIMS imaging for peptide content analysis


- Scanning mode, 5 μm /pixel
- Identify amino acid secondary ion signatures




$\text{C}_4\text{H}_8\text{N}^+$ 70D (P Immonium Ion)

Correlation of ion signatures to peptide array map

Phenylalanine map of peptide array:




TOF-SIMS scan for $C_8H_{10}N^+$ 120D
(F Immonium Ion)



Ion signature intensity correlates to amino acid position

TOF-SIMS scan for $C_8H_{10}N^+$ 120D
(F Immonium Ion)



	1	2	3	4	5	6	
	1	1	1	1	1	1	



Spot	peptide	TOF-SIMS	
1	LNWFGSGG	84.3	17
2	LPHFGSGG	65.0	13
3	LWFNGSGG	96.7	18
4	LHFKGSGG	74.6	15
5	LFNQSGGG	105	25
6	LFKEGSGG	134	36

Ion intensity varies with position, copy number, and environment

- Develop an *expectation score* for the ion intensity of an amino acid as a function of position within a peptide, peptide length, and copy number

$$S_{AA} = \frac{n_{AA}}{l_{pep}} \sum_{AA_i}^i (l_{pep} - p_{AA_i} + 1)$$



Spot	peptide	S _F	TOF-SIMS	
1	LDYAWFFYLRRHRGGSGG	2.71	73.4	13
2	LLVYVGNPRREQEGSGG	0	4.9	1.2
3	LSARWRPRRYEYNGSGG	0	7.5	1.5
4	LWWFKSNGFYDGSGG	2.53	69.8	20
5	LVGGHSLHLWRGSGG	0	7.3	1.5
6	LRHPNPFFKEHPNFGSGG	4.89	43.8	10
7	LWDYSVWPGAKGWGGSGG	0	4.0	1.6
8	LHAHGNNPHFNLHPQVWGSGG	0.60	11.4	2.1

Conclusions

- Peptide microarray immunosignaturing technology has immense promise for providing unique biodetection capabilities
- Additional characterization and optimization are necessary to further develop this technology and maximize impact
- Combined with complementary analytical techniques (e.g. genetic sequencing), new operational applications can be developed

