





Rapid Automated Portable Infection/exposure/disease Dx RAPIDx

SAND2008-5963P

*Greg J. Sommer
Biosystems Research Dept.*



Better, Faster, Cheaper, Smaller



Sandia is a multi-program laboratory operated by Sandia Corp., a Lockheed Martin Co., for the United States Department of Energy under Contract DE-AC04-94AL85000.



Acknowledgements

Funding Sources:

NIDCR (Saliva Diagnostics)

NIAID (Toxin Diagnostics)

LDRD (Biomarker Verification, Deployable Dx Pathogen Exposure)

Key Contributors:

Anson Hatch

Greg Sommer

Ying-Chih Wang

Anup Singh

Amy Herr (now at UC Berkeley)

Dan Throckmorton

Ron Renzi

Jim Brennan

Gabriela Chirica

Bal Ram Singh (UMASS Dartmouth)


Easwaran Ravichandran (UMASS Dartmouth)

Will Giannobile (U Michigan)

Steve Binder (Bio-Rad)

Exposure

City



Crowd



Field




Versatile Platform

Portable,
Rapid Triage



Low Cost, Non-invasive,
Convenient and Easy to use



8" x 8" x 4"

BSL Compatible



Screening

High Risk Traffic




Routine Clinic Visit



Research

Biomarkers of Toxin Exposure

Detection and treatment windows following exposure



BoNT A,B, E

Ricin

SEB

Shiga I,II

Epsilon toxin

IL-6

IL-2

TNF α



IFN γ

IgM

IgG

Deployable System


Powerful & Adaptable




Salient Features

- Microfluidic platform: gel shift assays
- Self-contained
- Easy to use
- Microliters of blood or saliva (non-invasive)
- Rapid (<10 min sample to result)
- Multiplexed (up to 32 assays simultaneously)
- Low cost per assay
- Low risk to operator
- Sensitive (0.1 pg/ml, 10 fM, Detection limit)


Microfluidic Immunoassays




*Electrophoretic
Separation*




Microfluidic Immunoassays



*Electrophoretic
Separation*




Microfluidic Immunoassays






- Solution phase technique
- No antibody immobilization
- No secondary antibody required
- Rapid, reproducible binding
- Simple electrokinetic manipulations on-chip
- Only specific Ab probes are labeled
- Clean baselines and easily interpreted electropherograms

*Electrophoretic
Separation*



Gel/Membrane Fabrication




Hatch et al., *Analytical Chemistry*, 2006

Meeting Needs for Protein Analysis

- Low abundance markers: need sub-pM detection
- Direct detection limit of LIF system is ~100 pM
- Solution preconcentration
 - Also simplifies on-chip mixing
 - Favorable kinetics for 10^{-9} or lower affinity constants




Hatch et al., *Analytical Chemistry*, 78, 4976-4984, 2006



Microfluidic Immunoassays

Load Specific Antibody (labeled)

- *Specific to toxin or biomarker*




Microfluidic Immunoassays

Load Specific Antibody (labeled)

- *Specific to toxin or biomarker*

Load and Mix Sample

- *Toxin concentrated at size exclusion membrane*
- *Binds to antibody (solution phase)*



Microfluidic Immunoassays

Load Specific Antibody (labeled)


- *Specific to toxin or biomarker*

Load and Mix Sample

- *Toxin concentrated at size exclusion membrane*
- *Binds to antibody (solution phase)*

Electrophoretic Separation

- *Detect level bound*



Microfluidic Immunoassays

Load Specific Antibody (labeled)


- *Specific to toxin or biomarker*

Load and Mix Sample

- *Toxin concentrated at size exclusion membrane*
- *Binds to antibody (solution phase)*

Electrophoretic Separation


- *Detect level bound*



Biomarkers of Disease





- Herr et al, *PNAS* 2007
- Non-invasive **saliva** sampling
- Statically validated MMP-8 as biomarker of periodontal disease
- **Presymptomatic** diagnosis, higher treatment efficacy
 - Traditional diagnosis based on pocket depth, bone loss






Diagnostics of Exposure

Preconcentration for fM detection





Hatch et al., *Analytical Chemistry*, 2006

Multiplexed Panel of Toxins, Host Response Markers – NIAID Funded U01





Meeting Needs for Protein Analysis

- Wide Channel Preconcentration:
 - Direct laser writing of membrane in channel
 - Seamless integration with separation channel
- Detection of 10 fM analyte (direct labeled)
 - >4 orders of magnitude preconcentration
 - <5 minute step
 - S/N ratio of 28





Wang et al., *microTAS* 2008

Botulinum Neurotoxin Enzyme Activity Assay



Botulinum Neurotoxin Enzyme Activity Assay



On-Chip Process


- Preconcentrate BoNT/A, mix probe, react, separate

Specific to active toxin

Signal amplified vs. direct-labeled immunoassay

- Detection limit ??

Spatial Multiplexing





8 Modules




1" Square
Chip


Simultaneous separations with all 8 modules x 4 analytes per



Spatial Multiplexing: Rotary Optical Scanner





Rotary Optical Scanner




Brennan et al.

Disposable Assay Cartridge



Bottom view




Top View

- isolated reagent storage
- valve-sealed gels on-chip until use
- simple loading of cartridge and sample
- concept amenable to low-cost mass fabrication

Approximate dimensions: 3 x 3 x 1.5 cm



Current Status



- TRL4 Device ~ 1.5 years
- On-going development of toxin and host response marker assays
- On-chip and off-chip sample pretreatment for complex samples
- Disposable cartridge engineering

Microscale Immobilized pH Gradients (μ IPG's)

Diffusion-controlled pH gradient immobilized within polyacrylamide matrix on-chip






pI Markers




AF647 Proteins

AF488 Proteins

Pore Limit Electrophoresis



Alternative Affinity Reagents



- Cytokines at size threshold for immuno-gel-shift assay
- Much improved with smaller affinity reagent such as affibody, aptamer
- Advantages in uniformity, production