


## Background


- Algal biofuel production cost can be reduced through elimination of pest-induced “pond crashes”<sup>1,2</sup>
- Current methods to eliminate pests include regular additions of pesticides, incurring a consistent cost<sup>3</sup>
- The cultivation of microalgae with unique bacterial consortia results in the death of their rotifer predators,<sup>4</sup> and potentially additional pests such as chytrids and phagotrophic algae
- Protective molecules produced in microalgae-bacterial consortia co-cultures can be identified with metabolomics tools



## Crop

## Pests

### Biofuel Microalga



*Microchloropsis salina*

### Rotifer




*Brachionus plicatilis*

### Chytrid



*Rhizophydium* sp.

### Phagotrophic Alga




*Poterioochromonas malhamensis*

- Hypothesis:** Select bacteria in our consortia produce small molecules that protect microalgae from predators and pathogens.
- Objective:** Identify protective small molecules produced by bacteria isolated from our unique consortia.


## Experimental Design

### Metabolomics



**Figure 1:** Workflow for identifying small molecules produced by protective bacteria in microalgae-bacterial co-cultures. Metabolites are extracted from filtered algal culture conditioned media and analyzed using spectroscopic methods (UHPLC-MS and <sup>1</sup>H NMR) prior to molecule identification.



### Bacterial Extractions for Bioassay Screening



**Figure 2:** Workflow for generating bacterial isolate chemical extracts. Extracts will be tested against chytrids and phagotrophic algae.


## Results

### Exometabolome Extracts of Different Algae-bacterial Co-cultures Have Different Chemical Profiles



**Figure 3:** Left) Algal culture growth curves for bacterial consortia experiment. Above left) Liquid chromatography mass spectrometry (LC-MS) chromatogram for an algal-bacterial (strongly protective) culture conditioned media extract. The mass spectrum for peak A (gold outline) shows m/z signals for distinct metabolites (red circles). Above Right) LC-MS chromatogram for an algal-bacterial (moderately protective) culture conditioned media extract. The mass spectrum for peak B (blue outline) shows an m/z signal for a distinct metabolite (red circle).


### Descriptive Statistics Indicate that Strata-X is the Best Extraction Method for Profiling Algal Culture Exometabolomes



**Figure 4:** Left) Algal culture growth curves for *Alteromonas macleodii* experiment. This bacterium was selected to determine whether it could protect algae from rotifers like the bacterial consortia. Middle) Principal component analysis of reversed phase C18 LC-MS data acquired for algae + rotifers (gold) and algae + bacteria + rotifers (blue) conditioned media extracts. Three different solid phase extraction methods are represented by triangles (C8), squares (PPL), or circles (Strata-X). Right) Loadings plot where each point represents a molecular mass with a corresponding retention time. These points are responsible for the separation and grouping of treatment groups.

## Conclusions

- Exometabolome extracts of algal cultures grown with moderately and very protective bacterial consortia have different chemical profiles
- Extracted molecules appear to be polar organic metabolites
- Metabolite concentrations are very low
- Strata-X is the best solid phase extraction method for extracting metabolites from algal culture conditioned media



## Future Work

- Determine chemical features in MS metabolomics spectra unique to microalgae-bacterial consortia co-cultures that did not crash in the presence of rotifers
- Bioassay-guided fractionation/isolation, spectroscopic characterization, and structure elucidation of protective molecules
- Acquisition of full biological profiles for protective metabolites
- Development of a bioassay and testing the bioactivity of bacterial isolate metabolites against chytrids and phagotrophic alga *P. malhamensis*

## References and Acknowledgements

- Sun A et al. (2011) Energy, 36, 5169-5179.
- Hamilton CE and Rossmeissl N. (2014) Department of Energy.
- Smith VH and Crew T. (2014) Algal Research, 4, 23-34.
- Fisher CL et al. (2019) 40, 101500.

I would like to thank Emily, Anne Marie, Bhuwan, and Sam for their advice, Dr. David Gaul for his assistance with LC-MS, and my undergraduate Hailey for her help with method development. This work was supported by LDRD funding from a collaboration with Sandia National Laboratories. SAND No. SAND2020

Sandia National Laboratories is a multi-mission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525. This poster describes objective technical results and analysis. Any subjective views or opinions that might be expressed in the presentation do not necessarily represent the views of the U.S. Department of Energy or the United States Government.