

Simultaneous 3D CH₂O LIF and tomographic PIV measurements at 10 kHz in lifted turbulent jet flames


B. Zhou¹, J. H. Frank¹, E. Huang¹, T. Li², B. Böhm², A. Dreizler²

¹Combustion Research Facility, Sandia National Laboratories, Livermore, California

²Reaktive Strömungen und Messtechnik, Technische Universität Darmstadt, Germany


Turbulent combustion involves thousands of chemical reactions that interact with turbulent flows over a wide spectrum of spatio-temporal scales. These interactions are highly non-linear, and inherently evolves in both spatial and temporal dimensions. Expanding diagnostic capabilities from conventional 10 Hz two-dimensional imaging to time-resolved three-dimensional (3D) visualization of key quantities (e.g. velocity and species) is a necessity to understand the dynamics of turbulent combustion phenomena. We demonstrate development of a new capability that extends species imaging (CH₂O) to time-resolved 3D measurements by employing a 100 kHz pulse-burst laser and an acoustic-optic deflector scanning system in a lifted turbulent jet flame. Temporal dynamics of the complex CH₂O topology and its interaction with the flow field in lifted turbulent jet flames are needed to better understand the flame propagation/stabilization mechanisms.

Tomo-PIV + 3D-LIF Imaging by Raster Scanning of Burst Laser




Jet Mixture	Jet Bulk Vel.	Jet Phi	Coflow Mixture	Coflow Bulk Vel.
DME/Air	13.5 m/s	9.5	Air	0.2 m/s

Acousto-Optic Deflector (AOD)




Synchronization of burst-pulse laser, AOD and TPIV laser




High-speed 3D Measurement of Flow Field and CH₂O LIF


CH₂O and Vel. Fluctuation v'




CH₂O and Vorticity Mag.



CH₂O and Strain Rate Norm



Auto-ignition Event



Summary

- Demonstrated feasibility of high-speed imaging of velocity field with 3D LIF imaging of formaldehyde using laser raster scanning
- Temporally correlated CH₂O and turbulence structures have been found. CH₂O pocket formation and propagation has been identified, suggesting occurrence of auto-ignition as a stabilization mechanism
- Next steps for Tomo-PIV/3D-LIF imaging data: Analyze dynamics of flame stabilization, propagation, and full 3D strain alignment

Acknowledgments

Sandia authors gratefully acknowledge the support of the US Department of Energy, Office of Basic Energy Sciences, Division of Chemical Sciences, Geosciences, and Biosciences. Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC., a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA-0003525.