

# Pulsed gradiometry in Earth's field

This research was developed with funding from the Defense Advanced Research Projects Agency (DARPA). The views, opinions and/or findings expressed are those of the authors and should not be interpreted as representing the official views or policies of the Department of Defense or the U.S. Government.

Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.


**Kaleb Campbell**

**Sandia National Labs, University of New Mexico  
Department of Physics and Astronomy**


# Applications of OPM's

## Magnetoencephalography (MEG)


- Magnetic fields generated from neuronal activity in the brain
- Fields in the femto-tesla to pico-tesla range - requires very sensitive magnetometers
- SQUID magnetometers are most sensitive - require magnetic shielding or field cancellation coils



5-sensor array




Partially covers the left hemisphere



Insert Person Here


# Basic Idea – Magnetic Gradiometer



\*Henry Tang, *Parametric Frequency Conversion of Resonance Radiation in Optically Pumped  $^{87}\text{Rb}$  Vapor*. Phys. Rev. A 7, 2010 (1973).


\*Vishal Shah, *System and Method for Measuring a Magnetic Gradient Field*. Patent. US10088535 (2018)

# How Sidebands are produced



- Atoms absorb angular momenta from  $\sigma^+$  light
- Make  $\Delta m_F = +1$  transitions to excited state

# How Sidebands are produced



- Atoms absorb angular momenta from  $\sigma^+$  light
- Make  $\Delta m_F = +1$  transitions to excited state

# How Sidebands are produced




$|F = 2, mf = 2 >$




$|F = 1, mf = 1 >$

# How Sidebands are produced



- $\mu W$  field Induces Rabi oscillations between levels
- Rabi frequency  $\Omega = \frac{\mu B}{\hbar}$
- $\Omega t = \pi/2$  known as  $\pi/2$  Pulse

# How Sidebands are produced




- Atoms are in a coherent superposition of both states

- Coherences are maximized

$$\rho = \begin{bmatrix} \rho_{11} & \rho_{12} \\ \rho_{21} & \rho_{22} \end{bmatrix} = \begin{bmatrix} 1/2 & 1/2 \\ 1/2 & 1/2 \end{bmatrix}$$

- Atoms oscillate at the ground state hyperfine frequency (6.834 GHz for  $\text{Rb}^{87}$ )
  - Nuclear spin and Electronic spin tumbling around each other in phase

# How Sidebands are produced



- Introduce Carrier (probe) light to the experiment

From Maxwell's equations (No current source, magnetic susceptibility, or spatial charge)

$$\nabla^2 E + \nabla(\nabla \cdot E) = \frac{-1}{c^2} \frac{\partial^2 D_e}{\partial^2 t}$$

where

$$D_e = E + 4\pi P$$


The electric displacement is related to the electric dipole moment of the atoms

$$D_e = E + 4\pi N \langle D \rangle$$

Or

$$D_e = (1 + 4\pi\chi)E$$

# How Sidebands are produced



$$\nabla^2 \tilde{\mathbf{E}} - \nabla(\nabla \cdot \tilde{\mathbf{E}}) = -\frac{\omega^2}{c^2} \left( \tilde{\mathbf{E}} + 4\pi n_a \langle \tilde{\mathbf{D}} \rangle \right)$$

- For a one-dimensional plane wave

$$\left( \frac{\partial}{\partial z} \pm i \frac{\mathbf{W}}{c} \cdot (1 + \chi/2) \right) \cdot \tilde{\mathbf{E}}_{\pm}$$

$$\langle \omega' | \chi | \omega \rangle = 4\pi n_a \sum_{\mu, \nu} \alpha_{\mu\nu} \tilde{\rho}_{\mu\nu} \delta_{(\omega' - \omega), \Omega_{\mu\nu}} e^{\frac{i}{c}(\omega' - \omega)z}.$$

- Sidebands are produced offset from the carrier beam by the hyperfine frequency (frequency of modulation)


$$\omega_+ = \omega_0 + \mu W$$

$\omega_0$  = Carrier frequency


$$\omega_- = \omega_0 - \mu W$$

- Similar to an Electro-optic modulator (EOM)

# Building a numerical model



Zoomed In On Sidebands






# Pulsed Experiment



Use an AOM as a switch on Pump

Pump is off when making measurement


- Pump doesn't degrade coherence
- Stronger pump can be used (better pumping)




# Building a numerical model



# QuSpin Physics Package Results



Gradient Noise Measurement (Earth's Field)




Photon Shot Noise:  $\sim 6 \text{ fT/cm/rt-Hz}$

Data and photo courtesy of QuSpin

# Compact single 795nm laser setup




- Potentially halve number of lasers needed
- Simpler experimental Setup



## How we do it


- Replace the 780nm probe with an attenuated pump
- Use an etalon to distinguish between carrier and sidebands
- Temperature/angle tune the etalon to the sideband mode



# Can we make the sensor dead-zone free?


## Case 1: Perpendicular field

- How can we optically pump to the end state?





## Case 2: Parallel field

- Only  $\sigma$  optical transitions allowed
- Requires  $\Delta m_F = 0$  microwave transition
- How to prepare in the  $|F=1, m_f = 1\rangle$  state?




# Case I: Ambient field perpendicular to the laser axis


We rotate the quantization axis, while maintaining the atomic population in the  $|2,2\rangle$  state.




Spin vector must follow torque vector to be adiabatic



To be adiabatic, the rotation rate of the field must be less than the Larmor precession frequency  $\frac{d\theta}{dt} \ll \omega_L$ .




Selection rules now allow  $\Delta m_F = 1$  transitions



## Case 2: Ambient field parallel to the laser axis

- We perform adiabatic rapid passage to transfer the population from the  $|2,2\rangle$  state to the  $|1,1\rangle$  state.
- We use a magnetic field ramp to simplify microwave and cover both cells.




Magnetic field must be ramped slow enough that the Spin vector follows the Torque vector

Selection rules allow  $\Delta m_F = 0$  transitions


$$F = 2$$

$$F = 1$$

$$m_f$$




# First Single Laser Experiment




Scan the probe frequency Between the F=1 and F=2 resonances





Scan the frequency between F=1 and F=2 during probe phase



# First Single Laser Experiment



Scan the probe frequency Between the F=1 and F=2 resonances



# Acknowledgements

**QUSPIN**



Vishal Shah-PI

Ying-Ju Wang




**Sandia  
National  
Laboratories**



Yuan-Yu Jau



Peter D.D. Schwindt



Igor Savukov

## Funding

**DARPA -Atomic Magnetometer for  
Biological Imaging In Earth's Native  
Terrain (AMBIENT)**

This research was developed with funding from the Defense Advanced Research Projects Agency (DARPA). The views, opinions and/or findings expressed are those of the authors and should not be interpreted as representing the official views or policies of the Department of Defense or the U.S. Government

Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.




# Thank You!

# How bad will this be?



## Beatnote Function

$$S(t) = \frac{E_1^2}{2} e^{-2t/T_1} + \frac{E_2^2}{2} e^{-2t/T_2} + E_1 E_2 e^{-t(\frac{1}{T_1} + \frac{1}{T_2})} \sin(2\pi f t + \varphi)$$



- $FOM = \sqrt{S}/w_2$
- Need to test this out with experimentally
- Square root dependence on amplitude
- This method is more efficient than ARP for small B-field gradients

# Dead-zone free operation


## Minimum field

- $f_{1 \rightarrow 2} - f_{1 \rightarrow 1} \gg \frac{1}{2\pi T_{\pi/2}}$
- For  $T_{\pi/2} = 0.1$  ms,  $B_{min} \gg 230$  nT


## How to switch between the two schemes

- Start up:
  1. Determine ambient field and direction using a field zeroing scheme
  2. Select scheme, microwave frequency, and direction of  $B_C$
  3. Begin operation.
- Continuous operation:
  1. Monitor signal size
  2. If signal size drops below threshold, switch scheme.
    1. If this fails, re-zero field.

## Quantization axis rotation



## ARP with longitudinal field



Need to understand better how to ramp  $B_C$  for ARP



## Atom shot noise limit

$$\delta B = \frac{1}{\gamma \sqrt{N T_2 \tau}}$$

$\gamma$  = gyromagnetic ratio  
 $T_2$  = transverse coherence time  
 $N$  = number of atoms  
 $\tau$  = measurement time

Decoherence limits noise to  $1/N^{1/2}$  scaling

Sensitivity improved by increasing  $T_2$  or  $N$

Most AMs do not operate at this limit.

## Photon shot noise limit

$$\delta B = \frac{\sqrt{2N_{ph}}}{dS(N_{ph})/dB} \propto \frac{1}{\sqrt{N_{ph}}}$$

Most AMs operate at or near the photon shot noise limit

Probe intensity contributes to  $T_2$  so  $N_{ph}$  cannot be made arbitrarily high.

Ground State Hamiltonian:

Spin angular momentum  $\mathbf{S}$

Nuclear angular momentum  $\mathbf{I}$

Total angular momentum  $\mathbf{F} = \mathbf{J} + \mathbf{I}$

$$H_{HF} = \mu_B g_S \mathbf{S} \cdot \mathbf{B} + \mu_B g_I \mathbf{I} \cdot \mathbf{B} + A_{HF} \mathbf{I} \cdot \mathbf{S}$$

$$g_S = 2.002, g_I = -0.000995 \text{ (} ^{87}\text{Rb), } -0.000399 \text{ (} ^{133}\text{Cs)}$$

$$E_{HF} = \frac{1}{2} A_{HF} (F(F+1) - I(I+1) - S(S+1))$$

$$^{87}\text{Rb} : \Delta E_{HF}/h = 6.835 \text{ GHz, } I = 3/2$$

$$^{133}\text{Cs} : \Delta E_{HF}/h = 9.192 \text{ GHz, } I = 7/2$$

At low field the nuclear and electron spins are combined,

$$H_B = \mu_B g_F \mathbf{F} \cdot \mathbf{B}, \Delta E_{mF} = g_F m_F \mu_B B$$

$$g_F \approx \pm \frac{1}{2F_{max}}, -F \leq m_F \leq F$$

$$\frac{g_S \mu_B}{h} = 2.8 \text{ MHz/G or } 28 \text{ Hz/nT (gyromagnetic ratio)}$$

$$^{87}\text{Rb} : g_F = \pm \frac{1}{4}, \frac{g_F \mu_B}{h} = 700 \text{ kHz/G or } 7 \text{ Hz/nT}$$

$$^{133}\text{Cs} : g_F = \pm \frac{1}{8}, \frac{g_F \mu_B}{h} = 350 \text{ kHz/G or } 3.5 \text{ Hz/nT}$$