

# A Dispersion Interferometer Diagnostic used for Low Electron Density Measurements in Magnetically Insulated Transmission Lines on Sandia's Z-Machine\*

Nathan Hines<sup>1</sup>, Sonal Patel<sup>2</sup>, Daniel J Scoglietti<sup>2</sup>, Robert H Dwyer<sup>1</sup>, Mark A Gilmore<sup>1</sup>, Darrell Armstrong<sup>2</sup>, David Bliss<sup>2</sup>, George Laity<sup>2</sup>, Michael E Cuneo<sup>2</sup>

1) The University of New Mexico, Albuquerque, NM 2) Sandia National Laboratories, Albuquerque, NM

\*Supported by the Laboratory Directed Research and Development program at Sandia National Laboratories, a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia LLC, a wholly owned subsidiary of Honeywell International Inc. for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525. This project was supported by LDRD project number 222428.

## Objective

- Reliable delivery of current to magnetically driven loads.
  - Electron flow can generate electrode plasmas and reduce coupling efficiency.
  - Increased electron flow further reduces coupling efficiency.
- A fiber-based Dispersion Interferometer (DI) will enable the first direct measurements of electron sheath flow on Z.
  - This will reduce the current lower limit for electrode plasma density measurements available by a factor of 100.
- This DI design operates with a Fundamental (F) wavelength at 1550 nm CW, with frequency-doubling to a Second-Harmonic (SH) wavelength at 775 nm.
  - The design will be fiber-coupled because of spatial limitations.




FIG 1: Diagram of basic DI design composed of the laser source, the first second-harmonic crystal (SHC), the plasma itself, the second second-harmonic crystal, a fundamental beam block, and a detector.

## UNM HelCat Experimental Setup




FIG 2.2: A 3-D render of UNM's HelCat device.

- (red) The DI path, traveling through plasma once.
- (black) Scanning double Langmuir probe.
- (blue) The mm-wave Interferometer path, traveling through plasma twice.

FIG 2.3: A 2-D axial cross section diagram of the three electron density diagnostics on UNM's HelCat device.

## DI Principles

- The measured change in phase of a wave with frequency  $\omega_\lambda$  will be related to the index of refraction of the plasma as:
$$\Delta\phi_\lambda = k_0 \int_0^l (N_\lambda - 1) dx \quad (1)$$
  - Where  $l$  is the plasma path length and...
$$N_\lambda = \sqrt{1 - \frac{p_e n_e e^2}{\epsilon_0 m_e \omega_\lambda^2}} \quad (2)$$
- For the DI, the total phase change is:
$$\Delta\phi_{total} = (2 * \Delta\phi_F) - \Delta\phi_{SH} \quad (3)$$
- Note: the first and second SH waves are perpendicularly polarized.
  - The reason for perpendicular polarizations of the 1<sup>st</sup> and 2<sup>nd</sup> SH waves is to split them into a differential detector while limiting the signal offset value.
- The information needed for an absolute electron density distribution is:
  - A pre-calibrated fringe sweep, defining the signal amplitude,  $A$ , and offset,  $DC$ .
  - A normalized electron density spatial profile of the plasma.
  - A measured phase change, from no plasma to plasma, between the two SH waves.
- With this, the total phase change is related to the measured signal as:
$$\Delta\phi_{total} = \sin^{-1}\left(\frac{V_0 + \Delta V - DC}{A}\right) - \sin^{-1}\left(\frac{V_0 - DC}{A}\right) \quad (4)$$
- This can be then compared to equation (3) to solve for the density.

## UNM HelCat DI Data




FIG 3.1: The "RF" signal sweep calibration data. Here we see the amplitude and DC-offset values we need.

Variables:  $A = 0.225 V$   $DC = 0.655 V$




FIG 3.2: The 15 shot average "RF" signal, with a red highlighted section representing the plasma operation period (220ms), a "higher RF" value (596.0-mV), a "lower RF" value (589.5-mV), and an orange mean curve.

Variables:  $\Delta V = 6.5 mV$   $V_0 = 0.59 V$




FIG 3.3: Derived absolute electron density of various radial profiles for the UNM HelCat machine, given the (above) signal data.

These absolute densities match expected results.

## Future Work

- Improve shot by shot variation and noise in the "RF" signal.
  - This is key for future experiments where shot count is limited.
- Collect HelCat profile density data with a scanning double Langmuir probe.
  - This will detail the normalized electron density radial profile shape inside the HelCat chamber.
  - Estimate the probe tip area to also get an absolute electron density distribution estimate, not just a normalized distribution shape.
- Collect HelCat mm-wave interferometer absolute line density data.
  - This will act as a verifying source against the DI diagnostic.
  - This will give an absolute electron density factor, with which to multiply into the scanning double Langmuir probe normalized distribution profile.
- Move the diagnostic to Sandia's Mykonos pulsed power machine and measure electron sheath flow in a moderate pulsed power environment.



- Testing the device capabilities at fast time scales (with the hope for time resolution) and in noisy environments.
- The device will shift to primarily fiber-coupled, from primarily open-beam on UNM's HelCat device.
- Employ the diagnostic at the Z machine.
  - Directly measure the electron densities in magnetically insulated flows.
- Develop a multi-channel DI design for spatial resolution capabilities.

## References

1. Lynn AG, "The HelCat dual-source plasma device," Rev Sci Instrum. Vol. 80, pp. 103501, 2009
2. V. Licht, "A sensitive dispersion interferometer with high temporal resolution for electron density measurements," Rev Sci Instrum, vol. 71, no. 7, pp. 2710-2715, 2000.
3. M. A. Van Zeeland, "Fiber optic two-color vibration compensated interferometer for plasma density Measurements," Rev Sci Instrum, vol. 77, pp. 10F325, 2006.