



The Development of an Optical Centrifuge for Chemical Dynamics Studies as Probed through Time-Resolved Coherent Raman Scattering

Christopher J. Klierer
SPIE Optics and Photonics
August 21, 2022



Optical centrifuge first demonstrated through rotational dissociation of Cl_2

VOLUME 85, NUMBER 3

PHYSICAL REVIEW LETTERS

17 JULY 2000

Forced Molecular Rotation in an Optical Centrifuge



D. M. Villeneuve,^{1,*} S. A. Aseyev,¹ P. Dietrich,^{1,2} M. Spanner,¹ M. Yu. Ivanov,¹ and P. B. Corkum¹

¹National Research Council of Canada, 100 Sussex Drive, Ottawa, Ontario, Canada K1A 0R6


²Freie Universität Berlin, Arnimallee 14, D-14195 Berlin, Germany

(Received 17 February 2000)

- First experimental demonstration by Corkum and coworkers in 2000.
- Rotational dissociation of Cl_2 observed through TOF-MS




Subsequent measurements demonstrated “super-rotor” reduction in collisional energy transfer rate






- Mullin and coworkers have made measurements using transient IR absorption for centrifuged CO₂, CO, and N₂O.
- Provides for the first time the ability to both centrifuge and study the energy transfer dynamics of molecules in high angular momentum states.

Two-Beam Hybrid Femtosecond/Picosecond RCARS




Phase-matching Condition



$$\mathbf{k}_{\text{CARS}} = \mathbf{k}_{\text{pump}} - \mathbf{k}_{\text{Stokes}} + \mathbf{k}_{\text{probe}}$$



Home-built Optical Centrifuge and Coherent Raman Probe




Home-built Optical Centrifuge and Coherent Raman Probe




Demonstrated Centrifuge of CO₂ to J=250, E_{rot} = 24000 cm⁻¹

Thermal Distribution at 300 K




- Chopping one of the centrifuge arms in the stretcher allows control over the final J-state distribution
- Can be limited to thermal distribution for comparison

Linewidths in the centrifuge for N₂ are significantly larger than at 294 K



- High-finesse coherence decays have been recorded in our lab for extreme accuracy in molecular thermometry
- Oriented angular momentum appears to modify collisional energy exchange rate


Comparison of collisional dephasing rates for CO₂ reveals significant difference for centrifuged molecules





- Dramatic enhancement of collisional dephasing rate observed for centrifuged CO₂ molecules.
- Comparison to earlier time-domain measured decoherence rates using fs/ps CARS reveals more than a factor of 2X faster dephasing.
- Comparison to thermally limited ensemble still results in enhanced decay rates.

Time evolution of thermally limited centrifuge pulse reveals growth of vibrational excitation

- Fits well to the rotational shift expected from $v=1$ CO_2 bending mode.
- Interesting because molecules are limited to thermal distribution
- Special collisional energy exchange due to oriented angular momentum?
- Resolution of vibrational excitation limited by standard grating-spectrometer solution




VIPA-based CARS spectrometer allows for spectral resolution down to ~ 100 MHz



- The vibrational population of N_2 during plasma discharge is assessed from the shift in the pure-rotational spectrum

VIPA-based CARS spectrometer isolates signal from ground and excited vibrational states



Future Directions

- Chemical dynamics studies coupled to velocity mapped imaging apparatus
- Energy exchange with surfaces where rotational energy is sufficient to drive surface reactions.
- Pulse shaping approaches to selectively spin up single rotational axes within polyatomic molecules.



Conclusions

- We have built and demonstrated a molecular optical centrifuge (3rd group in world to do so), coupled to a coherent Raman probe to assess the time-domain evolution of the excited rotational coherence.
- Surprising changes to coherence lifetimes observed for molecules for which the sense of rotation is controlled. Indicates modification to efficacy of collisional energy exchange.
- VIPA-based coherent Raman spectrometer allows high-spectral resolution analysis of pure-rotational scattering, and the separation of excited vibrational populations.

Thank you for your time!



Acknowledgements

Sandia Team:

Timothy Chen (Lead postdoc)

Brian Patterson

Scott Steinmetz

Dave Chandler

External

Collaborators/Discussions:

Ahren Jasper (Argonne)

Funding: BES GPCP

Happy to take your questions.

This material is based upon work supported by the Division of Chemical Sciences, Geosciences and Biosciences, Office of Basic Energy Sciences (BES), U.S. Department of Energy (USDOE). Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC., a wholly owned subsidiary of Honeywell International, Inc., for the USDOE's National Nuclear Security Administration under contract DE-NA0003525. This paper describes objective technical results and analysis. Any subjective views or opinions that might be expressed in the paper do not necessarily represent the views of the USDOE or the United States Government.

