



Electrochemical Aptamer Based Sensors for Measurements In Vivo and in Undiluted Blood

Alex Downs¹, Julian Gerson⁴, M. Nur Hossain⁵, Kaylyn Leung⁶, Kyle Ploense⁴, Michael Pham⁶,
Heinz-Bernhard Kratz⁵, Tod Kippin⁴, Phillipe Dauphin-Ducharme⁷, Kevin Plaxco^{3,6}

¹Sandia National Laboratories ²UC Santa Barbara Department of Mechanical Engineering, ³UC Santa Barbara Center for Bioengineering, ⁴UC Santa Barbara Department of Psychological and Brain Sciences, ⁵University of Toronto Scarborough Department of Physical and Environmental Sciences, ⁶UC Santa Barbara Department of Chemistry and Biochemistry, ⁷University of Sherbrooke Department of Chemistry.

Electrochemical Aptamer Based Sensors

DNA aptamers are artificially selected to bind to a target, such as a small molecule, protein, or metabolite





These sensors:

- **Function both in whole blood and in the living body**
- Enable real-time measurements of rising and falling target concentrations
- Are generalizable to a range of different targets


Sensor Miniaturization Using Nanoporous Gold²

Electrochemically-deposited nanoporous gold increases electrode surface area by 75-100X. This increase yields a 20-30X increase in sensor redox signal.



- **The nanoporous morphology enables miniaturization** of in-vein working electrode by 6-fold, while the dendrite morphology typically used in vivo does not.
- Here, we demonstrate in vivo dosing and measurement of the antibiotic, vancomycin, using square wave voltammetry

Sub-second Measurements Using Impedance¹






- When interrogating a sensor at its redox potential, a peak occurs in the phase component of impedance
- With target addition the peak position moves toward higher frequencies.
- Here, we use phase changes to measure the antibiotic, tobramycin.

Measuring phase at specific frequencies yields calibration curves.

Depending on frequency, the response:


- Increases in phase toward 0 degrees (Signal-on) or
- Decreases in phase angle toward -90 degrees (Signal-off)



Collecting phase at specific frequencies near the phase peak **achieves sub-second temporal resolution**


Phase interrogation yields stable measurements in whole blood for long time durations

Measurement Accuracy in Whole Blood³



- Temperature influences the response of electrochemical aptamer based sensors by changing the rate of electron transfer.
- Using **freshly collected, body temperature blood improves the sensor calibration** accuracy to greater than 10% in the clinical range.

Microneedle-based Measurements in Whole Blood



- The sensors, when implemented in a **microneedle format**, achieve measurement of vancomycin in body temperature whole blood

Citations and Acknowledgements

- [1] A. Downs, J. Gerson, K.L. Ploense, et al. *Analytical Chemistry* 92 (20), 14063-14068. 2020.
- [2] A. Downs, J. Gerson, M.N. Hossain, et al. *ACS Sensors* 6 (6) 2299-2306. 2021.
- [3] A. Downs, J. Gerson, K. Leung, et al. *Scientific Reports* 12, 5535, 2022.

This work was supported by the National Science Foundation Graduate Research Fellowship (1650114)¹⁻³, the Office of Naval Research (N00014-20-1-2164)^{2,3}, the National Institute of Health (R01AI145206, EB022015)^{2,3}, the Natural Sciences and Engineering Research Council of Canada (NSERC-RGPIN-2016-06122)³, and Sandia's Laboratory Directed Research and Development (LDRD) program. The content of the information does not necessarily reflect the position or the policy of the Government, and no official endorsement should be inferred.

Sandia National Laboratories is a multimission laboratory managed and operated by Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.