



Electrochemical aptamer-based sensors for measurements in vivo and in undiluted blood




Alex Downs

September 10th, 2022




Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.


My academic path




Faster Measurements



Smaller Sensors




More Accurate Quantification



Wearable sensors can revolutionize medical care, but we have them for very few biomarkers



Continuous glucose meter




Heo, Y., Kim, S. *Applied Sciences*. 2019.

Relies on enzymatic detection, which has not proven widely generalizable to other target

Electrochemical aptamer-based sensors ("EAB Sensors") offer a promising route for wearable sensing of many biomarkers




Aptamers: DNA, RNA, or peptide sequence artificially selected to bind specifically to a **small molecule, protein, or cell**






Wolter, Mayer. J. Neurosci.. 2017.


Redox
reporter



Target
binding






What makes this sensing platform unique?



- No washing steps - works directly at the point of detection
- Quantifies increasing and decreasing concentrations
- Works in undiluted bodily fluids, living body


Goal 1: Miniaturization of in vivo probes




Courtesy of
Kyle Ploense, Ph.D.



Approach: Increase surface area





Electrochemical alloying/dealloying





Downs, A., Gerson, J., Hossain, M., Ploense, K., Pham, M., Kraatz, H.B., Kippin, T., Plaxco, K., "Nanoporous Gold for the Miniaturization of In-Vivo Electrochemical Aptamer-Based Sensors." *ACS Sensors*. 2021, 6, 6, 2299–2306.

Nanoporous gold increases surface area and electrochemical signal



Nanoporous gold enables 6X miniaturization




Downs, A., Gerson, J., Hossain, M., Ploense, K., Pham, M., Kraatz, H.B., Kippin, T., Plaxco, K., "Nanoporous Gold for the Miniaturization of In-Vivo Electrochemical Aptamer-Based Sensors." *ACS Sensors*. 2021, 6, 6, 2299–2306.


Goal 2: Improve calibration accuracy




Temperature influences calibration curve



37°C data quantified with 25°C calibration curve





37°C data quantified with 37°C calibration curve






Downs, A., Gerson, J., Honeywell, K., Kippin, T., Plaxco, K., "Improved Calibration of Electrochemical Aptamer-Based Sensors."
Under submission to Scientific Reports.

Developing a calibration curve: media matters



Takeaways



- Nanoporous gold can enable miniaturization of EAB sensors for in vivo applications
- Calibrating measurements in the body, however, requires careful consideration of temperature and media
 - Matching temperature of calibration achieves most accurate measurements



Current goal at Sandia National Laboratories: Adapt aptamer sensors for microneedle-based measurements



30 μM dose



Thank you for your attention

Questions?



Acknowledgements:

- Plaxco, Kippin, and Kraatz group
- Sandia Chemical and Biological Sensors Department

Graduate School Funding: NSF GRFP, UCSB Graduate Division, Institute for Collaborative Biotechnologies

Postdoctoral Funding: Sandia National Laboratories Hruby Fellowship

If you're interested in working at a national lab or applying to postdoctoral fellowships there, reach out on LinkedIn or at amdowns@sandia.gov