

Application of Peridynamics to Predict Pharmaceutical Tablet Robustness

Sean Garner¹, William Ketterhagen¹, John Strong¹, Stewart Silling²

¹ AbbVie Development Sciences, Abbott Park, IL

² Sandia National Laboratories, Albuquerque, NM

AIChE Annual Meeting

November 15, 2022

SG, WK, and JS are employees of AbbVie and may own AbbVie stock. AbbVie sponsored and funded the study; contributed to the design; participated in the collection, analysis, and interpretation of data, and in writing, reviewing, and approval of the final presentation. SS is an employee of Sandia National Laboratories.

abbvie





Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia LLC, a wholly owned subsidiary of Honeywell International Inc. for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

Sandia National Laboratories is a multimission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

Introduction

- Tablets are a common oral solid dosage form
- Much is understood about the powder compaction process
 - Decades of manufacturing experience
 - Lab-based tools (e.g. compaction simulators)
 - Computational tools (e.g. FEM simulation)
- Despite the extensive experience with compaction, some issues are still observed in development and commercial tableting operations



Picking (Red) / Sticking (Blue)




Simmons and Gierer, (2011)

Capping



Wu et al., (2007)

Lamination




Anuar and Briscoe, (2010)

Chipping



Wu et al., (2007)

Erosion



Sinka et al., (2004)


Objectives

- Develop a peridynamics-based simulation approach to predict tablet robustness and fracture
 - Focus on unloading and ejection stages of the compaction process
- Predict tablet robustness (modes of tablet breakage) for
 - varying process conditions,
 - tooling geometries, and
 - material properties (as characterized by Drucker-Prager Cap model)
- Objectives aim to improve predictive capabilities for tablet failure, provide insight into breakage modes, and guide improvements to formulation/process.




Background – What is Peridynamics?

- A generalization of the standard theory of solid mechanics that allows fracture within its basic equations
 - Seamlessly models the transition from continuous deformation to fracture.



Peridynamic
simulation




Metallic glass crack tip
Images: Hofmann et al, 2008

First paper on peridynamics: Silling SA. Reformulation of elasticity theory for discontinuities and long-range forces. *Journal of the Mechanics and Physics of Solids*. 2000 Jan 1;48(1):175-209.

Reference book: Madenci E, Oterkus E. *Peridynamic theory and its applications* 2014. Springer, New York, NY.


Background – Peridynamic Approach

- The standard partial differential equations don't allow discontinuities such as cracks.
- Integral equations can be applied directly on discontinuities.




$$\begin{aligned}\frac{\partial \sigma_{xx}}{\partial x} + \frac{\partial \sigma_{xy}}{\partial y} + b_x &= 0 \\ \frac{\partial \sigma_{yx}}{\partial x} + \frac{\partial \sigma_{yy}}{\partial y} + b_y &= 0\end{aligned}$$

PDEs of standard solid mechanics



Peridynamic integral equation




$$\int_{\mathcal{H}} f(q, x) dq + b(x) = 0$$


- The material point x interacts directly with neighbors q within its horizon (cutoff distance).
- The bond between x and q has a bond force density $f(q, x)$ (force/volume²).
- $b(x)$ is the external body force density.
- \mathcal{H} is the family of x (material points within the horizon).
- The integral sums up all these forces acting on x .

Background – Peridynamic Material Model

- The material model determines $\mathbf{f}(\mathbf{q}, \mathbf{x})$ for every \mathbf{q} in \mathcal{H} , for every possible deformation of \mathcal{H} .
- Any material model from the standard theory can be used.





- Bonds can break irreversibly.
- After breakage, a bond cannot carry any tensile load.
- The criterion for bond breakage can be anything you can dream up.
- The simplest criterion is **critical bond strain**.
 - This criterion is used in the tablet simulations described later, but with dependence of the critical strain on pressure and solid volume fraction.




Autonomous Crack Growth

- Bonds break whenever they feel like it.
- When a bond breaks, it becomes more likely that a neighboring bond will also break.




Peridynamic Material Model for Tablet Breakage

- Bond force is linear in strain until failure.
- Bonds fail when their strain ϵ exceeds a critical strain ϵ_c .
- The critical strain depends on:
 - Hydrostatic pressure P .
 - Solid volume fraction (SDV).
- Parameters are calibrated from diametral compression & axial compression tests.




Bond failure dependence on solid volume fraction



Bond failure dependence on pressure

Future Workflow Envisioned




- Material-sparing, predictive approach to guide tablet development
- More responsive troubleshooting of compaction issues
 - The performance of each tooling iteration can be assessed *in silico*

Workflow for Model Development

Lab-based
Experiments

Computational

Calibration




Workflow for Model Development


Lab-based
Experiments

Computational

Calibration



Validation



Calibration Inputs

- Formulation

- 60% Acetaminophen (APAP)
- 39% Avicel PH102 (MCC)
- 1% Silicon Dioxide MP5

- Geometry

- 10 *mm* diameter flat-faced tooling
- Straight wall die

- Process Conditions – Compacts for Axial Compression Testing

- Fill weights: ~300 mg
- Compaction pressures: ~200 MPa
- Thickness: ~3.50 mm (fixed)

Resulting Solid Fractions:

In-die (SDV): 0.97
Out-of-die: 0.86

- Process Conditions – Compacts for Diametric Compression Testing

- Fill weights: ~750 mg
- Compaction pressures: ~230 MPa
- Thickness: ~7.7 mm (fixed)


Resulting Solid Fractions:

In-die (SDV): 0.97
Out-of-die: 0.89

Formulation
(composition,
material
properties),
Geometry
(tooling size/shape,
die taper, etc.)
Process
(compaction speed,
dwell time,
unloading speed)



Compaction Simulation




Compaction Simulation



Fette 2090i
Tablet Press

- Fette 2090i tablet press was simulated using a Huxley-Bertram compaction simulator
- All tablets produced in this study used a similar press profile with only changes in main compression gaps – turret speed and dwell time were the same for all tablets produced

Compaction Simulation for DPC Model Calibration



- Describes the limit to elastic deformation
- One curve per relative density (denser=stronger)
- Calibrated easily by press simulator experiments

1 Diametral and Simple Compression Strength Tests

- ❖ Two tests per density level





2 Fully Instrumented Die Compaction Experiments

- ❖ Only one experiment required for all densities







Compaction Simulation – DPC Shear Failure Surface Calibration

- Following Cunningham, et al., (2004) and Han, et al., (2008)



Compaction Simulation – DPC Cap Surface Calibration

- Following Cunningham, et al., (2004) and Han, et al., (2008)





Calibration Experiments


Axial &
Diametric
Compression
Testing



Axial Compression




Diametric Compression





FEM Simulation

- FEM simulations of powder compression
- The final compressed, in-die state becomes the input to the peridynamics simulation



FEM Simulation



Peridynamics Simulation – Overview

PD Simulation
of Ejection and
Axial/Diam.
Compression


- Inputs for PD simulation
 - DPC material model parameters from compaction simulation
 - In-die tablet stress and solid fraction distribution from FEM
 - Compression conditions as desired
- PD simulates same process as used in the experiments:
 - (1) Ejection of tablet from die
 - (2) Axial or diametric compression test



Peridynamics Simulation – Axial Compression


- Simulated mechanical testing ($SDV \approx 0.97$)

PD Simulation
of Ejection and
Axial/Diam.
Compression



Initial State



Start of failure



Peridynamics Simulation – Diametric Compression

PD Simulation
of Ejection and
Axial/Diam.
Compression

- Simulated mechanical testing ($SDV \approx 0.97$)




Drucker-Prager failure model during ejection



Tablet "capping" failure
Image: merlin-pc.com

Ejection of a pharmaceutical tablet from a rotary press




Workflow for Model Development


Lab-based
Experiments

Computational

Calibration




Validation



Validation Inputs

- Formulation [same as before]
 - 60% Acetaminophen (APAP)
 - 39% Avicel PH102 (MCC)
 - 1% Silicon Dioxide MP5
- Geometry [different tablet shape]
 - 10 mm x 17.3 mm oblong tooling
 - Straight wall die
- Process Conditions [different]
 - Fill weight: ~900 mg
 - Compaction pressure: ~200 MPa
 - Thickness: ~6.5 mm








Formulation
(60% APAP)
Geometry
(oblong tooling)
Process
(baseline
compaction
conditions)





FEM Simulation

- FEM simulations of the compaction (only) of the oblong tablet
- DPC material model parameters as before







SDV



Peridynamics Simulation of Ejection

Peridynamics
Simulation of
Ejection




- The PD simulation of ejection is initialized from an Abaqus simulation of compression.
- The initial pressure and solid fraction vary significantly with position.

Experimental Observations – Oblong Tablet

Tablet Press
(or Compaction
Simulation)




- Post-ejection μ CT
- Significant damage occurs during ejection



Peridynamics Simulation Prediction of Ejection



Peridynamics
Simulation of
Ejection





- For this tablet, the ejection simulation predicts a complete fracture and additional internal damage.
- Drucker-Prager condition is used as a nucleation condition for damage.
- Fracture surface is rough.


Experimental Observations – Oblong Tablet

- Mechanical Testing


Tablet Press
(or Compaction
Simulation)




Mechanical Testing Peridynamics Simulation Prediction



Start of diametral compression
Colors show damage



Post-failure
Colors show damage




- The post-ejection tablet model, including damage, is then subjected to compression.
- Failure load is ~52 N.

Summary & Future Work

- The use of peridynamics for prediction of tablet breakage looks promising
 - Material-sparing, predictive approach to guide tablet development and troubleshooting
- Presented a general workflow for calibration and validation
 - Quantitative agreement on breaking force
 - Qualitative agreement on damage / breakage pattern
 - Working on refinements for improved predictions

Peridynamics prediction of tablet capping



Good agreement with experimental observations of Wu et al. (2008).




abbvie

Drucker-Prager failure model during ejection




Tablet "capping" failure
Image: merlin-pc.com

Ejection of a pharmaceutical tablet from a rotary press




Fun applications


VIDEOS




Unstable crack path in a polyethylene membrane




Peeling of tape



Grinding
(includes friction in short-range forces)



EMU model of a balloon penetrated by a fragment




Early high speed photograph by Harold Edgerton
(MIT collection)


<http://mit.edu/6.933/www/Fall2000/edgerton/edgerton.ppt>




Mechanics Models for Powders in terms of *Continuum Mechanics*

Pressure Dependence of Mechanical Response



Tablet Geometry





10.0

10.0

10.0