






High-speed measurements in the Sandia free-piston shock tunnel



Cylinder startup NO PLIF



Sean P. Kearney
Engineering Sciences Center
Sandia National Laboratories
Albuquerque, NM 87175
spkearn@sandia.gov

Contributions from: Justin Wagner, Kyle Lynch, Elijah Jans, Josh Hargis, Kyle Daniel, Raj Bhakta, Charley Downing (Sandia)


Arlee Smith (AS Photonics)

Philip Varghese, Noel Clemens, John Murray, Spenser Stark, Dan Fries (Texas)



Sandia National Laboratories is a multimission laboratory managed and operated by National Technology & Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International Inc., for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-NA0003525.

Sandia Free-Piston Hypersonic Shock Tunnel (HST)



Tunnel Specifications


- Nozzle Exit Dia. = 0.36 m
- Test section diameter 0.5 m
- Run times of 1-2 milliseconds

U_{∞} (m/s)	H_0 (MJ/kg)	T_0 (K)	P_0 (MPa)
2850	4.6	3400	12
4060	9	6000	17

Target applications include high-temperature surface chemistry and hypersonic thermochemistry.



Survey of Upcoming Experiments in HST

- **HST introduces flow complexities**
 - Stagnation region gases react
 - Gas rapidly expanded through nozzle
 - Result: thermal non-eq., NO addition
- **Free-stream characterization necessary**
 - Temperature: CARS for heteronuclear molecules
 - Velocity: NO LIF
 - 100-kHz data with pulse-burst laser!
- **Examine boundary layer products**
 - Speciation/temperature of CO
 - Laser absorption
 - CARS (Coherent Anti-Raman Stokes Raman Scattering)




Pulse-burst laser for 100-kHz laser diagnostics


- Burst-mode lasers have allowed experimentalists to access high-speeds (10s to 100s of kHz)
- While powerful, these systems are not wavelength tunable—this prohibits application of **chemically specific** imaging and spectroscopic tools



Picosecond Optical Parametric Generation




Optical Parametric Oscillator (nanosecond)





UV Generation for NO PLIF


Coherent anti-Stokes Raman scattering (CARS)



Coherent anti-Stokes Raman scattering (CARS)




Calculated N₂ CARS Spectra





P. Danehy and S. Tedder
NASA Langley Research Center

Coherent Anti-Stokes Raman



Free-Stream Characterization: Temperature



- **Free-stream conditions**
 - Major source of uncertainty in shock tunnels
 - Temperature non-eq. in nozzle is expected
- **Simulation of nozzle temperatures**
 - Significant T_v differences between species
 - N_2 has highest degree of non-eq
- **Characterizing temperature non-eq. in HST**
 - Use CARS to measure T_{vib} , T_{rot} for N_2
 - Further improvement needed for T_{rot}
 - Next: O_2 CARS temp. measurements



100-kHz *nanosecond* CARS for HST free-stream




Noncolinear OPO produces wide spectral bandwidth for CARS detection

$$FWHM_{OPO} = \frac{c}{l_c |\nu_s - \nu_l|}$$




- 2X reduction in measurement uncertainty re: picosecond CARS
- Good sensitivity to thermodynamic nonequilibrium
- Improve sensitivity to T_{rot}
- O_2 measurements


CARS spectra can reveal nonequilibrium conditions in N_2 and O_2



100-kHz nitric oxide LIF imaging of transient flows




- $M = 9$ flight condition at ~ 130 kft
- Enthalpy ~ 5 MJ/kg
- NO PLIF Visualization over large, 70-mm field of view



Shock Tube Experiment

- $T = 3000$ K
- $U = 2.5$ km/s
- NO concentration $\sim 4\%$






100-kHz NO PLIF Imaging of Cylinder wake startup: $U = 2.5$ km/s, $T = 3000$ K

Free-Stream Characterization: Velocity

- NO is present in shock tunnel flow ($X_{\text{NO}} \sim 4\text{-}5\%$)
- Tracer for flow visualization
- Nitric Oxide Tagging Velocimetry
- Long fluorescence lifetime, $>100\text{ ns}$
 - $U_{\infty} = 3\text{ km/s} = 3\text{ }\mu\text{m/ns}$, $\Delta t \sim 100\text{ ns} \rightarrow \Delta x \sim 300\text{ }\mu\text{m}$
 - Track NO fluorescence at high image magnification

LIF beam tracks flow



TEMPERATURE/SPECIES MEASUREMENTS IN TPS BOUNDARY LAYERS

Shock Tunnel—Resistively Heated Models



- Impulse facility – insufficient test time
- Graphite coupons as TPS surrogate




- Absorption (QCL)
- Pulse-burst CARS
- CO PLIF

IC Plasma Torch Environment



- Long test duration (minutes)
- Low velocity but relevant temperature



Simultaneous T/CO measurement






50- μm \times 3-mm CARS volume



Summary and Conclusions



- Sandia relies heavily on laser diagnostics for high-speed measurements in its free-piston hypersonic facility
 - Short-duration, impulsive experiments
 - Thermal and chemical nonequilibrium
 - Free-stream and near-surface data





- Pulse-burst lasers can be adapted for CARS thermometry!
- Pulse-burst PLIF has been applied for flowfield imaging and velocimetry – 4-5% NO!

- Complementary measurement in ICP torch environment

- Additional laser-based methods include:
 - Laser absorption: QCL, VCSEL
 - Emission spectroscopy



100-kHz Pulse-Burst CARS in the Sandia Free-Piston Shock Tube



100-kHz Pulse-Burst CARS in the Sandia Free-Piston Shock Tube

