

6th Diesel Engine Emissions Reduction Workshop



Exhaust Aftertreatment, Gas Phase Measurements *in situ*

Neal W. Currier*, Michael J. Cunningham,
Cummins Inc.

William P. Partridge, John M. E. Storey, Samuel A. Lewis,
Robert W. Smithwick, Gerald L. DeVault
Oak Ridge National Laboratory

*neal.w.currier@cummins.com, (812)377-4625



Challenges of New Aftertreatment Technology



Increasing aftertreatment system complexity

- Driven by more stringent regulations
- Results in active aftertreatment systems
- Produce rapid dynamic events

Lack of fundamental understanding

- Need for engineering data
- Laboratory/on-engine differences



Measurement Needs



Critical Requirements:

- Catalyst surface temperature
- Catalyst surface composition
- **Exhaust gas composition**
- **Fast response times - $T_{90} \ll 1$ second**
- **Measure inside operating device**



Current Measurement Limitations



Available Gas Analysis Instrumentation

- Responds slowly, 1- 10 sec. T_{90}
- Samples large volume, 0.1- 5 l/min.
- Measures single component
- Has long sampling delays
- Requires sample pretreatment
- Do not measure absolute pressure



SpaciMS



Requirements


- Measure NO_x, O₂, CO₂, & HC simultaneously
- Sample at high speed - $T_{90} < 0.1$ seconds
- Sample inside catalyst monolith
- Disturb gas flow minimally

Spatially resolved, capillary inlet mass spectrometer






SpaciMS Measured NOx




DEER 2000 August 23, 2000

Combining Spatial and Temporal




DEER 2000 August 23, 2000

NOx Concentration Profiles



NOx Concentration Profiles



Advantages of SpaciMS



Resolve catalyst behavior spatially & temporally

- Measure anywhere in the catalyst device
- Follow fast events
- Provide temporal relationship between events

Investigate Chemistry

- Absorption, desorption/conversion rates
- concentration effects

Predict effects of catalyst size and shape

Diagnose failures

